

EASYPULSE

Portable Oxygen Concentrators

The Difference is Precision

The EasyPulse POCs are the perfect products for patients that demand quality and reliability.

VPSA technology keeps sieve beds cleaner for maximum product life.

Proven design and durability that supports your non-delivery business.

Lowest cost of ownership.

Best in industry warranty including sieve bed replacements.

When you place a Precision Medical POC with a patient, you can rely on our durable, proven design to support your business.

Precision Medical POCs

Precision Medical POCs can help your business by eliminating the high recurring costs of deliveries and the burden of managing cylinders. Other manufacturers use basic pressure only systems causing frequent sieve bed replacements. With our dependable POCs you can expect to minimize patient calls and enjoy lower cost of ownership.

Innovative VPSA technology with CMV

Only Precision Medical POCs use Vacuum Pressure Swing Adsorption technology ensuring each sieve bed is perfectly regenerated every cycle. This results in clean sieve beds that produce high oxygen purity for a long period of time.

With controlled minute volume (CMV), your oxygen patients receive therapeutic oxygen in the first third of the inspiratory phase for maximum efficiency, regardless of breath rate.

Easy to use

An easy to use push button operation allows patients to press one button to start their POC on the prescribed flow setting. Simple design reduces phone calls and troubleshooting.

Proven battery life

The internal and external batteries combine to provide up to 11 hours* of battery life. This duration allows patients to be mobile with the confidence that their oxygen needs will be met, at home as well as when traveling. CMV technology insures battery duration stays consistent on all flow settings.

Recognized design with longest product life

EasyPulse POCs' product design is recognized by providers not only in the USA, but around the world for being the most dependable POC in their fleet. With the lowest return rate of any POC, trust the company that has been a leader in hospital respiratory care since 1982.

*Battery duration times are based on new, fully charged batteries used with the EasyPulse POC system. Battery duration times will decline with battery age and environmental operating conditions.



Ordering Information

Device	3 L (PM4130)	5 L (PM4150)
EasyPulse POC Oxygen Concentrator includes: AC and DC power supplies, cords, carry bag and cannula	PM4130P	PM4150P
EasyPulse POC Oxygen Concentrator includes: AC and DC power supplies, cords, carry bag, cannula and external battery	PM4130PEB	PM4150PEB
EasyPulse POC Oxygen Concentrator includes: AC and DC power supplies, cords, backpack, cannula and external battery	PM4130PEBK	PM4150PEBK

Replacement accessories

Description	Part Number
Filter (Pack of 30)	505847-30
External Battery	506760
AC Power Supply	506629
AC Power Cord	506961
DC Car Charger	506750

Description	Part Number
Cannula	504833
Filter	505847
Carry Bag (PM4150 only)	506707
Back Pack (PM4150 only)	507303
Carry Bag (PM4130 only)	507947
Back Pack (PM4130 only)	508124

Specifications

Oxygen concentrations	Up to 95.5%
User interface	Push Buttons
Warranty	Device: 5 years Sieve Beds: 5 years Battery: 1 year
Size PM4130	PM4150
Height: 8.5 in (21.6 cm) Width: 6.5 in (16.5 cm Depth: 4.5 in (11.4 cm)	10.1 in (25.7 cm) 6.5 in (16.5 cm) 4.5 in (11.4 cm)

Specifications are subject to change without prior notice.

PM4130	PM4150	
Battery durations (at any breath rate)		
Setting 1: 5.5 hrs Setting 2: 3.8 hrs Setting 3: 2.8 hrs	Setting 1: 4.7 hrs Setting 2: 3.4 hrs Setting 3: 2.7 hrs Setting 4: 2.0 hrs Setting 5: 1.6 hrs	
2 batteries: up to 11 hours	2 batteries: 9.4 hrs	

Weight 4.9 lbs (2.22 kg) 6.6 lbs (2.98 kg)

The EasyPulse POCs features:

Lightweight, durable and proven designs

3L and 5L pulse settings utilizing VPSA technology

With CMV, patients receive therapeutic oxygen in the first third of inspiration

Simple, one touch operation

Fully functional on AC or DC power

FAA approved for use on commercial airlines

Optional backpack making it comfortable for patients to wear all day long

Best in industry 5-year warranty covering sieve bed replacements

Lowest cost of ownership

Battery packs with up to 9.4 hours with EasyPulse POC-5 and 11 hours for EasyPulse POC-3



For more information about the EasyPulse POC, call customer support at 1-800-272-7285

Email us at: info@precisionmedical.com





