



## FALLS MANAGEMENT AND PATIENT SAFETY CATALOG



The right products, best practices and safety essentials.

# Elevate patient safety.


Help reduce falls with targeted solutions.

Patients may be at risk of falling. The challenge is how do you identify and protect the patients most at risk while keeping everyone safe—and on their feet? Medline answers the call. We offer hundreds of **quality products and a unique falls management program** to help drive best practices and reduce patient falls.

**Strengthen your Falls Management with smart solutions** designed to help elevate patient safety and improve patient care. From falls risk recognition to safety essentials, we have everything you need to empower your team.

And your bottom line.

FALL RISK



## FALLS MANAGEMENT PROGRAM

4 Falls Management Program

## FALLS RISK RECOGNITION

6 Falls Prevention Slippers

8 Falls Prevention Gowns

9 Falls Prevention Kits

## TARGETED INTERVENTIONS

10 Confusion

14 Shower Furniture

15 Toileting

19 Bed Safety

## SAFETY ESSENTIALS

24 Alarms, Pads

28 Patient Monitoring System

30 Alarms

31 Gait Belts

32 Hand Protectors, Mitts, Splints


34 Limb Holders

37 Repositioning Belts

38 Safety Vests, Jackets, Wraps

39 Cones, Heel Protectors, Seat





“Patients deserve to be safe and clinicians deserve to be supported. We’re here to care for everyone.”

— Jay Hamilton-Roque BSN, RN

# Fight falls with our proven products. And powerful program.

Are you doing everything you can to protect your patients? The right products are critical to falls protection. But it's a complex issue that also demands expert support. That's why we've developed **Stand FIRM™—Falls Intervention and Risk Management**.

Developed by clinicians for clinicians, Stand FIRM is a simple, four-step program that makes falls management easy—and effective. We give you the practical tools and ongoing support you need to help protect patients and drive lasting change.



# STAND FIRM



Falls Intervention  
& Risk Management

## Q&A

**Q. Why aren't falls products enough to help protect patients from falls?**

**A.** Products are a key component of any falls management program. They are essential in supporting a strong clinical program to help keep patients safe.

**Q. We already have a falls management program. Why change?**

**A.** Stand FIRM is a comprehensive program and can complement your existing program, allowing you to assess patients from every angle and support your clinicians every step of the way. Your staff will be empowered to deliver the highest level of falls management care and you see lasting, meaningful change.

**Q. How easy is Stand FIRM to use?**

**A.** We've simplified our program into four easy-to-follow phases, carefully guided by us, every step of the way. We work with you to tailor Stand FIRM to your exact needs and seamlessly integrate it into your practice.

**Q. How does Medline continue to support us with Stand FIRM?**

**A.** Adopting the program is just the beginning. We partner with you for the long-term to evolve Stand FIRM as your needs change.

**Let's Stand FIRM against falls.**  
Contact your Medline Sales Representative  
today to schedule an onsite assessment.



## Falls Risk Recognition

**Best practices lead to better outcomes.**

Make it easy to do the right thing. Arm your team with the tools they need to quickly identify at-risk patients.



## Falls Prevention Slippers

- » Bright colored slippers help identify at-risk patients
- » Double-sided, non-slip treads help patients stay on their feet
- » Soft, breathable terrycloth interior keeps feet comfortable, dry and warm

### Yellow – Moderate Risk

| Item No.      | Size      | Tread Color | Pkg.     |
|---------------|-----------|-------------|----------|
| MDTPS4B01FPY  | Child     | White       | 48 pr/cs |
| MDTPS4B01FPYH | Child     | White       | ea       |
| MDTPS4B03FPY  | Medium    | White       | 48 pr/cs |
| MDTPS4B03FPYH | Medium    | White       | ea       |
| MDTPS4B04FPY  | Large     | White       | 48 pr/cs |
| MDTPS4B04FPYH | Large     | White       | ea       |
| MDTPS4B05FPY  | X-Large   | White       | 48 pr/cs |
| MDTPS4B05FPYH | X-Large   | White       | ea       |
| MDTPS4B06FPY  | 2X-Large  | White       | 48 pr/cs |
| MDTPS4B06FPYH | 2X-Large  | White       | ea       |
| MDTPS4B07FPY  | XXX-Large | White       | 48 pr/cs |
| MDTPS4B07FPYH | XXX-Large | White       | ea       |

### Red – High Risk

| Item No.      | Size               | Tread Color | Pkg.     |
|---------------|--------------------|-------------|----------|
| MDT211218R    | One size fits most | White       | 48 pr/cs |
| MDT211218RH   | One size fits most | White       | ea       |
| MDTPS4B01FPR  | Child              | Gray        | 48 pr/cs |
| MDTPS4B03FPR  | Medium             | Gray        | 48 pr/cs |
| MDTPS4B03FPRH | Medium             | Gray        | ea       |
| MDTPS4B04FPR  | Large              | Gray        | 48 pr/cs |
| MDTPS4B05FPR  | X-Large            | Gray        | 48 pr/cs |
| MDTPS4B06FPR  | XX-Large           | Gray        | 48 pr/cs |
| MDTPS4B07FPR  | XXX-Large          | Gray        | 48 pr/cs |

#### Identify The Risk

- » Yellow – Moderate Risk
- » Red – High Risk



Double-sided, non-slip treads help keep patients on their feet.



**Get the ongoing support you need to drive lasting change.**

Discover **STAND FIRM**, our Falls Management Program on page 4.



# Falls Prevention Gowns

## PerforMAX™ Falls Prevention IV Gowns

- » Quickly identify high-risk patients with this bright yellow gown
- » Waist ties are sewn into side seam and bar tacked for extra strength
- » Side tie eliminates the need for double gowning
- » Durable PerforMAX polyester fabric has superior color retention wash after wash

| Item No.     | Size      | Pkg.    |
|--------------|-----------|---------|
| MDTPG5ITSHRG | One Size  | 6 dz/cs |
| MDTPG7ITSHRG | One Size  | 6 dz/cs |
| MDTOG7ITSHRG | 3X-Large  | 1 dz/cs |
| MDTBG7ITSHRG | 10X-Large | 1 dz/cs |

### Identify The Risk

- Yellow — Moderate Risk
- Red — High Risk



MDTPG7ITSHRG



Falls Prevention Kits

» Packed with the essentials to identify moderate and high falls risk patients

| Item No.      | Description                        | Pkg.  |
|---------------|------------------------------------|-------|
| ➤ DYK100MFALL | Moderate Falls Risk Prevention Kit | 60/cs |
| ➤ DYK100HFALL | High Falls Risk Prevention Kit     | 60/cs |

Optional Component – Falls Risk Blankets

|                |                                       |       |
|----------------|---------------------------------------|-------|
| ➤ MDTSB8F14FPY | Moderate Falls Risk, 30 x 50", Yellow | 60/cs |
| ➤ MDTSB8F14FPR | High Falls Risk, 30 x 50", Red        | 60/cs |



DYK100MFALL



**Moderate Falls Risk Prevention Kit** includes a 3" x 3" yellow magnet for placement on a patient's doorframe to identify them as a moderate fall risk patient, a pair of extra-large slippers and a yellow fall risk ID band that's easily secured around the wrist.



DYK100HFALL



**High Falls Risk Prevention Kit** includes a 3" x 3" red magnet for placement on a patient's doorframe to identify them as a high fall risk patient, a pair of one size fits all red slippers and a yellow fall risk ID band that's easily secured around the wrist.



## Targeted Interventions

**Protect patients at  
their most vulnerable.**

Help prevent falls injuries with  
products designed to anticipate  
risks and enhance mobility.





### Activity Apron

- » Bright colors and texture variation provide visual and tactile stimulation for patients with decreased cognitive function
- » Includes clear picture frame for photos
- » Secures like an apron around waist
- » Machine washable

| Item No.    | Description    | Pkg. |
|-------------|----------------|------|
| MDTACTIVITY | Activity Apron | ea   |



MDTOVERLAY

### » Activity Overlay

- » Provides motor/sensory stimulation for patients with diminished cognitive function
- » Bright colors and texture provide visual and tactile stimulation
- » Attaches to lap or chair with sticky tape
- » Machine washable

| Item No.   | Description      | Pkg. |
|------------|------------------|------|
| MDTOVERLAY | Activity Overlay | ea   |

MDTFALLMATB



MDTFALLMATG

### Beveled Edge Fall Mats

- » Shock-absorbent mats help protect patients from fall-related injuries
- » Low profile, beveled edges help reduce trip hazard
- » Non-skid material helps prevent mat from sliding on floor

| Item No.    | Description         | Pkg. |
|-------------|---------------------|------|
| MDTFALLMATB | Brown, 23 x 68 x 1" | ea   |
| MDTFALLMATG | Gray, 23 x 68 x 1"  | ea   |

MDTFM24662B1



### Fall Mats

- » Our most economical fall mats help reduce injury due to a fall from bed
- » Absorbent, high-compression foam in bi-fold or non-folding style
- » Durable, easy-to-clean vinyl cover

#### Bi-Fold

| Item No.     | Size         | Pkg. |
|--------------|--------------|------|
| MDTFM24662B1 | 24 x 66 x 2" | ea   |
| MDTFM24663B1 | 24 x 66 x 3" | ea   |
| MDTFM36662B1 | 36 x 66 x 2" | ea   |
| MDTFM48662B1 | 48 x 66 x 2" | ea   |

#### Non-Folding

| Item No.   | Size         | Pkg. |
|------------|--------------|------|
| MDTFM24662 | 24 x 66 x 2" | ea   |
| MDTFM24663 | 24 x 66 x 3" | ea   |
| MDTFM36662 | 36 x 66 x 2" | ea   |
| MDTFM48662 | 48 x 66 x 2" | ea   |



HIPPREMIUMSH

## Men's Premium Hip Protector Briefs

- » Help prevent hip injury with shock-absorbing hip protector undergarments
- » Comfortable and discreet closed-cell foam pads are positioned directly over hip joints for maximum protection
- » Machine washable

| Item No.      | Size     | Waist  | Hips   | Pkg. |
|---------------|----------|--------|--------|------|
| HIPPREMIUMSH  | Small    | 30-38" | 35-37" | ea   |
| HIPPREMIUMMLH | Large    | 40-48" | 37-41" | ea   |
| HIPPREMIUMXLH | X-Large  | 44-54" | 41-45" | ea   |
| HIPPREMXXLH   | XX-Large | 48-60" | 45-50" | ea   |



PSY6031M

## Hipsters Men's Standard Briefs

- » Body-contouring, removable hip pads within the brief help provide protection throughout the critical hip fracture area
- » Hip pads are impact-absorbing, perforated EVA foam that allows skin to breathe, while reducing perspiration and resulting skin breakdown
- » Slip on easily over undergarments or wear as underwear
- » Virtually undetectable under clothing
- » Comfortable and washable polyester and cotton blend

| Item No. | Size   | Waist  | Hips    | Pkg. |
|----------|--------|--------|---------|------|
| PSY6031S | Small  | 28-30" | 35-37"  | ea   |
| PSY6031M | Medium | 30-34" | 37-41"  | ea   |
| PSY6031L | Large  | 34-38" | 41"-45" | ea   |



PSY6019M

## Hipsters Men's EZ-On Briefs

- » High impact absorbing removable pads help protect hip bones against injury from falls
- » Feature a crotchless design that allows patients to wear their own undergarments
- » Mesh fabric is water permeable for use during bathing
- » Virtually undetectable under clothing
- » Comfortable and washable polyester and cotton blend

| Item No. | Size   | Waist  | Hips   | Pkg. |
|----------|--------|--------|--------|------|
| PSY6019S | Small  | 28-30" | 35-37" | ea   |
| PSY6019M | Medium | 30-34" | 37-41" | ea   |
| PSY6019L | Large  | 34-38" | 41-45" | ea   |



### Hipsters Men's Incontinent Briefs

- » High impact-absorbing pads help protect hip bones against injury from falls
- » Discreet, low-profile standard foam pads are virtually undetectable under clothing to help improve compliance versus hard-shell products
- » Comfortable and washable polyester and cotton blend

| Item No.   | Size   | Waist  | Hips   | Pkg. |
|------------|--------|--------|--------|------|
| PSY6017SEA | Small  | 28–30" | 35–37" | ea   |
| PSY6017MEA | Medium | 30–34" | 37–41" | ea   |
| PSY6017LEA | Large  | 34–38" | 41–45" | ea   |



PSY6017SEA

### Hipsters Unisex Sweatpants

- » Removable hip pads provide protection throughout the critical fracture area to help reduce the risk of injury from a fall
- » Long pants provide additional comfort and warmth
- » Neutral navy color for women and men
- » Comfortable and washable polyester and cotton blend

| Item No.  | Size   | Waist  | Hips   | Pkg. |
|-----------|--------|--------|--------|------|
| PSY6009NM | Medium | 28–30" | 35–37" | ea   |
| PSY6009NL | Large  | 30–34" | 37–41" | ea   |



PSY6009NM

### Unisex Protective Caps

- » Absorbing, built-in soft foam pads help reduce the force of impact to the head from a fall
- » Discreet, low profile baseball style cap
- » Fleece interior offers comfort and reduced heat buildup
- » Comfortable and washable lightweight cotton blend

| Item No. | Size   | Circumference | Pkg. |
|----------|--------|---------------|------|
| PSY6040M | Medium | 22.5"         | ea   |
| PSY6040L | Large  | 24"           | ea   |



PSY6040M



MDS86960KDMB

## Knockdown Transfer Bench with Microban

- » Microban® antimicrobial product protection is built-in to this transfer bench, preventing the uncontrolled growth of bacteria, mold and mildew on the bath safety product and providing 24/7 cleanliness protection
- » Suction cup foot piece secures bench to tub floor
- » Lightweight, easy to assemble
- » Side arm provides additional stability

| Item No.      | Description    | Weight Capacity | Seat Height Range | Seat Depth | Width Between Arms | Overall Width | Pkg. |
|---------------|----------------|-----------------|-------------------|------------|--------------------|---------------|------|
| MDS86960KDMB  | Transfer Bench | 350 lbs         | 19.5–23.5"        | 16"        | N/A                | 26"           | 2/cs |
| MDS86960KDMBH | Transfer Bench | 350 lbs         | 19.5–23.5"        | 16"        | N/A                | 26"           | ea   |



MDS89745KDMB

## Knockdown Bath Bench with Microban

- » Microban antimicrobial product protection is built-in to make your bath safety product last longer and look cleaner
- » Helps prevent the uncontrolled growth of bacteria, mold and mildew on the bath safety product and provides 24/7 cleanliness protection
- » Bath benches help patients who have difficulty sitting in a standard bath tub or standing in the shower
- » Benches with a back give extra support and comfort; Benches without a back provide easy access and maneuverability when bathing
- » Non-slip, rubber feet resist slipping on the bathtub or shower floor

| Item No.      | Description | Weight Capacity | Seat Height Range | Seat Depth | Width Between Arms | Overall Width | Pkg. |
|---------------|-------------|-----------------|-------------------|------------|--------------------|---------------|------|
| MDS89745KDMB  | w/ Back     | 300 lbs         | 13.75–20.75"      | 12.25"     | N/A                | 17"           | 2/cs |
| MDS89745KDMBH | w/ Back     | 300 lbs         | 13.75–20.75"      | 12.25"     | N/A                | 17"           | ea   |
| MDS89740KDMB  | No Back     | 300 lbs         | 13.75–20.75"      | 12"        | N/A                | 17"           | 2/cs |
| MDS89740KDMBH | No Back     | 300 lbs         | 13.75–20.75"      | 12"        | N/A                | 17"           | ea   |



MDS89664KDMB

## Commode with Microban

- » Microban antimicrobial product protection is built-in to make your bath safety product last longer and look cleaner
- » Tool-free assembly
- » Powder-coated for easy cleaning
- » Splash guard included for over-the-toilet use
- » Features a convenient toilet paper holder
- » Hospital quality construction

| Item No.     | Description | Weight Capacity | Seat Height Range | Seat Depth | Width Between Arms | Overall Width | Pkg. |
|--------------|-------------|-----------------|-------------------|------------|--------------------|---------------|------|
| MDS89664KDMB | Steel       | 350 lbs         | 15.5–21.25"       | 19"        | N/A                | 25.5"         | 1/cs |



### Folding 3-In-1 Commodes

- » Folding frame design easily opens and folds, and will fold flat for convenient storage and transportation
- » Clip-on seat can be removed for easy cleaning
- » Includes pail, lid and splashguard
- » Elongated seat commode has a deeper seat to allow more room for cleaning

| Item No.   | Description           | Weight Capacity | Seat Height Range | Seat Depth | Width Between Arms | Overall Width | Pkg. |
|------------|-----------------------|-----------------|-------------------|------------|--------------------|---------------|------|
| G30213-4F  | Basic, Steel          | 350 lbs         | 16.5–22.5"        | 16"        | 19"                | 21.25"        | 4/cs |
| G30210-4F  | Premium, Aluminum     | 250 lbs         | 17.5–21.5"        | 16"        | 19"                | 21.25"        | 4/cs |
| MDS89664EL | Elongated Seat, Steel | 400 lbs         | 17–23"            | 16"        | 19"                | 22"           | 4/cs |



G30213-4F

### Non-Folding Commodes

- » Use in the bathroom as an elevated toilet seat or as toilet safety rails when combined with most standard toilets
- » Includes pail, lid and splashguard
- » Clip-on seat can be removed easily for cleaning
- » Gray powder-coated steel commode is rust-resistant and durable
- » Bariatric commodes are extra-wide for larger patients

| Item No.   | Description                 | Weight Capacity | Seat Height Range | Seat Depth | Width Between Arms | Overall Width | Pkg. |
|------------|-----------------------------|-----------------|-------------------|------------|--------------------|---------------|------|
| MDS89664   | Basic, Steel                | 350 lbs         | 15–21"            | 15"        | 19"                | 21.25"        | 4/cs |
| G30216B    | Bariatric, Steel            | 500 lbs         | 16.5–23.5"        | 18.5"      | 23.5"              | 27"           | 2/cs |
| MDS89664XW | Extra-Wide Bariatric, Steel | 650 lbs         | 18–23.5"          | 15"        | 24"                | 28"           | ea   |



MDS89664

### Drop-Arm Commodes

- » Arms swing out of the way for safe lateral transfers on and off the commode
- » Durable, rust-resistant, powder coated reinforced frame
- » Includes pail, lids and splashguard

| Item No.   | Description                | Weight Capacity | Seat Height Range | Seat Depth | Width Between Arms | Overall Width | Pkg. |
|------------|----------------------------|-----------------|-------------------|------------|--------------------|---------------|------|
| MDS89668   | Basic, Steel               | 250 lbs         | 17–21"            | 15"        | 17.5"              | 25.5"         | ea   |
| G98202     | Standard Non-Padded, Steel | 300 lbs         | 18.25–23.25"      | 18–20"     | 18.5"              | 22"           | ea   |
| G98204     | Premium Padded, Steel      | 350 lbs         | 20–25"            | 18–20"     | 18.5"              | 22"           | ea   |
| MDS89668XW | Bariatric, Steel           | 850 lbs         | 20–25"            | 19"        | 28"                | 29"           | ea   |



G98202



## Grab Bars


- » Help make bathrooms and stairways safer
- » Chrome bars have a non-slip knurled surface
- » Hanging hardware and instructions included

### Chrome

| Item No.    | Length | Pkg. |
|-------------|--------|------|
| MDS86012CHR | 12"    | 3/cs |
| MDS86016CHR | 16"    | 3/cs |
| MDS86018CHR | 18"    | 3/cs |
| MDS86024CHR | 24"    | 3/cs |
| MDS86032CHR | 32"    | 3/cs |

### White


| Item No. | Length | Pkg. |
|----------|--------|------|
| G87520M  | 12"    | 3/cs |
| G87530M  | 16"    | 3/cs |
| G87540M  | 18"    | 3/cs |
| G87550M  | 24"    | 3/cs |
| G87560M  | 32"    | 3/cs |
| G87570M  | 36"    | 3/cs |



MDTVSQUARE



MDTVBEVEL



MDTVCBEVEL

## Veriflex Shower Fall Mats

- » Build your own custom fall mat with tiles that snap together
- » Flexible, slip resistant surface and shock absorbing design helps impact and injury from a fall
- » Water flows through and around the mat to drain
- » Infused with nanosilver to limit bacterial and fungi growth
- » Beveled edge allows easy access for shower chairs and lifts

| Item No.   | Description            | Length         | Pkg. |
|------------|------------------------|----------------|------|
| MDTVSQUARE | Square Tile            | 11.75 x 11.75" | 4/pk |
| MDTVBEVEL  | Beveled Edge, Square   | 11.75 x 11.75" | 4/pk |
| MDTVCBEVEL | Corner Beveled Edge    | 11.75 x 11.75" | 4/pk |
| MDTVKIT    | Shower Mat Starter Kit | 11.75 x 11.75" | ea   |

Versatile modular design allows for a custom fit in most showers







G30270A

### Locking Raised Toilet Seats

- » Front clamping mechanism locks the seat onto the toilet, so it won't come loose or shift during use or transfer
- » Wide, contoured surface raises seat level by 5"
- » Fits most toilet bowls from 11" to 14"L
- » Optional padded armrests assist patient on and off the toilet

| Item No. | Armrests | Weight Capacity | Seat Elevation | Overall W x L | Opening W x L | Pkg. |
|----------|----------|-----------------|----------------|---------------|---------------|------|
| G30260   | No       | 350 lbs         | 5"             | 20.5 x 18"    | 8 x 9.5"      | 3/cs |
| G30270A  | Yes      | 350 lbs         | 5"             | 20.5 x 18"    | 8 x 9.5"      | 3/cs |



G30300

### Toilet Safety Rail

- » Assist patients with joint pain or at risk of falling
- » Handles adjust and rotate back for a wide range of positions
- » Foam armrests provide a comfortable, sure grip
- » Easy-to-clean aluminum frame mounts securely onto bowl
- » Adjustable bracket range of 18 to 24 inches

| Item No. | Description | Weight Capacity | Height Range | Width Range | Overall Depth Arms | Overall Depth Legs | Pkg. |
|----------|-------------|-----------------|--------------|-------------|--------------------|--------------------|------|
| G30300   | Aluminum    | 300 lbs         | 26 x 31"     | 18 x 24"    | 16.75"             | 12"                | 2/cs |



G30250

### Economy Raised Toilet Seat

- » Makes using the toilet easier for those with trouble bending
- » Simple design made of polypropylene resin is easy to clean
- » Non-slip pads and lip inside toilet offer additional stability

| Item No. | Weight Capacity | Overall W x L | Opening W x L | Pkg. |
|----------|-----------------|---------------|---------------|------|
| G30250   | 250 lbs         | 14.5 x 15.5"  | 8 x 9.5"      | 3/cs |

Compatible with Alarms sold on page 25.



MDT9130TL1

### Toilet Sensor Pads

- » Compact toilet sensor pad allows caregivers to monitor patients while maintaining their dignity and privacy
- » Position under the toilet seat
- » 30-day warranty
- » Alarm sold separately

| Item No.   | Pad Size   | Pkg. |
|------------|------------|------|
| MDT9130TL1 | 8 x 6 x 4" | ea   |
| MDT9130TL5 | 8 x 6 x 4" | 5/cs |



DYND80235DH

## Urinals

- » Patented splashproof lid and side storage help prevent accidental tipping and spilling for easier bedside care
- » Graduated in ounces and milliliters
- » 32-oz/1000ml clear plastic container

| Item No.     | Description                            | Pkg.  |
|--------------|--|-------|
| DYND80235DH  | Glow-in-the-dark Lid                   | ea    |
| DYND80234    | Mini Pack, Individually Shrink-Wrapped | 12/cs |
| DYND80235S   | Supreme Male Urinal                    | 48/cs |
| DYND80235SD  | Glow-in-the-dark Lid                   | 48/cs |
| DYND80235SNL | Supreme, No lid                        | 48/cs |
| PSY8270      | Urinal Bottle Holder Accessory         | ea    |



EKJ149B

## Spil-Pruf Male Urinals

- » Discreet and easy-to-use urinal helps increase patient confidence due to its spill and leak proof, odor-reducing features
- » Offer peace of mind to bed and wheelchair bound patients
- » Lightweight and durable polypropylene urinal holds up to 500cc of liquid

| Item No. | Description | Pkg.  |
|----------|-------------|-------|
| EKJ149B  | Black       | 12/cs |
| EKJ149CL | Clear       | 12/cs |



RCH33301

## UltraFlex® Silicone Male External Catheters

- » Built-in adhesive band offers ease of use and patient comfort
- » Avoid allergic reactions to latex proteins with 100% silicone
- » Non-sterile

| Item No. | Description       | Pkg.  |
|----------|-------------------|-------|
| RCH33301 | Small, 25mm       | 30/bx |
| RCH33302 | Medium, 29mm      | 30/bx |
| RCH33303 | Large, 32mm       | 30/bx |
| RCH33304 | Large, 36mm       | 30/bx |
| RCH33305 | Extra-Large, 41mm | 30/bx |

### Supra Low Air DPS Mattress Replacement System

- » Dual alternating pressure and low air-loss therapy in one mattress
- » Functions include alternating pressure, low air-loss, static float, auto firm and seat inflation
- » Low-pressure alarm with audio and visual indicators
- » 400-lb. weight capacity
- » 36 x 80 x 8"

| Item No.      | Description  | Pkg. |
|---------------|--|------|
| MDT24SUPRADPS | Supra DPS Mattress Replacement System                      | ea   |
| MDT24DPSRE    | Supra DPS Mattress Replacement System with 4" Raised Edges | ea   |



Also available with 4" raised air bolsters to help reduce the risk of patient falls and entrapment.



MDT24SUPRADPS

### Supra APL Mattress Replacement System

- » A versatile mattress replacement system suitable for any treatment protocol
- » Functions include alternating pressure, low air-loss, static float, easy-to-use weight dial for customized therapy
- » Cell-on-cell design prevents bottoming out during a power outage
- » Easy-to-clean, waterproof, quilted nylon cover reduces shear and friction
- » 400-lb. weight capacity, 36 x 80 x 8"

| Item No.      | Description                           | Pkg. |
|---------------|---------------------------------------|------|
| MDT24SUPRAAPL | Supra APL Mattress Replacement System | ea   |



MDT24SUPRAAPL



MDTE12C3576P

Pictured with optional pump.



Heel section features gently sloping air tubes which significantly reduce interface pressures in this sensitive area.

By adding the optional pump, EqualizeAire mattresses can be converted to an alternating pressure system featuring a 5-minute cycle time. Order systems with optional pump by adding a "P" at the end of the item number or MDTEAPUMP for the pump only.

## EqualizeAire 12000 Self-Adjusting Therapy Mattresses

- » Features twelve 5" air cells that dynamically adjust to patient's profile
- » Six anatomically designed pressure relief zones
- » Easy-to-use weight dial for customized therapy
- » Easy-to-clean, waterproof, quilted nylon cover reduces shear and friction
- » 500-lb. weight capacity

### With Pump

| Item No.     | Size         | Pkg. |
|--------------|--------------|------|
| MDTE12C3576P | 35 x 76 x 7" | ea   |
| MDTE12C3580P | 35 x 80 x 7" | ea   |
| MDTE12C3584P | 35 x 84 x 7" | ea   |

### Without Pump

| Item No.    | Size         | Pkg. |
|-------------|--------------|------|
| MDTE12C3576 | 35 x 76 x 7" | ea   |
| MDTE12C3580 | 35 x 80 x 7" | ea   |
| MDTE12C3584 | 35 x 84 x 7" | ea   |

## EqualizeAire 9000 Self-Adjusting Therapy Mattresses

- » Self-adjusting convertible, pressure-redistribution mattress features nine 5" air cells that dynamically adjust to patient's profile
- » Four anatomically designed pressure relief zones
- » Cell-on-cell design prevents bottoming out during a power outage
- » Easy-to-clean, waterproof, quilted nylon cover reduces shear and friction
- » 500-lb. weight capacity

### With Pump

| Item No.    | Size         | Pkg. |
|-------------|--------------|------|
| MDTE9C3576P | 35 x 76 x 7" | ea   |
| MDTE9C3580P | 35 x 80 x 7" | ea   |
| MDTE9C3584P | 35 x 84 x 7" | ea   |

### Without Pump

| Item No.   | Size         | Pkg. |
|------------|--------------|------|
| MDTE9C3576 | 35 x 76 x 7" | ea   |
| MDTE9C3580 | 35 x 80 x 7" | ea   |
| MDTE9C3584 | 35 x 84 x 7" | ea   |

Pictured with optional pump.


MDTE9C3576

Also available in non-powered four cell design. Contact your Medline sales representative for more information.

### TheraTech 3500 Memory Foam Mattresses

- » Tri-layered foam design features anatomically correct top layer for envelopment, sculpted middle layer for added pressure redistribution, and supportive base layer for comfort
- » Three anatomically designed pressure redistribution zones work in conjunction with the tri-layered foam design to mirror the body anatomically
- » Sloped heel section with hi-resiliency foam helps to redistribute heel pressure
- » Independent load bearing cells cradle bony prominences to relieve pressure and shear
- » 425-lb. weight capacity

| Item No. | Size           | Pkg. |
|----------|----------------|------|
| MSCT3576 | 36 x 76 x 6.5" | ea   |
| MSCT3580 | 36 x 80 x 6.5" | ea   |
| MSCT3584 | 36 x 84 x 6.5" | ea   |





TARGETED INTERVENTIONS

### TheraTech 5500 Memory Foam Mattresses

- » Tri-layered foam design features anatomically correct top layer for envelopment, sculpted middle layer for added pressure redistribution, and supportive base layer for comfort
- » Five anatomically designed pressure redistribution zones work in conjunction with the tri-layered foam design to mirror the body anatomically
- » Sloped heel section with a top layer of visco elastic memory foam helps redistribute heel pressure
- » Independent load bearing cells cradle bony prominences to relieve pressure and shear
- » 425-lb. weight capacity

| Item No. | Size           | Pkg. |
|----------|----------------|------|
| MSCT5576 | 36 x 76 x 6.5" | ea   |
| MSCT5580 | 36 x 80 x 6.5" | ea   |
| MSCT5584 | 36 x 84 x 6.5" | ea   |





MSCT7576

## TheraTech 7500 Memory Foam Mattresses

- » Tri-layered foam design features anatomically correct top layer for envelopment, sculpted middle layer for added pressure redistribution, and supportive base layer for comfort
- » Five anatomically designed pressure redistribution zones work in conjunction with the tri-layered foam design to mirror the body anatomically
- » Top layer features gel-infused visco memory foam for increased comfort and completely conforms to the contours of the body
- » Sloped heel section with two layers of visco elastic memory foam helps redistribute heel pressure
- » 450-lb. weight capacity

| Item No. | Description        | Pkg. |
|----------|--------------------|------|
| MSCT7576 | 36 x 76 x 7"       | ea   |
| MSCT7580 | 36 x 80 x 7"       | ea   |
| MSCT7584 | 36 x 84 x 7"       | ea   |
| MSCTVCR  | Versacare Mattress | ea   |
| MSCTCTR  | Totalcare Mattress | ea   |



PSY5725

## Seizure Side Rail Pads

- » Form fitted pads fit snugly over existing bedside rails
- » 1" layer of medium density foam helps safeguard against patient injury
- » Sheer, lightweight nylon outer shell helps prevent patient "rub" burns and allows unobstructed movement of the side rails to a lowered position
- » Open windows in upper section allow for easy access to controls on most bed models
- » Includes 2 upper and 2 lower pads

| Item No. | Description                                       | Pkg. |
|----------|---|------|
| PSY5742  | Fits CareAssist and Advanta 2 Beds                | ea   |
| PSY5746  | Fits TotalCare Med/Surg Beds                      | ea   |
| PSY5755  | Fits Secure I and II Beds with Rounded Side Rails | ea   |
| PSY5763  | Fits GO Bed II with Rectangular Side Rails        | ea   |
| PSY5725  | Fits Advanta Bed, Waterproof                      | 4/cs |



PSY6311L

## Torso Wedge

- » Help keep bed-ridden patients repositioned and comfortable
- » Ideal for use with alternating pressure mattresses to help maintain support without slipping like foam products
- » Easier to shape and position than regular foam pillows
- » Lightweight with polystyrene bead filling for reduced heat build up
- » Washable mesh cover

| Item No. | Description            | L x W x H     | Pkg. |
|----------|------------------------|---------------|------|
| PSY6311L | 30/60 degree angle, XL | 27 x 22 x 18" | ea   |



### Bed Assist Bar

- » Provides help getting in and out of bed
- » Features a soft grip for comfort, pocket for handy storage and straps for added safety
- » Built to last with quick and easy tool-free assembly
- » Includes a nylon bag for storage

| Item No.  | Description    | Pkg. |
|-----------|----------------|------|
| MDS6800BA | Bed Assist Bar | 3/cs |



MDS6800BA

### Side Rail Protector

- » Help reduce entrapment by closing up gaps within and between the bed side rails
- » Wipeable vinyl composition
- » Pads secure easily with multiple hook and loop attachments

| Item No.  | Description         | Pkg. |
|-----------|---------------------|------|
| MDT821570 | Side Rail Protector | ea   |

MDT821570



### Advantage Mattress Bolster

- » Transform any mattress to a raised edge mattress
- » 8" bolsters help prevent patient falls and lower the risk of entrapment
- » Quick installation and removal
- » Easy-to-clean waterproof nylon cover reduces shear and friction

| Item No.    | Description      | Pkg. |
|-------------|------------------|------|
| MDTMATBOLST | Mattress Bolster | ea   |

MDTMATBOLST



**Get the ongoing support you need to drive lasting change.**

Discover **STAND FIRM**, our Falls Management Program on page 4.

# Safety Essentials

**Handle high-risk patients with care.**

Support patients safely and securely with high-quality products at a great value.





MDT9100

### Quick Alert Nurse Call Alarms

- » For use with bed, chair, seatbelt, toilet and floor sensor mats
- » Distinguish between patients with 8 different tones and tunes
- » Includes 3 volume levels with low battery indicator, yellow protective cover, 2 AAA batteries
- » Basic alarm contains all the same features without delay option, tones, and tunes

| Item No. | Description                   | Pkg. |
|----------|-------------------------------|------|
| MDT9100  | Universal Alarm, Multi-tone** | ea   |
| MDT9105  | Universal Alarm**             | ea   |



### Recordable Alarm Monitor

- » Customizable with your desired message or alarm tone/sound
- » When patient sits up in bed or stands up from a chair, the monitor will alert you and get the attention of patient with your message
- » Protective boot for maximum protection
- » Auto reset functionality so alarm is always ready to use

| Item No.  | Description                | Pkg. |
|-----------|----------------------------|------|
| MDT800400 | Recordable Alarm Monitor** | ea   |



MDT8450Y

### Tamper-Resistant Personal Safety Alarm

- » Auto-sensing technology monitors patient movement and sounds an alarm if patient gets up from a bed or wheelchair
- » Prevent tampering and patient interference by pressing two buttons simultaneously to disengage the alarm
- » Low battery indicator light and volume control
- » Extra-durable wheelchair attachment clip included

| Item No.  | Description                                     | Pkg. |
|-----------|---|------|
| MDT8450Y  | Tamper-Resistant Alarm Monitor**                | ea   |
| MDT84BRKT | Optional Patient Alarm Bracket for MDT84 Series | ea   |



MDT8201

### Economy Alarm Monitor

- » Alarm alerts staff when patient gets out of bed or chair
- » Affordable monitor for use with most bed or chair sensor pads and floor mats
- » Alarm and low battery lights communicate functionality, so you're always powered up
- » Auto-sensing technology monitors patient movement and triggers the alarm

| Item No. | Description                    | Pkg. |
|----------|--------------------------------|------|
| MDT8201  | Universal Basic Alarm Monitor  | ea   |
| MDT8501  | Basic Economy Alarm, Wide Port | ea   |

\*All alarms are universally compatible with pads listed on page 26, excluding MDT85-series alarms.. \*\* Contact your Medline rep for cord options.





MDT5000Z

## Advantage® Magnetic Patient Alarm

- » Control patient wandering with an alarm that attaches to a bed, chair or patient clothing
- » Magnetic pull switch releases from any angle, helping to eliminate false alarms and clothing tears
- » Recessed on/off switch limits patient interference
- » Battery tester button ensures that alarm is always powered up
- » 2 AAA batteries included

| Item No.   | Description                                     | Pkg. |
|------------|---|------|
| MDT5000    | Alarm, Advantage, Magnetic Tether               | 5/cs |
| MDT5000Z   | Alarm, Advantage, Magnetic Tether               | ea   |
| MDT5000XCD | Advantage Magnetic Alarm Replacement Pull Cords | 3/bg |



CRE101968

## Crest Call Cords

- » Call cords specific for your custom call light system
- » Universal options in 7 ft and 12 ft.
- » Specialty Cords available upon request - please inquire with your Medline Representative

| Item No.    | Description   | Pkg.  |
|-------------|---|-------|
| CRE101968   | Nurse Call Cord, 7', 1/4" Both Ends                       | ea    |
| CRE10196812 | Nurse Call Cord, 12', 1/4" Both Ends                      | ea    |
| CRE113708   | Pendant Only, Red, Clean Cord                             | 15/pk |
| CRE113743   | Cord, Pad Call, 1/4", Plug, 10', Non-locking              | ea    |
| CRE114192   | Cord, Clean Cord, White, 100'                             | ea    |
| CRE114775   | Cord, Pull, 6', Glow, Dark White/Red                      | ea    |
| CRE7700DB12 | Cord, Call, 12', Beige, 1/4" Both Ends, Phone Plug        | ea    |
| CRE7700DB12 | Cord, Call, 12' Double, Beige, 1/4" Both Ends, Phone Plug | ea    |
| CRE7700DW10 | Cord, Call, Double, 10', White                            | ea    |
| CREB112C    | Call Station Cord, Bedside, Cornell, Double               | ea    |
| CREE1143C   | Station, Emergency Pull Cord                              | ea    |



MDT84NURSE

## Nurse Call Cable for MDT84 Series Alarm

- » Connects the patient alarms in the MDT84 series
- » 8-ft. long cable

| Item No.    | Description                           | Pkg. |
|-------------|---------------------------------------|------|
| MDT84NURSE  | Nurse Call Cable, Alarm, F/MDT84      | ea   |
| MDT84HRADPT | Y-Adapter, MDT8450 to Hil-ROM, 20 Pin | ea   |
| MDT84STRAP  | Chair Strap Alarm, F/MDT84            | ea   |



MDT91YELBT

## Miscellaneous Alarm Accessories

| Item No.   | Description                                 | Pkg. |
|------------|---|------|
| MDT91NURSE | Quick Alert Nurse Call Cord, 1/4" Both Ends | ea   |
| MDT9100AC  | Quick-Alert Series Alarm Adapter, AC        | ea   |
| MDT91YELBT | Yellow Cover, Quick-Alert Series            | ea   |
| MDT81STRAP | Nylon Strap with Quick-Release Buckle       | ea   |

\* Pads are only compatible with MDT8501 alarm series



MDTRUC1

## Reusable Bed/Chair Sensor Pads

- » When attached to an alarm, the sensor pad alerts staff that the patient has taken weight off of the pad and may be at risk of falling
- » Foldable, reusable and easy-to-store
- » High-grade, water-resistant vinyl material provides protection against incontinence

### Chair

| Item No. | Description     | H x W    | Pkg. |
|----------|-----------------|----------|------|
| MDTRUC1  | 30-Day Warranty | 12 x 15" | ea   |

### Bed

| Item No. | Description     | H x W    | Pkg. |
|----------|-----------------|----------|------|
| MDTRUB1  | 30-Day Warranty | 15 x 30" | ea   |



MDT821180C6

## Economy/Single Patient Use Sensor Pads

- » When pad is connected to a monitor alarm, it will trigger the alarm when weight is removed from the pad
- » Durable, high-grade vinyl pad is easy to store
- » For single patient use on a chair or bed (underneath the fitted sheet)
- » MDT82 series includes a full-size, non-slip grip sheet for maximum stabilization

### Chair

| Item No.    | Description      | H x W   | Pkg. |
|-------------|------------------|---------|------|
| MDT83130C6  | 30-Day Warranty  | 7 x 15" | 6/cs |
| MDT82145C6  | 45-Day Warranty  | 7 x 15" | 6/cs |
| MDT831180C6 | 180-Day Warranty | 7 x 15" | 6/cs |
| MDT821180C6 | 180-Day Warranty | 7 x 15" | 6/cs |
| MDT831365C6 | 1-Year Warranty  | 7 x 15" | 6/cs |
| MDT821365C6 | 1-Year Warranty  | 7 x 15" | 6/cs |

### Bed

| Item No.    | Description      | H x W    | Pkg. |
|-------------|------------------|----------|------|
| MDT83130B6  | 30-Day Warranty  | 10 x 30" | 6/cs |
| MDT82145B6  | 45-Day Warranty  | 10 x 30" | 6/cs |
| MDT831180B6 | 180-Day Warranty | 10 x 30" | 6/cs |
| MDT821180B6 | 180-Day Warranty | 10 x 30" | 6/cs |
| MDT831365B6 | 1-Year Warranty  | 10 x 30" | 6/cs |
| MDT821365B6 | 1-Year Warranty  | 10 x 30" | 6/cs |



MDT9190C5

## Water-Resistant Premium Sensor Pads

- » Weight-sensitive pad activates alarm if patient leaves a chair or bed
- » Water-resistant, high-grade vinyl
- » Foldable and easy to store
- » Use on a chair or bed (underneath the fitted sheet)

### Chair

| Item No.    | Description      | H x W    | Pkg. |
|-------------|------------------|----------|------|
| MDT85130C6* | 30-Day Warranty  | 13 x 13" | 6/cs |
| MDT9130C5H  | 30-Day Warranty  | 12 x 15" | ea   |
| MDT9130C5   | 30-Day Warranty  | 12 x 15" | 5/cs |
| MDT9190C5   | 90-Day Warranty  | 12 x 15" | 5/cs |
| MDT91180C5  | 180-Day Warranty | 12 x 15" | 5/cs |

### Bed

| Item No.    | Description      | H x W    | Pkg. |
|-------------|------------------|----------|------|
| MDT85130B6* | 30-Day Warranty  | 13 x 30" | 6/cs |
| MDT9130B5   | 30-Day Warranty  | 15 x 30" | 5/cs |
| MDT9190B5   | 90-Day Warranty  | 15 x 30" | 5/cs |
| MDT91180B5  | 180-Day Warranty | 15 x 30" | 5/cs |

# Sensible Care Proactive Patient Monitoring System

## Reduce bed falls with a high-tech solution

The only clinically proven technology available that has reduced bed falls and alarm fatigue.



Enhance patient safety and security with Sensible Care, a state-of-the art monitoring system that detects early patient movement patterns and alerts staff prior to a potential bed exit.

| Item No.     | Description     | Pkg. |
|--------------|-----------------|------|
| MDT9590B1    | Sensor Pad      | ea   |
| MDT9500      | Alarm Monitor   | ea   |
| MDT95DASH    | Dashboard       | ea   |
| MDT95SERVER  | Computer Server | ea   |
| MDT9590SHEET | Fitted Sheet    | ea   |

Sensor Pad  
MDT9590B1





## A. Sensor Pad

Lightweight, thin, reusable pad is inserted into the fitted sheet to monitor patient movement from the shoulders to the thighs, detect potential bed exits, and reduce bed falls. It also issues repositioning alerts.

## B. Alarm Monitor

When there's an unscheduled bed exit attempt, this bedside device flashes red, issues an audible message to remind the patient to wait for caregiver, and alerts staff.

## C-D. Dashboard and Computer Server

Placed at the nursing station, the nursing dashboard receives real time data from the alarm monitor and can simultaneously monitor up to 35 beds in real time. The computer server collects data from sensors and devices to generate unit-specific reports.

## E. Customizable Fitted Sheet

Poly sheet accommodates the sensor pad and is customized to your bed size.

## Free Mobile Device App

Nurses receive potential bed exit alerts and notifications to reposition patients via the free Sensible Care App available for Apple iOS® and Android™ devices.

C D





MDT8800

### Smart Caregiver® Wireless Pressure-Sensing Monitoring Systems

- » Wireless monitor helps reduce tripping hazards and eliminates broken or tangled pad cords
- » Cordless pads communicate with the monitor using wireless technology with both an audible and visual alert
- » Mount the monitor 15-20 feet away from the sensor pads to eliminate bedside alarm noise in the resident's room

#### Monitors

| Item No.  | Description                             | Pkg. |
|-----------|---|------|
| MDT8800   | Wireless Monitor (Up to 2 pads)         | ea   |
| MDT433CMU | Central Monitoring Unit (Up to 40 pads) | ea   |

#### Wireless Alarm Sensor Pads and Accessories

| Item No.    | Size                                    | Pkg. |
|-------------|---|------|
| MDT8830C1   | Cordless Chair Alarm Pad for MDT8800    | ea   |
| MDT8830B1   | Cordless Bed Alarm Pad for MDT8800      | ea   |
| MDT88365B1  | Bed Sensor Mat, 1-year warranty         | ea   |
| MDT88365C1  | Chair Sensor Mat, 1-year warranty       | ea   |
| SMRTIR1     | Infrared Patient Bed Alarm for MDT8800  | ea   |
| MDT8836XWB1 | Bed Sensor Mat, 1-year warranty, X-Wide | ea   |
| MDT8800FM   | Cordless Floor Mat                      | ea   |
| MDT433NC    | Nurse Call Button for MDT433CMU         | ea   |
| MDT433RB    | Reset Call Button for MDT433CMU         | ea   |
| MDT433PGD   | Pager with LCD for MDT433CMU            | ea   |
| MDT433PAGER | Pager without LCD for MDT433CMU         | ea   |



MDTIRM1

### Infrared Patient Bed Alarm

- » Alarm affixes to the headboard or bed frame, so patient can move freely within bed
- » If patient attempts to leave bed, the infrared beam detects motion and sends out your choice of alarm or chime tone
- » Alarm head swivels for precise beam alignment
- » Attach using the bracket (optional) or hook and loop strips (included)
- » Requires 9V battery (not included)

#### Alarm

| Item No. | Description    | Pkg. |
|----------|----------------|------|
| MDTIRM1  | Infrared Alarm | ea   |

#### Metal Mounting Bracket For Infrared Alarm

| Item No.    | Description            | Pkg. |
|-------------|------------------------|------|
| MDTIRM1BRKT | Metal Mounting Bracket | ea   |



MDT86100

### Anti-Wander Alarm Kits and Accessories

- » A great solution for patients with Alzheimers or dementia who may be prone to wandering
- » Kit includes three easy-to-use components: a door strip monitor, magnetic contact switch, and resident wristband
- » Easy installation with no hard wiring

| Item No.     | Description                  | Pkg. |
|--------------|------------------------------|------|
| MDT86100     | Anti-Wander Kit, One Door    | ea   |
| MDT86150     | Anti-Wander Kit, Two Doors   | ea   |
| MDT86100TC   | Caregiver Transmitter        | ea   |
| MDT86100TR   | Resident Transmitter         | ea   |
| MDT86100DT   | Door Tester                  | ea   |
| MDT86100W    | Reset Button                 | ea   |
| MDT86WRSTBND | Wristband (Only) Replacement | ea   |



MDT821203YS



MDT821203



MDT821207BLUL



MDT828203BWP



MDT828203PWP

### Cotton Gait Belts

- » Belt fits securely around the patient's waist and provides caregivers with a surface to grab when helping patients ambulate or transfer
- » Aid with lifting and transferring patients
- » Additional support helps reduce the risk of caregiver back injuries
- » Machine washable

| Item No.      | Design                              | Length | Buckle  | Pkg. |
|---------------|-------------------------------------|--------|---------|------|
| MDT821203YS   | Solid Yellow Fall Management        | 54"    | Metal   | ea   |
| MDT821203Y    | Solid Yellow Fall Management        | 60"    | Metal   | ea   |
| MDT821203YL   | Solid Yellow Fall Management        | 72"    | Metal   | ea   |
| MDT821203     | White with Blue and Red Pin Stripes | 60"    | Metal   | ea   |
| MDT828203     | White with Blue and Red Pin Stripes | 60"    | Metal   | 6/cs |
| MDT821203L    | White with Blue and Red Pin Stripes | 72"    | Metal   | ea   |
| MDT821203NL   | White with Black Stripes            | 100"   | Metal   | ea   |
| MDT821204     | Pastel Stripes                      | 60"    | Metal   | ea   |
| MDT828204C    | Pastel Stripes                      | 60"    | Metal   | 6/cs |
| MDT828208PASL | Pastel                              | 72"    | Metal   | ea   |
| MDT821203BS   | Solid Black                         | 54"    | Metal   | ea   |
| MDT821203B    | Solid Black                         | 60"    | Metal   | ea   |
| MDT821203BL   | Solid Black                         | 72"    | Metal   | ea   |
| MDT828203NAT  | Natural                             | 60"    | Metal   | ea   |
| MDT82803BLU   | Navy Blue                           | 60"    | Metal   | ea   |
| MDT821208NAT  | Natural                             | 60"    | Plastic | ea   |
| MDT821208PAS  | Pastel Stripes                      | 60"    | Plastic | ea   |
| MDT828209PASL | Pastel                              | 72"    | Plastic | ea   |

### Wipeable Gait Belts

- » Versatile belts provide help with mobilizing and transferring patients
- » Straps are durable and easy to keep clean
- » Save laundry time—wipe clean with standard PDI wipes/Hospital-grade disinfectant
- » Buckle provides easy adjustment

#### Polypropylene—Flexible, webbed material

| Item No.      | Design       | Length | Buckle  | Pkg. |
|---------------|--------------|--------|---------|------|
| MDT821206BLU  | Solid Blue   | 60"    | Plastic | ea   |
| MDT821206YEL  | Solid Yellow | 60"    | Plastic | ea   |
| MDT821206YELL | Solid Yellow | 72"    | Plastic | ea   |
| MDT821206B    | Solid Black  | 60"    | Plastic | ea   |
| MDT821206BL   | Solid Black  | 72"    | Plastic | ea   |
| MDT821206RAIN | Solid Black  | 60"    | Plastic | ea   |
| MDT821207BLU  | Solid Blue   | 60"    | Metal   | ea   |
| MDT821207BLUL | Solid Blue   | 72"    | Metal   | ea   |

#### Polyurethane—Smooth, durable material

|              |             |     |         |    |
|--------------|-------------|-----|---------|----|
| MDT828203PWP | Solid Pink  | 60" | Plastic | ea |
| MDT828203BWP | Solid Black | 60" | Plastic | ea |
| MDT821203WP  | Solid Blue  | 60" | Plastic | ea |
| MDT821203WPL | Solid Blue  | 70" | Plastic | ea |



MDT823256

## Hand Protectors

- » Provide a way for patients to exercise fingers during rehabilitation and to prevent scratching or picking and limit contracture in stroke patients
- » Secure to bed frame or wheelchair with included arm restraint
- » Machine washable polyester and cotton fabric is filled with resilient polyester pile
- » One size, for left or right hands
- » Safety check design has top flap for checking patients' hands
- » Finger-control shoelace style mitt separates fingers and features a breathable mesh backing

### Adult Hand Protectors

| Item No.  | Description    | Straps | Closure                        | Pkg. |
|-----------|----------------|--------|--------------------------------|------|
| MDT823256 | Flexible Palm  | Double | Metal Buckle with Teeth Slider | pr   |
| MDT823266 | Rigid Palm     | Double | Metal Buckle with Teeth Slider | pr   |
| MDT823243 | Safety Check   | Single | Hook-and-Loop Material         | pr   |
| MDT823242 | Finger-Control | Single | Shoelaces                      | pr   |

### Toddler and Youth Mitts

| Item No. | Description | L x W        | Wrist Circumference | Pkg. |
|----------|-------------|--------------|---------------------|------|
| PSY2911M | Toddler     | 3.93 x 1.75" | 5"                  | pr   |
| PSY2911L | Youth       | 6.5 x 1.75"  | 5.5"                | pr   |
| PSY2911S | Pediatric   | 3.54 x 1.75" | 4.75"               | pr   |

MDT823242

MDT823266



### Multi-Purpose Arm and Leg Splints

- » Keep patient from pulling at tubes, drains and lines with an arm splint
- » Easily applied with hook and loop closure
- » Vinyl mesh outer covering has a 100% cotton lining
- » Zippered opening enables beads to be added or removed for patient comfort

| Item No.   | Style               | Limb Length | Limb Circumference | Pkg. |
|------------|---------------------|-------------|--------------------|------|
| MDT828816P | Pediatric           | 7"          | 1.5–4"             | ea   |
| PSY8168XS  | Pediatric 2–7 years | 10"         | 1.5–4"             | ea   |
| MDT828816S | Small               | 11"         | 3–14"              | ea   |
| MDT828816M | Medium              | 13"         | 5–17"              | ea   |
| MDT828816L | Large               | 14"         | 8–21"              | ea   |
| PSY8169    | Adult Leg           | 18"         | 14–27"             | ea   |

MDT823243



MDT828816M



MDT823100

## Quick Release Limb Holders

- » Secure wrists or ankles to bed frame with hook and loop closure
- » Will not overtighten with patient movement
- » Non-detachable double straps offer extra range of motion restriction
- » 7 x 2.5" quilted cuff with double straps for range of motion

| Item No.    | Description                      | Cuff L x W | Pkg  |
|-------------|----------------------------------|------------|------|
| MDT823100   | Light Blue, Quick Release        | 7 x 2.5"   | 1 pr |
| MDT821301Q  | Black, Quick Release             | 7 x 4"     | 1 pr |
| MDT821301QL | Black, Long Strap, Quick Release | 7 x 4"     | 1 pr |



MDT823040

## Quick Control Limb Holders

- » Help keep patients safe who may be a danger to themselves and staff
- » 36" non-detachable straps offer additional range of motion restriction
- » 9 x 3" size features a slide harness strap and synthetic fur cuff for patient comfort
- » Machine washable

| Item No.   | Material                 | Straps                         | Cuff L x W | Pkg  |
|------------|--------------------------|--------------------------------|------------|------|
| MDT823040  | Orthopedic Felt, Cotton  | Metal Buckle and Hook-and-Loop | 10 x 3"    | 1 pr |
| MDT823070* | Synthetic Fur, Polyester | Metal Buckle with Teeth Slider | 9 x 3"     | 1 pr |

\*For a case of 6 pairs, order item MDT829070C.



MDT829030

## Disposable Limb Holder

- » Single-patient use for economy and efficiency
- » Secure wrists or ankles to a bed frame
- » Metal buckle with teeth slider cinch closure
- » 10 x 3" foam lined cuff with double 30" long non-stretch, non-detachable anchor straps

| Item No.  | Straps                         | Cuff L x W | Pkg     |
|-----------|--------------------------------|------------|---------|
| MDT829030 | Metal Buckle with Teeth Slider | 10 x 3"    | 6 pr/cs |



MDT829082C

## Reusable Padded Limb Holder

- » Secure wrists or ankles of at-risk patients to a bed frame
- » Will not tighten with patient movement
- » Slide harness minimizes chafing during wrist rotation
- » 10 x 4" fur lined cuff with 36" long single strap

| Item No.   | Straps                         | Cuff L x W | Pkg     |
|------------|--------------------------------|------------|---------|
| MDT829082C | Metal buckle with teeth slider | 10 x 4"    | 6 pr/cs |





### Foam Limb Holders

- » Prevent patients from pulling at medical devices, interrupting treatments and disrupting staff members
- » Padded foam cuff secures wrists or ankles to a bed frame
- » Will not tighten with patient movement
- » Adjustable 3/8" D-ring closure
- » 13.5 x 3" foam lined cuff with 55" long double straps

#### Adult

| Item No.  | Straps                                 | Cuff L x W | Pkg  |
|-----------|--|------------|------|
| MDT829520 | Double Strap with Metal D-ring closure | 13.5 x 3"  | 1 pr |
| MDT829528 | Single Strap with Metal D-ring closure | 13.5 x 3"  | 1 pr |

#### Infant

| Item No.   | Straps        | Cuff L x W | Pkg      |
|------------|---------------|------------|----------|
| MDT829473  | Hook-and-Loop | 5.5 x 1"   | 24 ea/bx |
| MDT829474P | Hook-and-Loop | 10 x 2"    | 24 ea/bx |



### Personal Safety Limb Holder

- » Cotton personal safety limb holders prevent patients from pulling out an IV line or catheter
- » For use on wrists or ankles, they won't over tighten with movement
- » Synthetic fur lined cotton cuff with adjustable hook-and-loop fastener and straps that tie to bed frame
- » 10 x 4" fur lined cuff with 44" long double straps

| Item No.  | Description           | Cuff L x W | Pkg     |
|-----------|-----------------------|------------|---------|
| MDT829052 | Synthetic Fur, Cotton | 10 x 4"    | 6 pr/cs |



### Disposable Personal Safety Limb Holder

- » Thick, breathable foam patient safety devices feature adjustable hook-and-loop closures to ensure a comfortable fit
- » Limits movement of wrists or ankles
- » Nylon D-ring provides added security
- » Anchor straps will not stretch
- » 8.75 x 3.25" foam lined cuff with 45" long double strap

| Item No.    | Description | Cuff L x W   | Pkg   |
|-------------|-------------|--------------|-------|
| NON243074   | Limb Holder | 8.75 x 3.25" | 48/cs |
| NON243074HH | Limb Holder | 8.75 x 3.25" | ea    |



### Cinch-Type Patient Safety Limb Holder

- » Foam limb holder limits movement of wrists or ankles
- » Hook-and-loop tab adheres anywhere on the foam backing
- » D-ring stress point prevents cuff from tightening when patient pulls on it
- » Adjusts for adult or pediatric use
- » 10 x 3" foam lined cuff with 40" long strap

| Item No.   | Description                    | Cuff L x W | Pkg     |
|------------|--------------------------------|------------|---------|
| MDT829020C | Cinch-Type with D-Ring Closure | 10 x 3"    | 6 pr/cs |



MDT829098W



MDT829090W



MDT829090WQ



PSY2214

## Wrist and Ankle Limb Holders

- » For patients who are a threat to themselves or others in locking and non-locking styles
- » Neoprene cuffs are gentle on skin and adjust to fit most limb sizes
- » Two universal keys included with locking cuffs
- » Hook and loop closure for easy adjustment
- » Machine washable

### Keylock Limb Holders

| Item No.   | Description | Cuff L x W  | Strap Length | Pkg. |
|------------|-------------|-------------|--------------|------|
| MDT829098W | Wrist       | 12 x 2.5"   | 46"          | pr   |
| MDT829099A | Ankle       | 14.5 x 2.5" | 44.5"        | pr   |

### Hook & Loop Ring Limb Holders

| Item No.   | Description | Cuff L x W  | Strap Length | Pkg. |
|------------|-------------|-------------|--------------|------|
| MDT829090W | Wrist       | 12 x 2.5"   | 48.5"        | pr   |
| MDT829091A | Ankle       | 14.5 x 2.5" | 50"L         | pr   |

### Quick-Release Limb Holders

| Item No.    | Description | Cuff L x W | Strap Length | Pkg. |
|-------------|-------------|------------|--------------|------|
| MDT829090WQ | Wrist       | 12 x 2.5"  | 45"          | pr   |
| MDT829091AQ | Ankle       | 12 x 2.5"  | 46.5"        | pr   |

### Stretcher Limb Holders

| Item No.   | Description | Cuff L x W  | Strap Length | Pkg. |
|------------|-------------|-------------|--------------|------|
| MDT8292750 | Wrist       | 12 x 2.5"   | 9L x 2"W     | pr   |
| MDT8292755 | Ankle       | 14.5 x 2.5" | 9L x 2"W     | pr   |

### Synthetic Leather Connected Locking Limb Holders

| Item No. | Description                             | L x W   | Strap Length | Pkg. |
|----------|---|---------|--------------|------|
| PSY2213  | Connecting Synthetic Cuff (Wrist/Ankle) | 38 x 3" | 32"          | pr   |
| PSY2214  | Wrist                                   | 38 x 3" | 32"          | pr   |
| PSY2215  | Keylock Synthetic Cuff (Ankle Only)     | 38 x 3" | 32"          | pr   |

## Repositioning Belts

- » Ideal for patients with minimal upper trunk control who require positioning support and are at risk of falling, or those who require a restraint for medical treatment
- » For use on a bed, stretcher or chair, wheelchair
- » For combative patients, utilize keylock restraints for maximum security
- » Flame-retardant cotton wraparound belt has a secure hook-and-loop closure
- » Roll belts are less restrictive in that they allow patients to sit up or roll from side to side in bed

### Locking Belts with Keylock

(For key replacements use item no.: MDT829074)

| Item No.   | Description                   | Pkg. |
|------------|-------------------------------|------|
| MDT829334C | Safety Belt w/Keylock, Cotton | ea   |
| MDT829337C | Three-Piece Swedish Belt      | ea   |

### Wrap-around Belts with Hook-and-Loop Closure

| Item No.   | Size                    | Max Length | Pkg. |
|------------|-------------------------|------------|------|
| MDT3658ML  | Med/Large, Beige Fabric | 77"        | ea   |
| MDT3658XL  | X-Large, Beige Fabric   | 96"        | ea   |
| MDT36582XL | X-Large, Beige Fabric   | 121"       | ea   |
| MDT36592XL | 2X-Large, Foam Fabric   | 121"       | ea   |

### Roll Belts with Tie Closure

| Item No.  | Description             | Pkg. |
|-----------|-------------------------|------|
| MDT821231 | Single Strap            | ea   |
| MDT821135 | Roll Belt, Double Strap | ea   |

### Roll Belts with Quick Release Buckles

| Item No.   | Description  | Pkg. |
|------------|--------------|------|
| MDT821231Q | Single Strap | ea   |
| MDT821135Q | Dual Strap   | ea   |

### Safety/Soft Belts

| Item No.  | Description        | Pkg. |
|-----------|--------------------|------|
| MDT822125 | Poly, 1pc          | ea   |
| MDT822126 | Poly, 2pc          | ea   |
| MDT822127 | Quick-Release, 2pc | ea   |

## Breezeline Pelvic Holders

- » Mesh holder helps prevent forward sliding in wheelchairs or geri-chairs
- » Straps have quick-release buckles
- » Machine washable

| Item No. | Description | L x W      | Upper Waist Strap | Lower Waist Strap | Pkg. |
|----------|-------------|------------|-------------------|-------------------|------|
| PSY4427S | Small       | 16.5 x 17" | 52.5"             | 50"               | ea   |
| PSY4427M | Medium      | 19 x 17"   | 56.5"             | 65"               | ea   |
| PSY4427L | Large       | 23 x 17"   | 60.5"             | 56"               | ea   |



MDT3658ML



MDT821231



MDT822125



PSY4427L





MDT828306L



MDT828002

## Quick Release Blended Safety Vests and Jackets

- » Keep patient securely in bed or in an upright position in a wheelchair
- » Vests are available with quick-release straps or ties and a zip back
- » Long sleeve jacket has an off-center zipper back to avoid spinal irritation
- » Machine washable cotton/polyester

### Vests—with quick-release straps, zip back

| Item No.   | Size                 | Pkg. |
|------------|----------------------|------|
| MDT828306M | Medium, 110–155 lbs. | ea   |
| MDT828306L | Large, 145–190 lbs.  | ea   |

### Vests—with tie straps, zip back

| Item No.    | Size                  | Pkg. |
|-------------|-----------------------|------|
| MDT828305S  | Small, 80–120 lbs.    | ea   |
| MDT828305M  | Medium, 110–155 lbs.  | ea   |
| MDT828305L  | Large, 145–190 lbs.   | ea   |
| MDT828305XL | X-Large, 180–230 lbs. | ea   |

### Jackets

| Item No.    | Size                 | Pkg. |
|-------------|----------------------|------|
| MDT8283360S | Small, 80–120 lbs.   | ea   |
| MDT8283360M | Medium, 110–155 lbs. | ea   |
| MDT8283360L | Large, 145–190 lbs.  | ea   |

## Tie-Back Safety Vests

- » Vests prevent patient from sitting up or sliding down in bed or wheelchair
- » Ties in back out of patient's reach with connecting straps that attach to bed frame or wheelchair tilt bars
- » Koolnit polyester mesh is comfortable and washable

| Item No.  | Size                 | Pkg. |
|-----------|----------------------|------|
| MDT828000 | Small, 80–120 lbs.   | 6/cs |
| MDT828002 | Medium, 110–155 lbs. | 6/cs |
| MDT828004 | Large, 145–190 lbs.  | 6/cs |

## Pulse Oximeter Probe Wraps

- » Help secure optical transducers and probe sensors without adhesive tape
- » Hold optical components firmly to help prevent misalignment
- » Universal design may be trimmed to fit neonatal foot, infant toe or pediatric finger

| Item No. | L x W         | Pkg.  |
|----------|---------------|-------|
| PSY6554  | 5.25" x 1.25" | 12/ea |
| PSY6554H | 5.25" x 1.25" | ea    |



PSY6554

## Palm Cones

- » Help provide anti-contracture hand therapy and prevent skin breakdown without restricting finger movement
- » Act as a diversionary device for patients prone to pulling at IV or feeding tubes
- » Use for exercise or to help flexion contractures of the hand
- » Highly absorbent terry cloth conforms to hand's natural shape
- » Adjustable foam strap

| Item No.  | Size          | Pkg. |
|-----------|---------------|------|
| MDT821300 | 2.5"W x 4.5"L | ea   |
| PSY5609   | 5"L           | ea   |



MDT821300

## HEELMEDIX® Heel Protectors

- » Off-load the heel and help redistribute pressure
- » Open heel design allows for airflow and easy monitoring
- » Absorbent inner lining helps keep skin dry, reduces friction and shear on heels
- » Help protect against foot drop
- » Optional wedge can be used to help prevent medial and lateral rotation

| Item No.    | Size   | Pkg. |
|-------------|--|------|
| MDT823330P  | Pediatric Petite, Slip-on, 5"-11" calf circumference | ea   |
| MDT823330   | One Size, Hook and Loop                              | ea   |
| MDT823330CS | One Size, Hook and Loop                              | 4/ea |
| MDT823330SW | One Size, Hook and Loop with Wedge                   | ea   |
| MDT823330XL | X-Large, Slip-on with Wedge, 19"-25" calf circum.    | ea   |
| MDT823330W  | Stabilizer Wedge                                     | 2/ea |



MDT823330

## Comfy Seat

- » Microfiber filling provides seating comfort
- » Covers gaps in wheelchair openings to help prevent pinching
- » Easily attaches to wheelchairs or geri-chairs with ties
- » Reinforced slits allow use of a positioning belt
- » Machine washable

| Item No. | Seat        | Back        | Side Panels | Pkg. |
|----------|-------------|-------------|-------------|------|
| PSY6526G | 22"W x 20"D | 22"W x 31"H | 20"D x 20"H | ea   |



PSY6526G



**MODERATE  
FALL  
RISK**



**Medline Industries, Inc.** Three Lakes Drive, Northfield, IL 60093

Medline United States  
1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

Medline Canada  
1-800-396-6996  
medline.ca | canada@medline.com

Medline México  
01-800-831-0898  
medlinemexico.com | mexico@medline.com

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this catalog. © 2018 Medline Industries Inc. © 2018 Stand FIRM Falls Management Program. Advantage, EqualizeAire, Medline, Safe-N-Sound, Sensor Pad, Smart Caregiver, Spil-Pruf, Supra, TheraTech, UltraFlex, Quick Alert are registered trademarks of Medline Industries, Inc. PerforMAX and Stand FIRM are trademarks of Medline Industries, Inc. Velcro is a registered trademark of Velcro Industries. Veriflex is a registered trademark of Boston Scientific Corporation. All rights reserved. MKT1899478 / LIT191R / 5M / PP / 39