

Venosan® USA

English

Indications: (20-30 mmHg)
Mild edema, slight varicose veins, to prevent varicose veins during pregnancy

Indications: (30-40 mmHg)
Varicose veins with tendency to edema, chronic venous insufficiency, after vein surgery, varicose veins during pregnancy

- Graduated compression
- Quality without compromise for the cost conscious patient


Español

Indicaciones: (20-30 mmHg)
Edema leve, venas varicosas pequeñas, prevención de venas varicosas durante el embarazo

Indicaciones: (30-40 mmHg)
Venas varicosas con tendencia al edema, insuficiencia venosa crónica, después de una intervención quirúrgica de venas, venas varicosas durante el embarazo

- Compresión graduada
- Para pacientes concientes de la economía y que no quieren renunciar a la calidad

MEASUREMENT SIZE GUIDE for custom fit to your exact measurements:



Circumference measurements	SMALL	MEDIUM	LARGE	X-LARGE
WAIST	up to 30" up to 76 cm	up to 33.5" up to 85 cm	up to 37" up to 94 cm	up to 40.5" up to 103 cm
HIP	up to 41" up to 104 cm	up to 45" up to 115 cm	up to 48.5" up to 123 cm	up to 52" up to 132 cm
MID-THIGH	16" - 20.5" 41 - 52 cm	20.5" - 24.5" 51 - 62 cm	24" - 29.5" 61 - 75 cm	29" - 33.5" 74 - 85 cm
BELOW KNEE	14" and over 36 cm and over	14" and over 36 cm and over	14" and over 36 cm and over	14" and over 36 cm and over
MID-CALF	11 - 13.5" 28 - 35 cm	13 - 15.5" 33 - 40 cm	15" - 17.5" 38 - 45 cm	17" - 19.5" 43 - 50 cm
ANKLE	6.25" - 7.25" 16 - 19 cm	7.5" - 9" 19 - 23 cm	9.25" - 10.75" 23 - 28 cm	11" - 12.5" 28 - 32 cm

STYLE	COLOR	SIZE	TOE	20-30 mmHg
BELOW KNEE	BEIGE	SMALL	OPEN	2030112
BELOW KNEE	BEIGE	MEDIUM	OPEN	2030113
BELOW KNEE	BEIGE	LARGE	OPEN	2030114
BELOW KNEE	BEIGE	X-LARGE	OPEN	2030115
BELOW KNEE	BEIGE	SMALL	CLOSED	2030122
BELOW KNEE	BEIGE	MEDIUM	CLOSED	2030123
BELOW KNEE	BEIGE	LARGE	CLOSED	2030124
BELOW KNEE	BEIGE	X-LARGE	CLOSED	2030125
BELOW KNEE	BLACK	SMALL	CLOSED	2030122BL
BELOW KNEE	BLACK	MEDIUM	CLOSED	2030123BL
BELOW KNEE	BLACK	LARGE	CLOSED	2030124BL
BELOW KNEE	BLACK	X-LARGE	CLOSED	2030125BL
MID-THIGH	BEIGE	SMALL	OPEN	2030212
MID-THIGH	BEIGE	MEDIUM	OPEN	2030213
MID-THIGH	BEIGE	LARGE	OPEN	2030214
MID-THIGH	BEIGE	X-LARGE	OPEN	2030215
MID-THIGH	BEIGE	SMALL	CLOSED	2030222
MID-THIGH	BEIGE	MEDIUM	CLOSED	2030223
MID-THIGH	BEIGE	LARGE	CLOSED	2030224
MID-THIGH	BEIGE	X-LARGE	CLOSED	2030225
MID-THIGH W/ SILICONE TOP	BEIGE	SMALL	CLOSED	2030222H
MID-THIGH W/ SILICONE TOP	BEIGE	MEDIUM	CLOSED	2030223H
MID-THIGH W/ SILICONE TOP	BEIGE	LARGE	CLOSED	2030224H
MID-THIGH W/ SILICONE TOP	BEIGE	X-LARGE	CLOSED	2030225H
PANTYHOSE	BEIGE	SMALL	CLOSED	2030322
PANTYHOSE	BEIGE	MEDIUM	CLOSED	2030323
PANTYHOSE	BEIGE	LARGE	CLOSED	2030324
PANTYHOSE	BEIGE	X-LARGE	CLOSED	2030325
PLUS SIZE PANTYHOSE	BEIGE	SMALL	CLOSED	2030422
PLUS SIZE PANTYHOSE	BEIGE	MEDIUM	CLOSED	2030423
PLUS SIZE PANTYHOSE	BEIGE	LARGE	CLOSED	2030424
PLUS SIZE PANTYHOSE	BEIGE	X-LARGE	CLOSED	2030425

STYLE	COLOR	SIZE	TOE	30-40 mmHg
BELOW KNEE	BEIGE	SMALL	OPEN	3040112
BELOW KNEE	BEIGE	MEDIUM	OPEN	3040113
BELOW KNEE	BEIGE	LARGE	OPEN	3040114
BELOW KNEE	BEIGE	X-LARGE	OPEN	3040115
BELOW KNEE	BEIGE	SMALL	CLOSED	3040122
BELOW KNEE	BEIGE	MEDIUM	CLOSED	3040123
BELOW KNEE	BEIGE	LARGE	CLOSED	3040124
BELOW KNEE	BEIGE	X-LARGE	CLOSED	3040125
MID-THIGH	BEIGE	SMALL	OPEN	3040212
MID-THIGH	BEIGE	MEDIUM	OPEN	3040213
MID-THIGH	BEIGE	LARGE	OPEN	3040214
MID-THIGH	BEIGE	X-LARGE	OPEN	3040215
MID-THIGH W/ SILICONE TOP	BEIGE	SMALL	OPEN	3040212H
MID-THIGH W/ SILICONE TOP	BEIGE	MEDIUM	OPEN	3040213H
MID-THIGH W/ SILICONE TOP	BEIGE	LARGE	OPEN	3040214H
MID-THIGH W/ SILICONE TOP	BEIGE	X-LARGE	OPEN	3040215H
PANTYHOSE	BEIGE	SMALL	OPEN	3040312
PANTYHOSE	BEIGE	MEDIUM	OPEN	3040313
PANTYHOSE	BEIGE	LARGE	OPEN	3040314
PANTYHOSE	BEIGE	X-LARGE	OPEN	3040315



Fiber Content 2030 75% Nylon, 25% Lycra® spandex, 3040 70% Nylon, 30% Lycra® spandex • Latex Free • Made in U.S.A.
©2003 SALZMANN AG. VENOSAN® is a registered trademark of SALZMANN AG St. Gallen, Switzerland