



| Re-Order Number | Description                    | Size             | Units |
|-----------------|--------------------------------|------------------|-------|
| 503150          | Heels-Off (Medium-Large)       | 24" x 16" x 7"   | Each  |
| 503151          | Heels-Off (Small)              | 24" x 13" x 5.5" | Each  |
| 503165          | Heels-Off w/Gel (Medium-Large) | 24" x 16" x 7"   | Each  |
| 503166          | Heels-Off w/Gel (Small)        | 24" x 13" x 5.5" | Each  |

### Purpose

Heels-Off is designed to off-load both heels to prevent pressure sores from forming and can aid in the healing of existing ulcers. It elevates the legs to relax the spine and prevents leg and lower back pain.

In addition to the benefits of the standard Heels Off, the new Heels Off w/Gel has a surface layer of gel to provide enhanced pressure relief and therapeutic cooling for the lower leg.

### Easy to Use

Place the higher end of the Heels-Off device under the patient's knees while suspending the heels

### Easy Care

The Low-Shear II cover could be wiped clean with a spray and wipe cleanser. If disinfected, avoid using chlorine (bleach) chemicals. Do not launder.

[Click here to view detailed product descriptions and prices.](#)

Call 800-397-5899 to order today!