



Kangaroo™ Connect Enteral Feeding Pump

The Kangaroo Connect enteral feeding pump can enhance quality of life for patients through best-in-class pump performance

The Kangaroo Connect pump is the...

- Smallest
- Lightest
- Most Waterproof

...enteral pump available on the market.

The Kangaroo Connect pump has an:

- Accurate alarm system
- Ergonomic design

Feature Comparison Kangaroo™: Always Improving

40 years of research resulted in the features of the Kangaroo Connect enteral feeding pump.

	Kangaroo™ Joey Enteral Feeding Pump	Kangaroo™ Connect Enteral Feeding Pump	Zevex Infinity Orange™* Enteral Feeding Pump
Accuracy +/-	7%	5%	5%
Weight	1.7 lbs	0.7 lbs	0.9 lbs
IPX Rating	4	6	5
Battery Life	8 hrs	24 hrs	24 hrs
72 Hr. History	•	•	
Auto-Prime	•	•	



Imagine a pump with:

- No fragile door
- A color screen
- 2 sensors for improved alarm management
- A full service sales team and dedicated repair facility

For more information please visit
Covidien.com/KangarooConnect

A pump that could improve your patients' quality of life:
This is Kangaroo™.

Kangaroo™ Connect Enteral Feeding Pump

FEEDING SET

- Cassette attachment allows quick and easy feeding bag changes.
- “No drip chamber” design helps to reduce clogging
- Simple, one-handed loading of pump cassette
- ENFit connection



ORDERING INFORMATION

Item Code	Description	Ship Case
384400	Kangaroo™ Connect Enteral Feeding Pump	1
384491	Kangaroo™ Connect Power Cord with Adapter US	1
384492	Kangaroo™ Connect Pole Clamp	1
384494	Kangaroo™ Connect Car Charger	1

Pump Sets

77000FD	Kangaroo™ Connect Safety Screw Spike Set Non-Sterile	30
77100FD	Kangaroo™ Connect Feeding Set 1000 mL Non-Sterile	30
77150FD	Kangaroo™ Connect Feeding Set 1600 mL Non-Sterile	30
77500FD	Kangaroo™ Connect Feeding Set 500 mL Non-Sterile	30

Backpacks

770035S	Kangaroo™ Connect Backpack, Adjustable Straps, Purple, Small	1
770035M	Kangaroo™ Connect Backpack, Adjustable Straps, Purple, Medium	1
770035L	Kangaroo™ Connect Backpack, Adjustable Straps, Purple, Large	1
770037S	Kangaroo™ Connect Backpack, Adjustable Straps, Black, Small	1
770037M	Kangaroo™ Connect Backpack, Adjustable Straps, Black, Medium	1
770037L	Kangaroo™ Connect Backpack, Adjustable Straps, Black, Large	1
770036S	Premium Kangaroo™ Connect Backpack, Black, Small	1
770036M	Premium Kangaroo™ Connect Backpack w/ Tabletop Stand, Black, Medium	1
770036L	Premium Kangaroo™ Connect Backpack w/ Tabletop Stand, Black, Large	1

For additional information on Nutritional Delivery Solutions please visit:

Covidien.com/kangarooconnect

Covidien.com/patientcare-edu

Covidien.com/medicalsupplies

Pending 510(k) clearance, these products are not currently available for sale in the U.S.

COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG.

™* Trademark of its respective owner.

Other brands are trademarks of a Covidien company. © 2015 Covidien.

H8873 – 200 – 0515

