

Mask Quick Reference Guide



Mask Quick Reference Guide



Nasal Masks and Spare Parts

F&P Eson™ 2

F&P Eson™

F&P Zest™ Q

Zest™

FlexiFit™ 407

FlexiFit™ 406 Petite

FlexiFit™ 405

Full Face Masks and Spare Parts

F&P Simplus™

Forma™

FlexiFit™ 431

FlexiFit™ 432

Nasal Pillows Masks and Spare Parts

F&P Brevida™

F&P Pilairo™ Q

Opus™ 360

Oral Mask and Spare Parts

Oracle™

F&P Eson™ 2



Confidence comes easy with F&P Eson 2. To improve on the F&P Eson – designed for performance and comfort – our development team focused on ease of use.

- **RollFit™ Seal**

The F&P Eson 2 retains the RollFit Seal, designed to deliver an effective seal and lessen the likelihood of leaks. Available in three seal sizes.

- **Intuitive Headgear**

This high quality headgear provides multiple fitting options for patients, including the ability to stretch the mask directly over the head.

- **Easy Frame**

The Easy Frame – fits all three seal sizes. Now has a very low profile which provides a clear line of sight while retaining seal stability.

- **VisiBlue™**

Blue highlights incorporated into key mask components assist patients with mask use, orientation, disassembly and reassembly.








- **Washable Diffuser**

The innovative diffuser, which reduces noise and draft, is now washable and more durable.

Product Code	Description	Quantity	HCPCS Code
ESN2SA	Eson 2 Small	1	A7034 + A7035
ESN2MA	Eson 2 Medium	1	
ESN2LA	Eson 2 Large	1	

F&P Eson™ 2

The F&P Eson 2 Nasal Mask spare parts are listed below.

Product Code	Description	HCPCS Code	
400ESN221* 400ESN222*	Eson 2 Headgear Small Eson 2 Headgear Medium-Large	A7035	
400ESN214* 400ESN215* 400ESN216*	Eson 2 Seal Small Eson 2 Seal Medium Eson 2 Seal Large	A7032	
400FPH241	Female Swivel 2		
400ESN251	Eson 2 Clips		
400ESN261	Washable Diffuser		
400ESN271	Eson 2 Elbow		
400ESN231* 400ESN232* 400ESN233*	Mask Kit - Small (no headgear) Mask Kit - Medium (no headgear) Mask Kit - Large (no headgear)	A7034	

* Reimbursable replacement parts

F&P Eson™



The F&P Eson is designed to perform in tune with you and your needs. The three key components, the *RollFit Seal*, *ErgoFit Headgear* **and** *Easy Frame*, work in harmony. These components have been designed to deliver the comfort, seal and easy use that Fisher & Paykel Healthcare masks are known for.

- **RollFit Seal**

As the name suggests, the one-piece seal ‘rolls’ back and forth on the bridge of the nose to adjust automatically. In doing so, the RollFit Technology reduces pressure on the bridge of the nose.

- **ErgoFit Headgear**

This breathable, ergonomically designed headgear has been designed to self-locate high on the rear of the head, allowing for head movement (sideways, up and down).

- **Easy Frame**

This low-profile frame is stable, durable and small, and ensures a clear line of sight. The one frame fits all three seal sizes and has an Easy-Clip frame attachment to aid assembly after cleaning.

Product Code	Description	Quantity	HCPSC Code
400449	Eson Mask Small		A7034 +
400450	Eson Mask Medium		A7035
400451	Eson Mask Large		
	Complete system includes:		
	• RollFit Seal	1	
	• ErgoFit Headgear	1	
	• Diffuser Material plus spares	6	
	• Q Cover	1	


F&P Eson™

The F&P Eson Nasal Mask spare parts are listed below.

Product Code	Description	HPCPS Code	
400HC567* 400HC568*	Headgear – Small Headgear – Medium/Large	A7035	
400HC126* 400HC127* 400HC128*	Seal – Small Seal – Medium Seal – Large	A7032	
400HC228	Diffuser – 10 pack + 1 Q Cover		
400HC569	Buckle and Clips		
400HC570	Swivel		
400HC571* 400HC572* 400HC573*	Mask Kit – Small (no headgear) Mask Kit – Medium (no headgear) Mask Kit – Large (no headgear)	A7034	







* Reimbursable replacement parts

F&P Zest™ Q

Product Code	Description	HCPCS Code	
400444 400445 400446	F&P Zest Q Petite F&P Zest Q F&P Zest Q Plus Complete system includes: <ul style="list-style-type: none"> • Easy-Clip Silicone Seal • FlexiFoam™ Cushion • Stretchgear™ Headgear • Glider™ Strap • Exhaust Flow Diffuser Cover • Exhaust Flow Diffuser Material – 5 pack + 1	A7034 + A7035	
400HC557* 400HC542* 400HC558*	FlexiFoam Cushion and Silicone Seal – Petite FlexiFoam Cushion and Silicone Seal FlexiFoam Cushion and Silicone Seal – Plus	A7032	
400HC314* 400HC316*	Headgear Headgear – Petite	A7035	
400HC225	Exhaust Flow Diffuser Cover		
400HC226	Exhaust Flow Diffuser Material – 10 pack		
400HC227	Elbow (includes Swivel, Cover and Diffuser Material)		
900HC418	Glider Strap		
400HC576* 400HC577* 400HC578*	Mask Kit – Petite (no headgear) Mask Kit – Standard (no headgear) Mask Kit – Plus (no headgear)	A7034	

* Reimbursable replacement parts

Zest™

Product Code	Description	HCPSC Code	
400439 400440 400441	Zest Petite Zest Zest Plus Complete system includes: <ul style="list-style-type: none"> • Easy-Clip Silicone Seal • FlexiFoam Cushion • Stretchgear Headgear • Glider Strap 	A7034 + A7035	
400HC542* 400HC557* 400HC558*	FlexiFoam Cushion and Silicone Seal FlexiFoam Cushion and Silicone Seal – Petite FlexiFoam Cushion and Silicone Seal – Plus	A7032	
400HC314* 400HC316*	Headgear Headgear – Petite	A7035	
400HC221	Elbow plus Swivel		
900HC418	Glider Strap		
400HC559* 400HC543* 400HC560*	Mask Kit – Petite (no headgear) Mask Kit – Standard (no headgear) Mask Kit – Plus (no headgear)	A7034	



* Reimbursable replacement parts

FlexiFit™ 407

Product Code	Description	HCPCS Code	
HC407	FlexiFit 407 Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (one size) • Foam Cushion • Headgear • Glider Strap 	A7034 + A7035	
400HC501*	Foam Cushion and Silicone Seal	A7032	
400HC301*	Headgear (one size)	A7035	
400HC202	Exhaust Flow Elbow		
900HC418	Glider Strap		
400HC502*	Mask Kit (no headgear)	A7034	









* Reimbursable replacement parts

FlexiFit™ 406

Product Code	Description	HCPCS Code	
HC406	FlexiFit 406 Petite Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (one size) • Foam Cushion • Headgear • Glider Strap 	A7034 + A7035	
400HC509*	Foam Cushion and Silicone Seal	A7032	
400HC307*	Headgear (one size)	A7035	
400HC202	Exhaust Flow Elbow		
900HC418	Glider Strap		
400HC510*	Mask Kit (no headgear)	A7034	

* Reimbursable replacement parts

FlexiFit™ 405

Product Code	Description	HCPCS Code	
HC405	FlexiFit 405 Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (Small) • Silicone Seal (Large) • Foam Cushion • Headgear (one size) • Exhaust Flow Diffuser Cap plus Spare • Exhaust Flow Diffuser Material plus Spares • Glider Strap 	A7034 + A7035	
900HC427* 900HC428*	Foam Cushion and Silicone Seal – Small Foam Cushion and Silicone Seal – Large	A7032	
900HC429*	Headgear (one size)	A7035	
900HC439	Exhaust Flow Diffuser Material – 10 pack		
900HC440	Exhaust Flow Diffuser Cap plus Spare		
900HC441	Exhaust Flow Elbow and Diffuser Material with Spares		
900HC418	Glider Strap		
900HC406*	Mask Kit (no headgear)	A7034	

* Reimbursable replacement parts

F&P Simplus™








F&P Simplus

When engineering the F&P Simplus our designers set out to create a mask that revolutionized Full Face comfort, seal and easy use. The F&P Simplus is the result of intensive research and development, many prototypes and thousands of design hours. We consulted patients, equipment providers and clinicians from around the world to develop a mask that simply fits + performs.

- RollFit Seal**
As the name suggests, the one-piece seal “rolls” back and forth on the bridge of the nose to adjust automatically. In doing so, the **RollFit Technology** reduces pressure on the bridge of the nose.
- ErgoForm Headgear**
This breathable, ergonomically designed headgear self-locates high on the rear of the head, allowing for unrestricted head movement (sideways, up and down).
- Easy Frame**
This low-profile frame is stable, durable and small, and ensures a clear line of sight. The one frame fits all three seal sizes and has an **Easy-Clip Frame Attachment** to aid assembly after cleaning.






Product Code	Description	Quantity	HCPSC Code
400475	F&P Simplus Small		A7030 +
400476	F&P Simplus Medium		A7035
400477	F&P Simplus Large		
	Complete system includes:		
	• RollFit™ Seal	1	
	• ErgoForm™ Headgear	1	
	• Easy Frame	1	
	• Swivel	1	

F&P Simplus™

Product Code	Description	HCPCS Code	
400HC579* 400HC580* 400HC581*	F&P Simplus RollFit Seal - Small F&P Simplus RollFit Seal - Medium F&P Simplus RollFit Seal - Large	A7031	
400HC582* 400HC583*	ErgoForm Headgear - Small ErgoForm Headgear - Medium/Large	A7035	
400HC587	F&P Simplus Headgear Clips + Buckle		
400HC570	Female Swivel		
400HC584* 400HC585* 400HC586*	F&P Simplus Mask Kit - No Headgear Small F&P Simplus Mask Kit - No Headgear Medium F&P Simplus Mask Kit - No Headgear Large	A7030	






* Reimbursable replacement parts

Forma™

Product Code	Description	HCPCS Code	
400470 400471 400473	Forma – Small Forma – Medium & Large Forma – X-Large Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (Small or Medium & Large or X-Large) • FlexiFoam Cushion (Small or Medium & Large or X-Large) • Headgear with Easy-Release Cord • Glider Strap 	A7030 + A7035	
400HC121* 400HC122* 400HC123* 400HC124*	FlexiFoam Cushion and Silicone Seal – Small FlexiFoam Cushion and Silicone Seal – Medium FlexiFoam Cushion and Silicone Seal – Large FlexiFoam Cushion and Silicone Seal – X-Large	A7031	
400HC315*	Headgear with Crown Strap (one size)	A7035	
400HC203	Elbow and Vented Non-Rebreathing Valve		
400HC552* 400HC553* 400HC554* 400HC555*	Small Mask Kit (no headgear) Medium Mask Kit (no headgear) Large Mask Kit (no headgear) X-Large Mask Kit (no headgear)	A7030	






* Reimbursable replacement parts

FlexiFit™ 431

Product Code	Description	HCPCS Code	
HC431	FlexiFit 431 Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (Small) • Silicone Seal (Medium) • Silicone Seal (Large) • Headgear with Easy-Release Cord • Glider Strap 	A7030 + A7035	
400HC112* 400HC104* 400HC103*	Silicone Seal – Small Silicone Seal – Medium Silicone Seal – Large	A7031	
400HC302*	Headgear with Crown Strap (one size)	A7035	
400HC203	Elbow and Vented Non-Rebreathing Valve		
400HC503*	Mask Kit (no headgear)	A7030	

* Reimbursable replacement parts

FlexiFit™ 432

Product Code	Description	HCPCS Code	
HC432S HC432M HC432L HC432XL	FlexiFit 432 – Small FlexiFit 432 – Medium FlexiFit 432 – Large FlexiFit 432 – X-Large Complete system includes: <ul style="list-style-type: none"> • Silicone Seal (Small or Medium or Large or X-Large) • Foam Cushion (Small or Medium or Large or X-Large) • Headgear with Easy-Release Cord • Glider Strap 	A7030 + A7035	
400HC113* 400HC114* 400HC115* 400HC119*	Foam Cushion and Silicone Seal – Small Foam Cushion and Silicone Seal – Medium Foam Cushion and Silicone Seal – Large Foam Cushion and Silicone Seal – X-Large	A7031	
400HC309*	Headgear with Crown Strap (one size)	A7035	
400HC203	Elbow and Vented Non-Rebreathing Valve		
400HC519* 400HC521* 400HC522* 400HC541*	Small Mask Kit (no headgear) Medium Mask Kit (no headgear) Large Mask Kit (no headgear) X-Large Mask Kit (no headgear)	A7030	

* Reimbursable replacement parts

F&P Brevida™



Designed for Confidence. Built for Comfort.

Gaining patient confidence is key to successful CPAP therapy and patient confidence begins with a mask that fits and is comfortable. Developed using extensive patient centred research, F&P Brevida features:

- **Adjustable Headgear**
Simple, adjustable headgear to deliver an individual fit and designed to stay in place during sleep.
- **Innovative AirPillow™ Seal**
The AirPillow seal inflates to form a ‘pillow’ of air in and around the nose for a gentle, effective seal. The sealing region is made of extra thin silicone which is gentle on the nose and includes side supports to reduce the likelihood of the seal dislodging during sleep.
- **VisiBlue™**
VisiBlue highlights are blue colour cues, incorporated into the key mask components, which assist patients with mask use, orientation, disassembly and re-assembly. VisiBlue helps to make mask education simpler and more intuitive.
- **Washable Diffuser**
The diffuser – which reduces noise and draft – is now washable and reusable.

All seal sizes (XS-S and M-L) fit on one mask frame.

Product Code	Description	HCPCS Code
BREISMA	Brevida Nasal Pillows Mask (XS-S and M-L seals)	A7034 + A7035

F&P Brevida™

Product Code	Description	HCPCS Code	
400BRE113*	Brevida Seal Extra Small-Small	A7033	
400BRE114*	Brevida Seal Medium-Large	A7033	
400BRE120*	Brevida Headgear	A7035	
400BRE151	Brevida Headgear Clips		
400BRE161	Brevida Diffuser		
400BRE171	Brevida Elbow		
400FPH241	Female Swivel 2		
400BRE131* 400BRE132*	Brevida Mask Kit (no headgear) XS-S Brevida Mask Kit (no headgear) M-L	A7034	

* Reimbursable replacement parts

F&P Pilairo™ Q





Weighing in as one of the lightest nasal pillows masks (1.94 ounces/55 grams*), the F&P Pilairo Q has been further developed and re-engineered to improve patient satisfaction. The F&P Pilairo Q has just three key components and two minimalist headgear options.

- **AirPillow™ Seal**
The AirPillow seal is 3D contoured with fine Variable Thickness Silicone (VTS™) and uses CPAP pressure to self-inflate to create an effective yet gentle double seal into and around the nose.
- **Easy Frame**
This low-profile frame is stable, durable and small, and ensures a clear line of sight. The delivery tube is lightweight, very flexible and designed to minimize pressure or pull on the face.
- **Headgear Options: Adjustable or Stretchwise™ Headgear**
 - a. **Adjustable Headgear** – This breathable, ergonomically designed headgear fits a greater range of head sizes and enables sizing adjustment to the patient's preference.
 - b. **Stretchwise Headgear** – This minimal headgear automatically adjusts to the wearer for a quick and easy fit.

Product Code	Description	Quantity	HCPCS Code
400421	F&P Pilairo Q Nasal Pillows Mask Complete system includes: <ul style="list-style-type: none">• Adjustable Headgear• Stretchwise Headgear• AirPillow Seal• Flexible Tube	1 1 1 1	A7034 + A7035







F&P Pilairo™ Q

Product Code	Description	HCPCS Code	
400HC125*	AirPillow Seal	A7033	
400HC328†	Adjustable and Stretchwise Headgear – Combo pack	A7035	
400HC229	Elbow Cover and Diffuser – 10 pack		
400HC224	Female Swivel		
400HC565*	Pilairo Q Nasal Pillows Mask Kit (no headgear)	A7034	








† F&P Pilairo Q with Stretchwise Headgear

* Reimbursable replacement parts

Opus™ 360

Product Code	Description	HCPCS Code	
HC482	Opus 360 Complete system includes: <ul style="list-style-type: none"> • Mask Base and Flexitube Assembly • Headgear • Silicone Seal (Small, Medium and Large) • Male Swivel • Spare Tube-Anchoring Strap 	A7034 + A7035	
400HC116* 400HC117* 400HC118*	Silicone Pillow – Small Silicone Pillow – Medium Silicone Pillow – Large	A7033	
400HC310*	Headgear with Tube-Anchoring Strap (one size)	A7035	
400HC311	Tube-Anchoring Strap		
400HC312	Headgear – Tri-Glide Buckle		
400HC526*	Opus 360 Nasal Pillows Mask Kit (no headgear)	A7034	

* Reimbursable replacement parts

Product Code	Description	HCPCS Code	
HC452	<p>Oracle</p> <p>Complete system includes:</p> <ul style="list-style-type: none"> • Flexitube Assembly • Headgear with Easy-Release Cord • Silicone Seal (Small and Large) • Storage Case 	A7044 + A9999	
400HC105 400HC106	<p>Silicone Seal – Small</p> <p>Silicone Seal – Large</p>		
400HC300	Headgear Accessory plus Clip	A9999	
900HC439	Exhaust Flow Diffuser Material – 10 pack		
900HC440	Exhaust Flow Diffuser Cap plus Spare		
900HC441	Elbow Kit		
400HC204	Nasal Plugs – Small, Medium and Large		

Notes



