

Tubular Elastic Dressing Retainers

McKesson Tubular Elastic Dressing Retainers are made from an elastic net material. They securely hold various types of dressings with no adhesive needed. They come in a variety of sizes to fit any shape (sizes 1-11) and hold 30 yards per box.

Features and Benefits

- Designed to serve as a secondary dressing
- Applies gentle pressure to keep bandages and medical devices securely in place without adhesives
- Easy to customize and conforms to body contours
- Allows flexibility and maximum air flow
- Available in a wide range of sizes to suit a variety of patient needs
- Not made with natural rubber latex



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Tubular Elastic Dressing Retainer, Size 1, 30 yds.	1016	1ea/bx, 10bx/cs
Tubular Elastic Dressing Retainer, Size 2, 30 yds.	1017	1ea/bx, 10bx/cs
Tubular Elastic Dressing Retainer, Size 3, 30 yds.	1018	1ea/bx, 10bx/cs
Tubular Elastic Dressing Retainer, Size 4, 30 yds.	1019	1ea/bx, 10bx/cs
Tubular Elastic Dressing Retainer, Size 5, 30 yds.	1020	1ea/bx, 10bx/cs
Tubular Elastic Dressing Retainer, Size 6, 30 yds.	1021	1ea/bx, 10bx/cs
Tubular Elastic Dressing Retainer, Size 7, 30 yds.	1022	1ea/bx, 10bx/cs
Tubular Elastic Dressing Retainer, Size 8, 30 yds.	1023	1ea/bx, 10bx/cs
Tubular Elastic Dressing Retainer, Size 9, 30 yds.	1024	1ea/bx, 10bx/cs
Tubular Elastic Dressing Retainer, Size 10, 30 yds.	1025	1ea/bx, 10bx/cs
Tubular Elastic Dressing Retainer, Size 11, 30 yds.	1026	1ea/bx, 10bx/cs