



**DRIVE
COMPLIANCE.**

Isolation Gowns
AAMI Level Protection

PROTECT YOUR PEOPLE AND YOUR BOTTOM LINE.

You want to know you're protected. Our protective apparel goes through rigorous testing and meets strict FDA guidelines. We follow the highest quality standards to make sure you get the protection you need.

As a leader in protective apparel, we also have the most comprehensive product line and low acquisition costs. You can find the exact products you need—and maintain cost efficiencies.

We Follow High Quality Standards

Our isolation gowns meet strict guidelines developed by the Association for the Advancement of Medical Instrumentation (AAMI). The chart below describes the criteria for AAMI Levels 1, 2 and 3. Each level is assigned a different color. Choose the right gown by matching the color on the neck binding to the AAMI level on the chart.

Levels of Barrier Protection—AAMI PB70:2012¹

	LEVEL 1	LEVEL 2		LEVEL 3	
Test	AATCC 42:2000	AATCC 42:2000	AATCC 127:1998	AATCC 42:2000	AATCC 127:1998
Requirements at 4% AQL	Water Impact ≤ 4.5 g	Spray Impact ≤ 1.0 g	Hydrostatic Pressure ≥ 20 cm	Spray Impact ≤ 1.0 g	Hydrostatic Pressure ≥ 50 cm
Anticipated Risk of Exposure	Low	Moderate		Moderate to High	



Color-coded neck binding for easy identification



NONLV100

NONLV200

NONLV325

AAMI Level 1 Isolation Gown

Recommended for: Med/Surg, Visitors, Contact Precautions

- » Made from multi-ply material
- » Side ties eliminate the difficulties of tying behind the back
- » Choose between over-the-head or neck tie styles

Item No.	Description	Pkg.
Medium Weight, Neck Tie, Elastic Wrists		
NONLV100	Regular/Large, Yellow	100/cs
NONLV100XL	X-Large, Yellow	100/cs
Medium Weight, Over-the-Head, Thumb Loop Wrists		
NONLV140	Regular/Large, Blue	100/cs
NONLV140XL	X-Large, Yellow	100/cs

AAMI Level 2 Isolation Gown

Recommended for: Med/Surg, ICU, Emergency, Airborne and Contact Precautions

- » Made from multi-ply material with side ties to eliminate the difficulties of tying behind the back
- » Choice of over-the-head and tape tab neck styles
- » Available with elastic wrists or thumb loops

Item No.	Description	Pkg.
Medium Weight, Tape Tab Neck, Elastic Wrists		
NONLV200M	Medium, Yellow	100/cs
NONLV200	Regular/Large, Yellow	100/cs
NONLV200XL	X-Large, Yellow	100/cs
NONLV200XXL	XX-Large, Yellow	100/cs
NONLV200XXXL	XXX-Large, Yellow	100/cs
Medium Weight, Tape Tab Neck, Thumb Loop Wrists		
NONLV200TL	Regular/Large, Yellow	40/cs
NONLV200TLXL	X-Large, Yellow	40/cs
Medium Weight, Over-the-Head, Thumb Loop Wrists		
NONLV240	Regular/Large, Yellow	100/cs
NONLV240XL	X-Large, Yellow	100/cs
Lightweight, Over-the-Head, Thumb Loop Wrists		
NONLV242	Regular/Large, Yellow	100/cs
NONLV242XL	X-Large, Yellow	100/cs

AAMI Level 3 Isolation Gown

Recommended for: Med/Surg, ICU, Emergency/Trauma, Airborne, Contact and Droplet Precautions

- » Comfortable and breathable materials with classic tie waist
- » Choice of tape tab or hook and loop neck
- » Available with elastic wrists or knit cuffs

Item No.	Description	Pkg.
Heavyweight Multi-Ply, Tape Tab Neck, Elastic Wrists		
NONLV315	Regular/Large, Yellow	100/cs
NONLV315XL	X-Large, Yellow	100/cs
Premium Multi-Ply, Hook and Loop Neck, Knit Cuffs		
NONLV325	Regular/Large, Blue	50/cs
NONLV325XL	X-Large, Blue	50/cs
Microporous Breathable, Tape Tab Neck, Knit Cuffs		
NONLV350	Regular/Large, White	50/cs
NONLV350XL	X-Large, White	50/cs

Gown Features



Side ties eliminate the difficulties of tying behind the back
Available on AAMI Levels 1 and 2 gowns



Neck tie style
Available on AAMI Level 1 gowns



Thumb loop wrists
Available on AAMI Level 2 gowns



Hook and loop neck style
Available on AAMI Level 3 gowns

MEDLINE PPE COMPANION™

Protection That's Always By Your Side

Streamline caregiving. Move PPE closer to the point of care. This unique isolation cart has all you need to protect yourself—and your patients.

Isolation Cart

A Precaution Identification

Multilingual instructions show the clinically appropriate level of personal protection.

B Organized Storage

Multilingual labels on each drawer.

C Hand Protection

Three glove sizes. More storage behind.

D Face Protection

At-a-glance supply of face shields and earloop masks, ample storage for N95 respirators and more.

E Body Protection

Gowns dispense one at a time with a touch and a pull.

F Gown Storage

Additional storage for a minimum of 50 gowns.

G Quick Disinfection

Sanitize hands. Disinfect surfaces.



Know you're protected. Contact your Medline Representative for more information.







Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US    

Reference

1. Considerations for selecting protective clothing used in healthcare for protection against microorganisms in blood and body fluids. Centers for Disease Control and Prevention website. Table 3: ANSI/AAMI PB 70:12 Classification of barrier performance of surgical gowns, other protective apparel, surgical drapes and drape accessories. Available at: <http://www.cdc.gov/niosh/nppt/topics/protectiveclothing/>. Accessed May 3, 2016.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.

© 2016 Medline Industries, Inc. Companion is a trademark and Medline is a registered trademark of Medline Industries, Inc. MKT1678080 / LIT489R / 500 / SG / 5