

Product Information: Nipples

For more information, contact your Abbott Nutrition Representative or visit
www.abbottnutrition.com

Nipples

- 1 Provides the flexibility to choose a ready-for-use nipple and ring that match each infant's sucking strength or need.
- 1 **Infant** (latex-free)—for infants with developed sucking strength.
- 1 **Premature** (latex-free)—for preterm infants with less developed sucking strength.
- 1 **Orthodontic** (contains latex)—for newborn infants, to stimulate oral exercise while feeding.
- 1 **Slow Flow** (latex-free)-to support the underdeveloped "suck-swallow-breath" coordination of preterm infants. Provides a more controlled flow rate vs faster-flow nipples for increased physiological stability and a more efficient sucking pattern.¹
- 1 Easy to use.
- 1 Nipples and rings individually prepackaged and ready for use.
- 1 Easy, effective venting system.
- 1 Facilitates latching and consistent flow.
- 1 Disposable—not for reuse.



Safety Precautions

- 1 DO NOT use as a pacifier, as it presents a choking hazard.
- 1 DO NOT reuse, as reuse can cause deterioration or separation of the rubber, presenting a choking hazard.
- 1 Keep out of reach of children.

¹ Chang YJ, et al. *J Nurs Res* 2007;150(3):215-223.

Availability

List Number	Item
00079	Similac Infant Nipple & Ring / 250 ct
00094	Similac Premature Nipple & Ring / 250 ct
50512	Similac Orthodontic Nipple & Ring / 250 ct
53894	Similac Slow Flow Nipple and Ring / 250 ct

For more information, contact your Abbott Nutrition Representative or visit
www.abbottnutrition.com

Nutrition Information - Nipple Not Applicable

	Value

For more information, contact your Abbott Nutrition Representative or visit
www.abbottnutrition.com