

NIDEK *Medical*

**Proudly Presents
The Most Powerful
Member of our
Nuvo Family of
Oxygen Concentrators**



.....Mark 5
Nuvo 8

- Liter Flow -- 2 to 8 LPM
- Ideal for High Flow Patients
- Low Maintenance Costs
- Made in the USA

CE
0413



NIDEK Medical

Additional Features:

- Built on the Reliable Nuvo Platform
- Easy Access to filter, battery, and hourmeter
- Superior Grade of Molecular Sieve
- Spring Mounted High Efficiency Thomas Nexus Compressor

Specifications*

Nuvo 8 Oxygen Concentrator

**Nuvo 8 Oxygen Concentrator
with Oxygen Concentration Status Indicator (OCSI)**

115 Volt - 60 Hz

580

585

230 Volt - 50/60 Hz

980

985

Average Power Consumption
Oxygen Concentration (at 8 LPM)
Liter Flow
Outlet Pressure
Weight
Dimensions

500 Watts
90% (+5% / -3%)
2 to 8 LPM
17 PSIG (1.2 BAR)
53.5 lbs / 24.3 kg
15.5"W x 15.6"D x 27.8"H

490/625 Watts
90% (+5% / -3%)
2 to 8 LPM
17 PSIG (1.2 BAR)
55.5 lbs / 25.2 kg
39.4 cm x 39.6 cm x 70.6 cm

Accessories and Re-order Items

Compressor Intake Filter
External Cabinet Filter
Humidifier
Cannula with 7' tubing
Nuvo 8 User's Guide

9250-1180
9250-1025
9012-8774
9012-8780
2010-2238STD

9250-1180
9250-1025
9012-8774
9012-8780
2010-2238CE

* NOM - Specifications are subject to change without notice

2010-2238