

**MARK 5 Nuvo [M5C5] Serial No.** \_\_\_\_\_

**Date first used:** \_\_\_\_\_

\_\_\_\_\_

**Maintained by:** \_\_\_\_\_

\_\_\_\_\_

**Your distributor:** \_\_\_\_\_

**Address :** \_\_\_\_\_

\_\_\_\_\_

**Telephone :** \_\_\_\_\_

**PREVENTIVE MAINTENANCE:**

- a. Wash cabinet dust filter weekly.
- b. Inspect inlet air filter at each patient visit.  
Replace filter annually, or more often depending on environment.
- c. Check oxygen concentration every 5,000 hours or one year.

The manufacturer's instructions for the **preventive maintenance** of the devices are defined in the maintenance manual. Check with your service provider for any updates to the recommended schedules.

The work must be carried out by suitably trained technicians certified by the manufacturer.

**Use original spare parts only (see Pg. 7).**

Upon request, the supplier can provide circuit diagrams, spare parts lists, technical details or any other information of use to qualified technical personnel for parts of the device which are designated as being the manufacturer's responsibility or by the manufacturer as repairable.

**Medical Device Regulations require users and service providers to report to the manufacturer any incident that could, if repeated, result in injury to**