

Durability That Outlasts the Competition.

The Oxlife Independence by O2 Concepts is the only concentrator on the market with all light-weight metal construction and standard 5 year warranty. It was specifically designed to handle the rigors of 24/7 use. With continuous flow up to 3LPM and pulse settings up to 6 (96 mL), the Oxlife Independence is the ideal non-delivery solution for your most demanding oxygen patients.



- Patented ESA technology enables the Oxlife Independence to have the longest battery life in its class
- Can be used on all settings while plugged in a car (DC mode)
- Integrated DC power supply no external brick
- Durable exterior made of all light-weight metal to survive your most demanding oxygen patients

- FAA Approved
- Convenient access to preventative maintenance door for fast and easy pump filter change
- Integrated cart with 6" wheels offers "best in class" mobility and lighter overall weight
- Custom split battery configuration allows for "hot swapping" of battery. No loss of oxygen when switching battery







- Flow settings in .5 increments for better patient titration
- Large LCD screen with multi-color background (blue, red, yellow) — easy for visual and hearing impaired patients.

| SETTING | DOSAGE (ML) |
|---------|-------------|
| .5 | 8 |
| 1 | 16 |
| 1.5 | 24 |
| 2 | 32 |
| 2.5 | 40 |
| 3 | 48 |
| 3.5 | 56 |
| 4 | 64 |
| 4.5 | 72 |
| 5 | 80 |
| 5.5 | 88 |
| 6 | 96 |
| | |

Product Specifications*

Dimensions

Weight**

Oxygen Flow Rate

Battery Duration***

Sound Level @ 2 Pulse

Battery Recharge Time****

Oxygen Concentration

Oxygen Output Pressure

AC Power

DC Power

Operating Temperature

Altitude Operating Range

Transportation/Storage Temperature

Transportation/Storage Humidity

Warranty

Height: 20.2" • Width: 10.9" • Depth: 8"

Device: 16.7 lbs • Battery: 1.35 lbs/each

Pulse Dose Settings 0.5-6 in 0.5 increments (8mL-96mL)

0.5-3 LPM in Continuous Mode in 0.5 LMP increments

5.75 hours @ 2 pulse at 20 bpm • 2.5 hours @ 2LPM continuous

40 dBA

1.5 hours - 1 battery • 2.5 hours - 2 batteries

87%-95.6% at Sea Level at all settings

5psig +/- 1 psig (34 +/- 7 kPa)

Universal Input: 100-240 VAC • 50-60 Hz • 2 Amps

Output 24 Volts • 6.3 Amps

12 -15 VDC • 10 Amps

50°F (10°C) to 104°F (40°C)

0 - 13,123 Feet (0 - 4,000 Meters)

-4°F (-20°C) to 140°F (60°C)

0-95% Non-Condensing

5 Year Limited • 1 Year on Accessories/Battery



^{**}Production weight does not include cart kit



^{***}Battery duration shown is for extended battery system

^{****}Battery recharge time will increase if unit is running while charging

Ordering Information

| Part Description | Part Number | | |
|---|-------------|--|--|
| Systems | | | |
| Oxlife Independence POC Standard System The Oxlife Independence POC Standard System includes oxygen concentrator, 1 battery, pull handle and wheels, AC and DC power supplies, accessory bag and user manual. | 800-0001-1 | | |
| Oxlife Independence POC Extended Battery System The Oxlife Independence POC Extended Battery System includes oxygen concentrator, 2 batteries, pull handle and wheels, AC and DC power supplies, accessory bag and user manual. | 800-0001-2 | | |
| Accessories | | | |
| AC Power Cord/Power Supply | 800-1022 | | |
| DC Power Cord | 800-1001 | | |
| Battery - Lithium Ion Rechargeable | 800-1002S | | |
| Accessory Carry Bag | 800-1003 | | |
| Battery Charger - Desk Top Single Bay | 800-1004 | | |
| Battery Charger - Multiple (6 Bay) | 800-1005 | | |
| Inlet Air Filter (Pack of 6) | 800-1006 | | |
| Wheels (2) | 800-1007 | | |
| User Manual | 800-1008 | | |
| Cart Kit (Includes Wheels and Pull Handle) | 800-1010 | | |
| Pull Handle Kit | 800-1011 | | |
| Pump Filter Kit | 800-1013 | | |
| Humidifier Kit | 800-1015 | | |

*US Federal Law Restricts these devices to sale by or on the order of a physician



914 N. Broadway, Suite 100 Oklahoma City, OK 73102

T: 877-867-4008 F: 877-867-1314 www.o2-Concepts.com



