

one solution



OXLIFE INDEPENDENCE






THE PIECE THAT COMPLETES YOUR NEEDS.

The OxLife Independence by O2 Concepts is a portable oxygen concentrator designed with both patients and providers in mind. Boasting both continuous flow up to 3L/min and pulse settings up to 6 (96 mL), the OxLife Independence gives patients total freedom with 24/7 oxygen on the go while offering providers a one solution approach to maximize profits in today's turbulent reimbursement climate.

-  Patent-pending ESA technology enables the OxLife Independence to be the most efficient portable oxygen concentrator in the market

- Only POC in class that can be used on all settings while charging in car



- Integrated DC power supply -no external brick

- Custom split battery configuration gives patients the flexibility of lighter weight or extended battery life

- Integrated cart kit design offers lighter overall weight

- Easy access filter door -no tools necessary

- 6" wheels allow easier patient transport



- Longest battery life in class

- Medicare HCPCS codes E1390/E1392 (Pending)

- FAA Approved 

PULSE FLOW SETTINGS

SETTING	DOSAGE (ML)
.5	8
1	16
1.5	24
2	32
2.5	40
3	48
3.5	56
4	64
4.5	72
5	80
5.5	88
6	96



PULSE FLOW SETTINGS IN .5 INCREMENTS FOR BETTER PATIENT TITRATION

PRODUCT SPECIFICATIONS*

DIMENSIONS	Height: 20.2" . Width: 10.9" . Length: 8"
OXYGEN FLOW RATE	Pulse Dose Settings 1-6 in 0.5 Increments (8mL-96mL) 1-3 LPM in Continuous Mode in 0.5 LPM Increments
OXYGEN CONCENTRATION	90% +/- 3% at Sea Level
OXYGEN OUTPUT PRESSURE	5 psig +/- 1 psig (34 +/- 7 kPa)
AC POWER	Universal Input: 100-240 VAC . 50-60 Hz . 2.0 Amps Output: 24 Volts . 6.3 Amps
DC POWER	12 - 15 VDC . 10 Amps
OPERATING TEMPERATURE	50°F (10°C) to 104°F (40°C)
ALTITUDE OPERATING RANGE	0 - 13,123 Feet (0 - 4,000 Meters)
TRANSPORTATION/STORAGE TEMPERATURE	-4°F (-20°C) to 140°F (60°C)
TRANSPORTATION/STORAGE HUMIDITY	0 - 95% Non-Condensing
WARRANTY	3 Year Limited . 1 Year on Accessories/Battery

*O2 CONCEPTS RESERVES THE RIGHT TO CHANGE SPECIFICATIONS OF THIS PRODUCT AT ANYTIME WITHOUT NOTICE

ORDERING INFORMATION*

PART DESCRIPTION	PART NUMBER
SYSTEMS	
OXLIFE INDEPENDENCE POC SYSTEM (The Oxlife Independence POC System includes oxygen concentrator, 1 battery, pull handle and wheels, AC and DC power cords, accessory bag and user manual)	800-0001
ACCESSORIES	
ACCESSORY CARRY BAG	800-1003
BATTERY - LITHIUM ION RECHARGEABLE	800-1002
BATTERY CHARGER - DESK TOP SINGLE BAY	800-1004
BATTERY CHARGER - MULTIPLE (6 BAY)	800-1005
CART KIT (INCLUDES WHEELS AND PULL HANDLE)	800-1010
INLET AIR FILTER (PACK OF 6)	800-1006
AC POWER CORD/POWER SUPPLY	800-1000
DC POWER CORD	800-1001
PROTECTIVE COVER	800-1009
PULL HANDLE KIT	800-1011
USER MANUAL	800-1008
WHEELS (2)	800-1007

*US FEDERAL LAW RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN



14001 MCAULEY BLVD, SUITE 170 . OKLAHOMA CITY, OK 73134
 T: 877-867-4008 . F: 877-867-1314 . WWW.O2-CONCEPTS.COM



PROUDLY MADE IN THE USA

MKT-0001 Rev: A