

SPECIFICATIONS

DIMENSIONS	
Device Dimensions with Handle/Wheels	H: 20.29 in H: 51.5 cm W: 10.85 in W: 27.55 cm D: 9.45 in D: 24 cm
WEIGHTS	
Device Weight	16.7 lbs.
Cart Kit (Wheels & Pull Handle)	2.17 lbs.
Battery	1.4 lbs.
AC Charger	1.9 lbs.
DC Charger	.37 lbs.
MODES OF OPERATION	
Continuous	0.5 to 3 LPM in 0.5 LPM increments
Pulse	0.5 to 6.0 (8mL-96mL) setting increments
Battery Specification	14.4V Lithium Ion Battery
External Power Supply Power Input	AC power: 120 VAC, 50/60 Hz @ 2.5 amps DC power: 12 VDC
Oxygen Purity	87 - 95.6% in all settings
Operating Temperature	50°F (10°) to 104°F (40°C)
Operating Humidity	10% - 95% @ 82.4°F (28°C)
Operating Environment	Free of smoke, pollutants, and fumes.
Transport/ Storage Temperature	-4°F (-20°) to 140°F (60°C)
Transport/ Storage Humidity	0 -95% non-condensing
Altitude	0-13,123ft or (0-4000m)
Operating Time	24 hours per day when connected to an external AC or DC power source.
Cannula Specification	DO NOT use cannula tubing longer than 7ft (2.13m) when using Pulse Flow Mode Do NOT use cannula tubing longer than 50ft (15.25m) when using Continuous Flow Mode.