

Impregnated Gauze

*All Impregnated Gauzes are Sterile.



Medicated Petrolatum Gauze Non-Adherent Dressing

Non-adherent fine mesh gauze impregnated with 3% Bismuth Tribromophenate in a Petrolatum blend. Not made with natural rubber latex.

Cat. No.	Description
C-WMD18P	1 x 8 in.
C-WMD22	2 x 2 in.
C-WMD44	4 x 4 in.
C-WMD59	5 x 9 in.

Oil Emulsion Non-Adherent Dressing

Non-adherent cellulose acetate fabric impregnated with USP Mineral Oil. This dressing is used directly on the wound and allows exudate to pass through to a secondary dressing. Not made with natural rubber latex.

Cat. No.	Description
C-WNM33	3 x 3 in.
C-WNM38	3 x 8 in.
C-WNM383	3 x 8 in., 3 ea/pk
C-WNM316	3 x 16 in.
C-WNM59	5 x 9 in.

Petrolatum Gauze Non-Adherent Dressing

Non-adherent fine mesh gauze impregnated with USP White Petrolatum. Helps promote wound healing through a moist environment and are highly conformable. Not made with natural rubber latex.

Cat. No.	Description
C-WNM18	1 x 8 in.
C-WNM39	3 x 9 in.
C-WNM318	3 x 18 in.
C-WNM336	3 x 36 in.
C-WNM636	6 x 36 in.

Packing strips

We offer two types of packing strips: plain or impregnated with Iodoform. All are made with 100% cotton. Used for drainage of infected or open wounds. Not made with natural rubber latex.

Cat. No.	Description
Plain	
C-PG145P	1/4 in. x 5 yd
C-PG125P	1/2 in. x 5 yd
C-PG15P	1 in. x 5 yd
C-PG25P	2 in. x 5 yd

Iodoform

C-PG145I	1/4 in. x 5 yd
C-PG125I	1/2 in. x 5 yd
C-PG15I	1 in. x 5 yd
C-PG25I	2 in. x 5 yd