

A close-up photograph of two hands touching a blue and pink fabric sling. The hands are positioned to highlight the texture and quality of the material. The blue fabric has a fine, woven texture, and the pink fabric is a smooth, solid color. The hands are shown from a top-down perspective, with fingers gently touching the fabric.

**STRONG HOLD,
SOFT TOUCH.**

Safe Patient Handling Slings
Comfort and confidence for patients and staff

CHOOSE THE RIGHT SLINGS IN THE RIGHT FABRIC.

Your lifting equipment helps prevent injuries for both patients and caregivers. However, not all slings are created equal.

Our product portfolio includes eight innovative fabric choices with comfort-enhancing details to meet every patient need.



FREE

Schedule Your Safe Patient Handling Assessment

Know if your facility is compliant with Occupational Safety and Health Administration (OSHA) guidelines. We can help.

At no cost to you, we'll assess your entire facility and provide you with a concise report with everything you need to know to comply with OSHA guidelines.¹ For more information, or to schedule your assessment, **contact your Medline representative today.**

At-A-Glance Sling Selection

Decide which sling style is most appropriate for your intended use.

	Bathing and Toileting	Transferring	Repositioning
Polyester Mesh (pg. 8)	✓		
Toileting (pg. 13)	✓		
U-shaped (pg. 7, 8, 10, 11, 14, 15)	✓	✓	
Stand-assist (pg. 5, 12, 14)	✓	✓	
Knitted fabric (pg. 10, 11)		✓	
Spacer fabric (pg. 6, 7)		✓	
Single-patient use (pg. 14, 15, 17)		✓	
ComfortRise™ Sling (pg. 6)		✓	✓
Comfort Glide® (pg. 4)		✓	✓
Repositioning (pg. 6, 7, 17)			✓

Sizing Tool

- 1 Choose your preferred sling style.
- 2 Weigh the patient to determine the needed capacity.
- 3 Measure the patient for a proper fit using the chart below.

General Transfer Slings

Measure the patient's shoulder circumference (A):

Small	Medium	Large	X-Large
28"-45" (71-114 cm)	37"-57" (94-145 cm)	47"-67" (119-170 cm)	57"-80" (145-203 cm)

Toileting Slings

Measure the patient's chest circumference (B):


Small	Medium	Large	X-Large
26"-42" (66-107 cm)	30"-48" (76-122 cm)	38"-60" (97-152 cm)	58"-73" (147-185 cm)


Stand-Assist Slings

Measure the patient's chest circumference (B):

Small	Medium	Large	X-Large
28"-42" (71-107 cm)	30"-48" (76-122 cm)	36"-60" (91-152 cm)	50"-79" (127-201 cm)

Sizing tool is a guideline only.





Safe Patient Handling:

Hybrid Microfiber Nylon

Soft, breathable microfiber top,
smooth nylon bottom.

Comfort Glide® Sling

Microfiber/Nylon

A revolutionary innovation. All-in-one glide sheet and sling helps improve safety, efficiency and patient care. The skin-friendly, breathable material can remain in bed with the patient throughout their length of stay (single-patient use).


Glide sheet—the smooth underside minimizes friction for easy sliding, and the handles along the perimeter enable repositioning by multiple caregivers. Proven Comfort Glide Sheet technology is compatible with our Comfort Glide wedges and UltraSorbs® drypads.

Sling—the twelve strong straps offer attachment options to maximize patient comfort and safety.

- » Reduce staff lifting tasks during transfers
- » Save staff time with fewer steps to use a lift
- » Eliminate SKUs with multi-use products
- » Have more time to spend with patients

RECOMMENDED USE

Repositioning tasks and
patient transfers.



MDSC40018


USE WITH
6-PT.
CRADLES

Item No.	Style	Dimensions	Wt. Cap.	Pkg.
MDSC40018	Flat, Over Mattress	42" x 78" (106.7 x 198.1 cm)	1,000-lb. (453.6 kg)	12/cs

Luxurious and Resilient:

Textured Grip Fabric

Strong and flexible, soft to the touch, neoprene adds advanced functionality to a standard sling.



Grip Sling for Stand-Assist Lifts

Neoprene

The Grip Sling adheres to clothing and is less likely to slide up under the arms and create pressure during a stand-assist transfer.

- » Uniquely comfortable for patient satisfaction
- » Specially textured neoprene material
- » Soft padding in the underarm and lumbar areas

RECOMMENDED USE


Patient toileting and transfers with stand-assist lifts for extra support for standing.



MDSGRIPSLING3

USE WITH
STAND
ASSIST
LIFTS

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSGRIPSLING2	Medium	30" x 48" (76.2 x 121.9 cm)	500-lb. (226.8 kg)	ea
MDSGRIPSLING3	Large	36" x 60" (91.4 x 152.4 cm)	500-lb. (226.8 kg)	ea



Pressure Ulcer Prevention and Treatment:

Spacer Fabric

Warp-knitted polyester “spacer” fabric.

- » **Redistributes pressure.** Studies show that the spacer fabric improves surface-area contact.²
- » **Breathable.** Air and moisture permeable, and compatible with low air-loss mattresses.
- » **Reduces shear.** The fabric stretches in all directions to help reduce shear.
- » **Comfortable.** Studies show that patients prefer the soft feel over traditional sling fabrics, polyester and slipfit.²

ComfortRise™ Repositioning Sheet and Sling in One


Spacer Fabric

Dual functionality reduces patient disruptions, opportunities for staff back injuries and SKUs for facilities to purchase, store and put into service.

- » Smooth nylon underside glides easily over beds
- » Twelve straps for redistributing weight by multiple caregivers
- » Eliminates a glide sheet purchase and costs similar to a sling


RECOMMENDED USE

Repositioning and routine lifting, especially for patients with pressure ulcers or limited upper-body function and tone.



USE WITH
6-PT.
CRADLES

Item No.	Style	Dimensions	Wt. Cap.	Pkg.
CLIFT	Flat, Over Mattress	42" x 78" (106.7 x 198.1 cm)	700-lb. (317.5 kg)	ea




U-Shaped Sling

Spacer Fabric

RECOMMENDED USE

Toileting and routine lifting, especially for patients with pressure ulcers or limited upper-body function and tone.

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSPUNVS	Small	33.5" x 41" (85.1 x 104.1 cm)	450-lb. (204 kg)	ea
MDSPUNVM	Medium	40" x 48" (101.6 x 121.9 cm)	450-lb. (204 kg)	ea
MDSPUNVL	Large	44" x 53" (111.8 x 134.6 cm)	450-lb. (204 kg)	ea
MDSPUNVXL	X-Large	47" x 60" (119.4 x 152.4 cm)	450-lb. (204 kg)	ea



USE WITH
6-PT.
CRADLES


Hammock Sling

Spacer Fabric

RECOMMENDED USE

Toileting and routine lifting, especially for patients with pressure ulcers or limited upper-body function and tone.

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSPHAMS	Small	33.5" x 37.5" (85.1 x 95.3 cm)	450-lb. (204 kg)	ea
MDSPHAMM	Medium	38" x 42" (96.5 x 106.7 cm)	450-lb. (204 kg)	ea
MDSPHAML	Large	41" x 43" (104.1 x 109.2 cm)	450-lb. (204 kg)	ea
MDSPHAMXL	X-Large	46.5" x 46.5" (118.1 x 118.1 cm)	450-lb. (204 kg)	ea



USE WITH
6-PT.
CRADLES


Repositioning Sling

Spacer Fabric

RECOMMENDED USE

Turning or repositioning patients in bed, especially for those with pressure ulcers.

Item No.	Style	Dimensions	Wt. Cap.	Pkg.
MDSPPPOS	Flat	42" x 78" (106.7 x 198.1 cm)	440-lb. (200 kg)	ea
MDSPPPOSE	Fitted, with Elastic	42" x 78" (106.7 x 198.1 cm)	440-lb. (200 kg)	ea




USE WITH
6-PT.
CRADLES



Soft, Sturdy and Colorfast:

Polyester Mesh

Our polyester mesh is ultra-strong and soft to the touch. The synthetic material retains its brilliant color, wash after wash. Skillful workmanship and reinforced stitching support high weight capacities.



USE WITH
6-PT.
CRADLES

MDSMMHS1

U-Shaped Sling


Polyester Mesh Fabric

RECOMMENDED USE

Patient bathing and toileting.

With Head Support

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMMHS1	Small	33" x 49" (83.8 x 124.5 cm)	450-lb. (204 kg)	ea
MDSMMHS2	Medium	38.5" x 54" (97.8 x 137.2 cm)	450-lb. (204 kg)	ea
MDSMMHS3	Large	42.5" x 59" (108 x 149.9 cm)	450-lb. (204 kg)	ea




USE WITH
6-PT.
CRADLES


MDSMM2

Without Head Support


Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMM2	Medium	39" x 47" (99 x 119.4 cm)	450-lb. (204 kg)	ea
MDSMM3	Large	43" x 48" (109.2 x 121.9 cm)	450-lb. (204 kg)	ea



MDSMR110



MDSMR114



MDSPM111



Full-Body Sling

Polyester Mesh Fabric

RECOMMENDED USE

Routine lifting, especially for amputee patients.

Without Commode Opening

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMR110	Medium	30" x 47" (76.2 x 119.4 cm)	450-lb. (204 kg)	ea
MDSMR111	Large	34" x 52" (86.4 x 132.1 cm)	450-lb. (204 kg)	ea
MDSMR140	XX-Large	42" x 52" (106.7 x 132.1 cm)	700-lb. (317.5 kg)	ea

With Commode Opening

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMR114	Medium	32" x 48" (81.3 x 48")	450-lb. (204 kg)	ea
MDSMR115	Large	35" x 50" (88.9 x 127 cm)	450-lb. (204 kg)	ea
MDSMR116	X-Large	37" x 56" (94 x 142.2 cm)	450-lb. (204 kg)	ea
MDSMR141	XX-Large	42" x 50" (106.7 x 127 cm)	700-lb. (317.5 kg)	ea


Pediatric Sling

Polyester Mesh Fabric

RECOMMENDED USE

Bathing and toileting of pediatric patients over 35-lb. (16 kg).


Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSPM111	One Size Fits Most	21" x 33" (53.3 x 83.8 cm)	100-lb. (45.4 kg)	ea



Strong and Comfortable:

Knitted Fabric

This fabric has just the right amount of “give” for a secure, comfortable and easy fit. But don’t let it fool you. This heavy-duty choice is designed to withstand high weight capacities and rigorous healthcare use.



USE WITH
6-PT.
CRADLES

MDSMPHS1

U-Shaped Padded Sling


Knitted Fabric

RECOMMENDED USE

Patient transfers.

With Head Support

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMPHS1	Small	31.5" x 50" (80 x 127 cm)	450-lb. (204 kg)	ea
MDSMPHS2	Medium	38" x 54" (96.5 x 137.2 cm)	450-lb. (204 kg)	ea
MDSMPHS3	Large	44" x 64" (111.8 x 162.6 cm)	450-lb. (204 kg)	ea
MDSMPHS4	X-Large	48" x 73" (121.9 x 185.4 cm)	700-lb. (317.5 kg)	ea




USE WITH
6-PT.
CRADLES


MDSMP2

Without Head Support

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMP2	Medium	38" x 50" (96.5 x 127 cm)	450-lb. (204 kg)	ea
MDSMP3	Large	44.5" x 51" (113 x 129.5 cm)	450-lb. (204 kg)	ea
MDSMP4	X-Large	46" x 59" (116.8 x 149.9 cm)	700-lb. (317.5 kg)	ea



USE WITH
4-PT.
CRADLES



USE WITH
4-PT.
CRADLES

U-Shaped Padded Sling

Knitted Fabric

RECOMMENDED USE


Patient transfers.

With Head Support

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMHC70012	Medium	37" x 47" (94 x 119.4 cm)	500-lb. (204 kg)	ea
MDSMHC70011	Large	41" x 57" (104.1 x 144.8 cm)	500-lb. (226.8 kg)	ea
MDSMHC70010	X-Large	46" x 60" (116.8 x 152.4 cm)	700-lb. (317.5 kg)	ea

Without Head Support

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMHC70002	Medium	34" x 36" (86.4 x 91.4 cm)	500-lb. (226.8 kg)	ea
MDSMHC70001	Large	41" x 57" (104.1 x 144.8 cm)	500-lb. (226.8 kg)	ea
MDSMHC70000	X-Large	46" x 46" (116.8 x 116.8 cm)	700-lb. (317.5 kg)	ea



USE WITH
4-PT.
AND
6-PT.
CRADLES

Pediatric Padded Sling

Knitted Fabric

RECOMMENDED USE

Transfers of pediatric patients over 35-lb. (16 kg).

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSPP121	One Size Fits Most	21" x 32" (53.3 x 81.3 cm)	100-lb. (45.4 kg)	ea



USE WITH
STAND
ASSIST
LIFTS


Stand-Assist Padded Sling

Knitted Fabric/Flannel

RECOMMENDED USE

Patient toileting and transfers with stand-assist lifts for extra support for standing.

Item No.	Size	Chest Circ.	Wt. Cap.	Pkg.
MDSMSA1	Small	28"-42" (71-107 cm)	500-lb. (226.8 kg)	ea
MDSMSA2	Medium	30"-48" (76-122 cm)	500-lb. (226.8 kg)	ea
MDSMSA3	Large	36"-60" (91-152 cm)	500-lb. (226.8 kg)	ea
MDSMSA4	X-Large	50"-79" (127-201 cm)	600-lb. (272.2 kg)	ea



USE WITH
STAND
ASSIST
LIFTS

Stand-Assist Buttock Strap

Knitted Fabric

RECOMMENDED USE



An accessory for extra support with padded stand-assist slings.

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMBSS	Standard	34" x 7" (86.4 x 17.8 cm)	400-lb. (181 kg)	ea
MDSMBSB	X-Large	58" x 11" (147.3 x 28cm)	700-lb. (317.5 kg)	ea

Luxurious and Resilient:


Padded Flannel

Soft flannel provides a secure and easy fit as comfortable as your favorite shirt. Yet strong enough to support patients up to 450-lb.



USE WITH
4-PT.
AND
6-PT.
CRADLES

MDSMR112



USE WITH
4-PT.
CRADLES

USE WITH
STAND
ASSIST
LIFTS

MDSMT02

Woven Padded Full-Body Sling

Flannel/Cordura

RECOMMENDED USE

Routine lifting, especially for amputee patients.

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMR112	Medium	30" x 47" (76.2 x 119.4 cm)	450-lb. (204 kg)	ea
MDSMR113	Large	34" x 52" (86.4 x 132.1 cm)	450-lb. (204 kg)	ea
MDSMR117	X-Large	36" x 57" (91.4 x 144.8 cm)	450-lb. (204 kg)	ea

Padded Toilet Sling

Flannel/Knitted Fabric

RECOMMENDED USE

Patient toileting. Can also be used with stand-assist lifts.


Item No.	Size	Chest Circ.	Wt. Cap.	Pkg.
MDSMT01	Small	26"-42" (66-107 cm)	450-lb. (204 kg)	ea
MDSMT02	Medium	30"-48" (76-122 cm)	450-lb. (204 kg)	ea
MDSMT03	Large	38"-60" (97-152 cm)	450-lb. (204 kg)	ea
MDSMT04	X-Large	58"-73" (147-185 cm)	450-lb. (204 kg)	ea



Single Patient Use:


Polypropylene

Cost-effective polypropylene slings are designed for single-patient use. The right choice for facilities concerned about slings not returning from the laundry.




USE WITH
6" PT.
CRADLES

MDSMDHS1



USE WITH
6" PT.
CRADLES

MDSMD3



USE WITH
STAND
ASSIST
LIFTS

MDSMDSA663

U-Shaped Sling, Single-Patient Use

Polypropylene

RECOMMENDED USE

Patient bathing and toileting.

With Head Support

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMDHS1	Small	33" x 49" (83.8 x 124.5 cm)	700-lb. (317.5 kg)	10/cs
MDSMDHS2	Medium	38.5" x 54" (97.8 x 137.2 cm)	700-lb. (317.5 kg)	10/cs
MDSMDHS3	Large	42.5" x 59" (108 x 149.9 cm)	700-lb. (317.5 kg)	10/cs
MDSMDHS4	X-Large	47" x 73" (119.4 x 185.4 cm)	700-lb. (317.5 kg)	10/cs

Without Head Support

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSMD2	Medium	39" x 47" (99 x 119.4 cm)	700-lb. (317.5 kg)	10/cs
MDSMD3	Large	43" x 48" (109.2 x 121.9 cm)	700-lb. (317.5 kg)	10/cs

Stand-Assist Sling, Single-Patient Use

Polypropylene

RECOMMENDED USE

Patient toileting and transfers.


Item No.	Size	Chest Circ.	Wt. Cap.	Pkg.
MDSMDSA662	Medium	30"-48" (76-122 cm)	500-lb. (226.8 kg)	10/cs
MDSMDSA663	Large	36"-60" (91-152 cm)	500-lb. (226.8 kg)	10/cs

Limb Lift Sling, Single-Patient Use

Polypropylene

RECOMMENDED USE

For lifting to provide care, or for ROM tasks.



MDSLIMBLIFT

USE WITH
2-PT.
CRADLES

USE WITH
4-PT.
AND
6-PT.
CRADLES

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSLIMBLIFT	One Size Fits Most	40" x 6.5" (101.6 x 16.5 cm)	450-lb. (204 kg)	10/cs


Low-Back Sling, Single-Patient Use

Polypropylene

RECOMMENDED USE

More active patients with lower acuity.

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSGD432	Medium	38" x 39" (96.5 x 99 cm)	600-lb. (272.2 kg)	10/cs
MDSGD433	Large	40" x 44" (101.6 x 111.8 cm)	600-lb. (272.2 kg)	10/cs
MDSGD434	X-Large	41" x 47" (104.1 x 119.4 cm)	600-lb. (272.2 kg)	10/cs



MDSGD432

USE WITH
2-PT.
CRADLES


High-Back Sling, Single-Patient Use

Polypropylene

RECOMMENDED USE

More active patients with lower acuity.

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSGDH442	Medium	38" x 48" (96.5 x 121.9 cm)	600-lb. (272.2 kg)	10/cs
MDSGDH443	Large	39.5" x 54" (100.3 x 137.2 cm)	600-lb. (272.2 kg)	10/cs
MDSGDH444	X-Large	42" x 56" (106.7 x 142.2 cm)	600-lb. (272.2 kg)	10/cs



MDSGD442

USE WITH
2-PT.
CRADLES


U-Shaped Clip Sling, Single-Patient Use

Polypropylene

RECOMMENDED USE

Patient bathing and toileting.

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
BESCLIPSLNGD2	Medium	35" x 47" (88.9 x 119.4 cm)	600-lb. (272.2 kg)	10/cs
BESCLIPSLNGD3	Large	38" x 51" (96.5 x 129.5 cm)	600-lb. (272.2 kg)	10/cs
BESCLIPSLNGD4	X-Large	42" x 55" (106.7 x 139.7 cm)	600-lb. (272.2 kg)	10/cs



BESCLIPSLNGD2

USE WITH
CLIP
CRADLES


Durable and Fade-Resistant:

Synthetic Polyester

Support high weight capacities with polyester slings. Retains its color longer than traditional fabrics.

U-Shaped Clip Sling

Polyester Fabric




USE WITH
CLIP
CRADLES

RECOMMENDED USE

Patient bathing and toileting with clip-style positioning cradles.

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
BESCLIPSLING2	Medium	35" x 47" (88.9 x 119.4 cm)	600-lb. (272.2 kg)	ea
BESCLIPSLING3	Large	38" x 51" (96.5 x 129.5 cm)	600-lb. (272.2 kg)	ea
BESCLIPSLING4	X-Large	42" x 55" (106.7 x 139.7 cm)	600-lb. (272.2 kg)	ea



USE WITH
2-PT.
CRADLES


Low-Back Sling

Polyester Fabric

RECOMMENDED USE

With more active patients with lower acuity.

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSGP472	Medium	38" x 39" (96.5 x 99 cm)	600-lb. (272.2 kg)	ea
MDSGP473	Large	40" x 44" (101.6 x 111.8 cm)	600-lb. (272.2 kg)	ea
MDSGP474	X-Large	41" x 47" (104.1 x 119.4 cm)	600-lb. (272.2 kg)	ea



USE WITH
2-PT.
CRADLES

High-Back Sling

Polyester Fabric

RECOMMENDED USE

With more active patients with lower acuity.

Item No.	Size	Dimensions	Wt. Cap.	Pkg.
MDSGPH482	Medium	38" x 48" (96.5 x 121.9 cm)	600-lb. (272.2 kg)	ea
MDSGPH483	Large	39.5" x 54" (100.3 x 137.2 cm)	600-lb. (272.2 kg)	ea
MDSGPH484	X-Large	42" x 56" (106.7 x 142.2 cm)	600-lb. (272.2 kg)	ea


Gait Training Sling

Polyester Fabric

RECOMMENDED USE

Therapeutic treatments for regaining the ability to walk.

Item No.	Chest Circ.	Wt. Cap.	Pkg.
MDSMGAIT	40"-54" (101.6-137.1 cm)	300-lb. (136 kg)	ea



MDSMGAIT

USE WITH
2-PT.
CRADLES


Single-Patient Use Repositioning Sling

Polyester/Cotton Fabric

RECOMMENDED USE

Turning or repositioning patients in bed.

Item No.	Dimensions	Wt. Cap.	Pkg.
MDSTRIDR011	42" x 78" (106.7 x 198.7 cm)	1,100-lb. (499 kg)	10/cs



MDSTRIDR011

USE WITH
4-PT.
AND
6-PT.
CRADLES


Reusable Repositioning Sling

Polyester/Cotton Fabric

RECOMMENDED USE

Turning or repositioning patients in bed.

Item No.	Dimensions	Wt. Cap.	Pkg.
MDSTRIRR01	42" x 78" (106.7 x 198.7 cm)	1,100-lb. (499 kg)	ea



MDSTRIRR01

USE WITH
4-PT.
AND
6-PT.
CRADLES

Sling Maintenance Best Practices

- » Check condition before each use. If there is any fraying or visible wear and tear, do not use.
- » Single-patient use slings should be thrown away when soiled or upon patient discharge.
- » Reusable slings should be replaced every six months.
- » Follow care instructions on wash tag. If illegible, do not use.
- » Keep at least two reusable slings per patient on hand—one available and one in the laundry.

For more solutions for safe patient handling, contact your Medline Representative today.



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US    

References: 1. Inspection Guidance for Inpatient Healthcare Settings. Occupational Safety and Health Organization (OSHA). Available at: www.osha.gov/dep/enforcement/inpatient_insp_06252015.html. Accessed May 3, 2017. 2. Mellson J, Richardson N. The Impact of Hoist Sling Fabrics on Gluteal Interface Pressure While Sitting in Healthy Individuals: A Controlled Pre-Post Test Study. American Journal of SPHM (Safe Patient Handling and Movement). 2012;2(3):79-85. Available at: www.americanjournalofsphm.com/shop_ajsphm/index.php?route=product/product&product_id=88. Accessed May 3, 2017.

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure. © 2017 Medline Industries, Inc. All rights reserved. ComfortRise is a trademark, and Medline, Comfort Glide and UltraSorbs are registered trademarks of Medline Industries, Inc. MKT1788408 / LIT486 / 43