

## In-Room™ Sharps Containers with Always-Open Lid



851201

### In-Room™ System with Always-Open Lid

- Always-Open lids provide convenient, quick and easy point-of-use sharps disposal.
- Accommodates larger syringes.
- Economical horizontal drop maximizes container volume.
- Containers easily lock for final disposal.
- 8.5" wide opening will accept up to 60mL syringe.
- Wall enclosures, brackets and adhesive-backed holders are available.



85221

Code	Description	Dimensions	Ship Case
851201	5 Quart, Clear	11"H x 4.75"D x 10.75"W	20
851301	5 Quart, Transparent Red	11"H x 4.75"D x 10.75"W	20
85321	2 Gallon, Clear	11.5"H x 6"D x 13.75"W	10
85321R	2 Gallon, Transparent Red	11.5"H x 6"D x 13.75"W	10
85221	3 Gallon, Clear	16.25"H x 6"D x 13.75"W	10
85221R	3 Gallon, Transparent Red	16.25"H x 6"D x 13.75"W	10

## In-Room™ Sharps Containers with Mailbox-Style Lid



85021

### In-Room™ System with Mailbox-Style Lid

- Mailbox-style lid can be kept closed when container is not in use.
- Economical horizontal drop maximizes container volume.
- Containers lock for final disposal.
- Wall enclosures, brackets and adhesive-backed holders are available.



85121

Code	Description	Dimensions	Ship Case
85021	2 Quart, Clear	6.25"H x 4.75"D x 10.75"W	20
85031	2 Quart, Transparent Red	6.25"H x 4.75"D x 10.75"W	20
85121	5 Quart, Clear	11"H x 4.75"D x 10.75"W	20
85131	5 Quart, Transparent Red	11"H x 4.75"D x 10.75"W	20