

INVACARE®

PLATINUM™ 10

OXYGEN CONCENTRATOR



The **Invacare® Platinum™ 10 Oxygen Concentrator** designed to provide an alternative to liquid systems and is known for reliability and performance. The 10 liter unit features quiet operation, a top handle and bottom recesses for easy handling and a side filter access door. In addition, the concentrator is compatible with the Invacare® HomeFill® Oxygen System and is available with an **Invacare® SensO₂® Oxygen Purity Sensor**.

[Click here to view detailed product information and pricing.](#)

Call 800-397-5899 to order today!

Features

Pressure-based system is superior to timing-based systems

Self diagnostic electronics designed to minimize troubleshooting time

Top handle for ease of mobility

Battery-free power loss alarm, Invacare® Sieve-GARD™ Monitor protects sieve from degradation, compressor 35 psi pressure-relief valve and low flow alarm

Provider Benefits

Designed to provide continuous flow up to 10 liters per minute as an alternative to liquid systems

All Platinum 10 stationary concentrators include a HomeFill system port for inventory flexibility

Patient-accessible reset button helps minimize emergency service calls



Patient Benefits

User friendly design with easy to see control panel, humidifier bottle holder and filter access

HomeFill system compatibility allows patients with flow rates up to 5LPM to enjoy life's experiences with small, refillable cylinders

Reliable concentrator provides unlimited, worry-free oxygen supply

HCPCS Code: E1390



1. Invacare® HomeFill® Oxygen System connection
2. Easy to maintain filter
3. Also available with Invacare® SensO₂® Oxygen Purity Sensor



Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2017

Invacare Corporation.
All rights reserved.
Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification.
Form No. 11-190
Rev. 06/17
170582



Model #	Dimensions	Weight	Shipping Weight	Quantity	Oxygen Flow	Flow Range	Maximum Outlet Pressure	Pressure Relief Mechanism
IRC10LX	26 3/8" H	54 lb.	60 lb.	1 Per Case	10 LPM	2-10 LPM	9+/- 0.5 PSI	35 PSI
IRC10LX02	18 3/8" W 14 3/8" D			18 Per Pallet				+/- 3.5 PSI
	O ₂	Filter	Electrical Requirements	Rated Current Input	Average Power Consumption	Sound Level	Limited Warranty	
	94%-87% At All Flow Rates	Cabinet Compressor Inlet, Outlet HEPA	120 VAC +/- 10%, 60 Hz (132VAC/108VAC)	5.2 A	550 W	58 db (A)	3 Years	

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.