



# Aqua Creek Products

The Leaders of Recreation, Fitness and Ability

Pool Access Chairs are an ideal solution for swimming pools with ramps or zero depth entries. These chairs can be used in showers, locker rooms and pool environments. The Pool Access Chairs feature an anti-tip design with non-rusting parts constructed from medical grade PVC. Swing out armrests and retractable footrest make lateral transfers easy for any individual.

## Features:

- 250 lb (113 kg) weight capacity - Pool Access Chair
- 400 lb (181 kg) weight capacity - Heavy Duty Chair
- Durable, non-corrosive PVC construction
- Stainless steel hardware
- Plastic rear wheels with rubber tread
- Anti-tip design
- Swing out armrests for lateral transfers
- One year warranty

## Accessories:

- Chest strap - *Item #F-426CS*
- Pool Access Chair with optional mesh seat  
*Item #F-011PAC-MS*
- Seat belt - *Item #F-11PACSB*



**Pool Access Chair**  
*Item #F-011PAC*



**Heavy Duty Pool Access Chair**  
*Item #F-400PWC*

**Now available in knockdown version**  
*Item #F-400PWC-KD*

**PROUDLY SOLD BY:**