

Compartment Style Compressor Nebulizer

HCPCS: E0570-1, E0570, A7003, A7013



FIGURE A



FIGURE B



FIGURE C

- Includes disposable neb kit and tubing
- Attractive, Sleek Design
- Durable, Quiet System
- Piston-Driven Compressor
- Convenient Storage Compartment

PRODUCTS

Description

MQ6000	Compartment Style Compressor Nebulizer, EA
--------	--

SPECIFICATIONS

Description

Liter Flow	, 8 lpm
Maximum Pressure	, 30 psi
Nebulizer Medication Capacity	5ml
Noise Level	55 dBA
Operating Humidity Range	10% to 95% RH
Operating Pressure	, 10 psi
Operating Temperature Range	10° to 40°C (50° to 104°F)
Particle Size	0.5 µm to 5.0 µm
Power Consumption	80 VA
Power Requirments	120V, 50-60 Hz
Storage Humidity Range	10% to 95% RH
Storage Temperature Range	-20° to 70°C (-4° to 158° F)
Warranty	Limited 5 Year
Weight	4.3 lb

ACCESSORIES

Description

18254	Reusable Neb Kit 10/case, CS
MASK 001A	Adult Aerosol mask CS/50, CS
NEB KIT 500	Neb Kit, Inc Mouth pc&Tub50/cs, CS