



BEDSIDE CARE ESSENTIALS



**A Comprehensive Range of Plastic
and Stainless Steel Medical Utensils**

GET REAL SAVINGS EVERY DAY. CHOOSE A PROVIDER WITH END-TO-END EXPERTISE.

You're under pressure to cut costs. But not at the expense of quality. We understand. And we've got your back. Because helping you win on quality and price is what Medline—the top provider of bedside essentials in the world—does best.

We answer the call as your single supply source with bedside essentials available both individually and in pre-assembled kits. But we don't stop there. As the only provider who controls the entire supply chain, we're continually innovating operational performance to match your need for ongoing value.

Manufacturer-Direct Pricing

We leverage economies of scale. Eliminate third-party expenses. And pass the savings on to you.

- » Largest, privately held manufacturer and distributor of healthcare and surgical supplies in the U.S.
- » Growing network of ISO 13485-certified and FDA QSR-compliant manufacturing facilities

Increased Productivity

Partnering with Medline as your single supply source is not only a smart savings strategy, it's also a big time saver.

- » Simplified ordering and invoicing process
- » No need to track and manage multiple vendors

Quality Assurance

Rigorous quality control and testing mean our patient utensils are efficient and tough enough to resist cracking and leaking.

- » Engineered to stack, store and ship as economically as possible
- » Six quality control and test labs certified to ISO 17025
- » 800+ highly trained quality professionals worldwide

Speedy Delivery

To minimize shipping lead times and reduce freight costs, our warehouses are located close to our customers.


- » 43 distribution centers in North America
- » 425+ trucks nationwide
- » 1-day ship to 96% of the U.S.

Hands-On Help

Experience the proactive support of a dedicated Medline Representative and a team of experts focused on your success.


- » Customer service on call 24/7 for overnight shipping and emergency orders
- » Continuously exceeding a 99% service level
- » Easy ordering by phone, fax, internet or through your Representative

greensmart.

For pigment-free products, simply look for the greensmart™ leaf. 

Talk to a cost savings expert. Contact your Medline Representative for more product information and samples today.



	Bedpans	4-7
	Bedside Containment	8-11
	Stainless Steel Components	12-14
	Autoclavable Plastic Products	15
	Bedside Basins	16-17
	Bedside Hydration	18-21
	Bedpan and Urinal Holders	22-23

Bedpans

Our durable bedpans come in a variety of shapes and sizes to aid pediatric through bariatric patients. They have a non-stick surface for optimal comfort and are not made with natural rubber latex.





Stack-A-Pan Bedpans

These bedpans nest together, saving storage space.

- » Fits over standard toilets when seat is raised
- » Graduated in ounces and milliliters to 40-oz./1200 ml
- » Measured total volume of 55-oz./1600 ml
- » 14 L x 11 W x 2.75" H (35.6 x 27.9 x 7 cm)


Item No.	Description	Pkg.
DYND80245	Graphite	50/cs
DYND80216	Gold	50/cs
DYNEC80245C 	Pigment Free	50/cs



Pontoon Bedpans

This stackable design maximizes storage space.

- » Saddle shaped with a high-rolled front
- » Graduated in ounces and milliliters to 40-oz./1200 ml
- » Measured total volume of 55-oz./1600 ml
- » Tested to support up to 350-lbs.
- » 14 L x 11 W x 2.75" H (35.6 x 27.9 x 7 cm)


Item No.	Description	Pkg.
DYND80217	Graphite	20/cs
DYND80211	Gold	20/cs
DYNEC80217C 	Pigment Free	20/cs



Contour Bedpans

Lightweight and saddle-shaped, these bedpans have a coccyx groove for extra comfort.

- » A high, rolled front and undercut rear helps prevent spills
- » Grooved bottom separates liquid from solid
- » Graduated in ounces and milliliters to 40-oz./1200 ml
- » Measured total volume of 85-oz./2500 ml
- » Tested to support up to 400-lbs.
- » 13 L x 11 W x 4.25" H (33 x 27.9 x 10.8 cm)

Item No.	Description	Pkg.
DYND80281	Graphite	8/cs
DYND80282	Gold	8/cs
DYNEC80281C 	Pigment Free	8/cs




Durafit® Contour Bedpans

These bedpans are extra large and have a 450-lb. weight capacity.

- » Measured volume of 65-oz./1900 ml
- » Grooved bottom helps keep liquid and solids separate
- » 15 L x 12 W x 3.5" H (38.1 x 30.5 x 8.9 cm)

Item No.	Description	Pkg.
DYND80290	Graphite	20/cs



Fracture Bedpans

Ideal for immobile patients, these bedpans have a tapered front for easy placement.

- » Loop handle for trouble-free maneuvering
- » Measured total volume of 26-oz./770 ml
- » Tested to support up to 250-lbs.
- » Suitable as a female urinal
- » 12 L x 9 W x 1" H (30.5 x 22.9 x 2.54 cm)


Item No.	Description	Pkg.
DYNC8522	Graphite	24/cs
DYNC8521	Gold	24/cs
DYNEC8522C	Pigment Free	24/cs

Deluxe Fracture Bedpans

These bedpans have a larger weight capacity and lower tapered sides than the standard version.

- » Loop handle for easy maneuvering
- » Measured total volume of 52-oz./1500 ml
- » Tested to support up to 275-lbs.
- » Suitable as a pediatric or female urinal
- » 10.5 L x 10 W x 1.88" H (26.7 x 25.4 x 4.8 cm)

Item No.	Description	Pkg.
DYND80230A	Graphite	24/cs
DYNEC80230C	Pigment Free	24/cs



Large Fracture Bedpans

These bedpans bridge the size gap between our Deluxe and Mega Fracture options.

- » Loop handle facilitates maneuvering
- » Measured volume of 75-oz./2200 ml
- » 16 L x 10.5 W x 3" H (40.6 x 26.6 x 7.6 cm)

Item No.	Description	Pkg.
DYNC8555	Graphite	12/cs
DYNC8551	Gold	12/cs

Mega Fracture Bedpans

These extra large bedpans support up to 900-lbs.

- » Tapered front slides easily under patient
- » Loop handle facilitates maneuvering
- » Measured volume of 80-oz./2400 ml
- » Suitable as a bariatric female urinal
- » 17.75 L x 16.25 W x 5" H (45.2 x 41.4 x 12.7 cm)

Item No.	Description	Pkg.
DYNC8552	Graphite	4/cs
DYNEC8552C 	Pigment Free	4/cs

Bedside Containment

Make bedside care easier with our wide assortment of plastic containers.



Urinals

Our urinals have a patented splash-proof lid and store on their side, preventing accidental tipping and spilling.

- » Graduated in ounces and milliliters
- » Measured volume of 32-oz./1000 ml
- » Available with a glow-in-the-dark lid

Item No.	Description	Pkg.
DYND80235S	Supreme Male Urinal	48/cs
DYND80235SD	Glow and Let Flow, Deluxe Male Urinal, Glow-in-the-Dark Cover	48/cs
DYND80234	Male Urinal, Mini-Pack, Individually Shrink-Wrapped	12/cs
DYND80235R	Male Urinal, Retail Packaged	6/cs
DYND88235	Autoclavable Urinal w/Blue Lid	12/cs



DYND80235S



DYND80235SD

24-Hour Urine Collection Bottle

Designed for overnight collection.

- » Graduated in ounces and milliliters
- » Measured volume of 96-oz./3000 ml
- » Wide, easy-grip handle
- » Can be used with calculi strainer, sold separately

Item No.	Description	Pkg.
DYND80024	24-Hour Urine Container, Amber	20/cs
DYND4712	Calculi Strainer	100/cs



Specimen Pans

Provides efficient collection of stool and urine samples.

- » Graduated in ounces and milliliters
- » Reinforced ribbing for durability
- » Fits any standard size commode
- » Available in one-piece and two-piece options

Item No.	Description	Pkg.
DYND36600	Standard Specimen Pan, 30-oz./900 ml, Translucent w/Black Graduations	100/cs
DYND36608	Standard Specimen Pan, 30-oz./900 ml, White w/Black Graduations	100/cs
DYND36605	Deluxe Specimen Pan, 40-oz./1200 ml, White w/Black Graduations	100/cs
DYND36500	Two-Piece Specimen Collection Container, 40-oz./1200 ml, Light Blue	100/cs




DYND36605



DYND36600



DYND80416



DYND80419



DYND80417

Triangular Intake/Output Containers

Choose from black or molded graduations in both ounces and cubic centimeters for accurate measurements.

- » Measured volume of 32-oz./1000cc
- » Flexible and rigid options

Item No.	Description	Pkg.
DYND80417	Flexible, Translucent w/Molded Graduations	200/cs
DYND80419	Flexible, Translucent w/Black Graduations	200/cs
DYND80416	Rigid, Clear w/Molded Graduations	200/cs



Denture Cups/Containers

These containers are designed with larger lips so cups separate with minimal effort and sticking.

- » Clear lids make contents easy to identify
- » Lid imprinted for patient-specific labeling: Name, room, date and contents

Item No.	Description	Pkg.
DYND70292	Clear Lid, 8-oz., Translucent	250/cs
DYND70293	Clear Lid, 8-oz., Aqua	250/cs
DYND70295	Clear Lid, 8-oz., Orange	250/cs
DYND70294B	Hinged Lid, 8-oz., Hot Pink	200/cs



Bulk Bowls

Our durable plastic bowls are available in a variety of sizes with molded graduations.

- » Light blue with graduations in ounces and cubic centimeters
- » Sterile bowls are packaged individually in easy-open peel pouches
- » Non-sterile bowls are packed in bulk

Item No.	Description	Pkg.
DYND50310	Sterile 8-oz. Bowl	50/cs
DYND50315	Sterile 16-oz. Bowl	100/cs
DYND50320	Sterile 32-oz. Bowl	50/cs
DN5C30016	Non-Sterile 16-oz. Bowl	250/cs
DYND50350	Non-Sterile 32-oz. Bowl	250/cs



DYNDX02763

DYND80000

DYND90000

Medicine Cups

These translucent cups accurately dispense liquid and dry medications.

- » Molded graduations
- » Two cups comply with 2015 ISMP metric-only safety standards

	DYND70000	DYND80000	DYND90000	DYNDX02763	DYNJ05195
Sterile	No	No	No	No	Yes
Volume: 1-oz.	•	•	•		
Volume: 2-oz.				•	•
Cubic Centimeters		•	•	•	•
Drams		•			
Fluid Ounces		•		•	•
Milliliters	•	•	•		
Tablespoons		•			
2015 ISMP Safety Standards	•		•		

Item No.	Description	Pkg.
DYND70000	1-oz. Non-Sterile, Graduated in ml only	5,000/cs
DYND80000	1-oz. Non-Sterile, Graduated in cc, drams, fl. oz., ml and tablespoons	5,000/cs
DYND90000	1-oz. Non-Sterile, Graduated in cc and ml	5,000/cs
DYNDX02763	2-oz. Non-Sterile, Graduated in cc and fl. oz.	1,840/cs
DYNJ05195	2-oz. Sterile, Graduated in cc and fl. oz., Pre-Packaged	100/cs

1.6 million

Number of medicine cups Medline produces every 24 hours using the world's most advanced manufacturing system of its kind.



Stainless Steel Components

Our stainless steel utensils offer a heavy-duty, long-lasting and easy-to-clean alternative to plastic.





Medicine Cups

This 2-oz. cup is graduated in cubic centimeters, fluid ounces, tablespoons and teaspoons for accurate and easy measuring.

» Non-sterile

Item No.	Description	Pkg.
DYND052Z	2-oz./60 ml	1/ea
DYND052A	2-oz./60 ml	12/cs

Iodine Cups

Easy to hold during procedures, these cups are ideal for soaking sponges.

» Choose from 6-oz. and 14-oz.

Item No.	Description	Pkg.
DYND056GZ	6-oz./180 ml	1/ea
DYND056G	6-oz./180 ml	24/cs
DYND0514GZ	14-oz./420 ml	1/ea

Graduated Measures

These measures have a large, smooth pouring spout and easy grip handle.

» Graduated in cubic centimeters and fluid ounces

Item No.	Description	Pkg.
DYND05T1062Z	16-oz./500 ml	1/ea
DYND05T1062	16-oz./500 ml	12/cs
DYND05T1063Z	32-oz./1000 ml	1/ea
DYND05T1063	32-oz./1000 ml	12/cs

Bedside Tray

The large, open bead on this tray's edges ensures trouble-free stacking, cleaning and holding.

» Covered corners for easy sterilization

» Seamless, oblong design

» 13.7 L x 9.8 W x 0.75" H (34.8 x 24.9 x 1.9 cm)

Item No.	Description	Pkg.
DYND0513FZ	Large Bedside Tray	1/ea

Griffin Style Beakers

The closed handle and seamless construction ensure safe transfer of liquids.

» Available in three sizes

Item No.	Description	Pkg.
DYND05125BZ	4-oz./125 ml	1/ea
DYND05250BZ	8-oz./250 ml	1/ea
DYND052WZ	64-oz./2000 ml	1/ea

Tumblers

These sturdy tumblers are virtually unbreakable and easy to clean.

» Round edges for comfortable drinking

» Easy to hold

Item No.	Description	Pkg.
DYND05T2107Z	7-oz./200 ml	1/ea
DYND05T2107	7-oz./200 ml	12/cs

Urinals

Seamless construction and an open flange make this urinal easy to clean and handle.

» Measured volume of 32-oz./1000 ml

Item No.	Description	Pkg.
DYND05T5GZ	Male Urinal, Stainless	1/ea



DYND05T5GZ



DYND05107Z



DYND0530Z



DYND05130Z



DYND05134Z

Solution Basins

These easy-clean basins feature a no-spin recessed bottom for stability on flat surfaces.

- » Large, rounded-radius corners for easy pouring
- » Available in three convenient sizes

Item No.	Description	Pkg.
DYND05130Z	5-qt./5.5 L	1/ea
DYND05130	5-qt./5.5 L	12/cs
DYND05134Z	7-qt./7.7 L	1/ea
DYND05134	7-qt./7.7 L	12/cs
DYND05136Z	9.5-qt./10.45 L	1/ea
DYND05136	9.5-qt./10.45 L	12/cs

Sponge Bowls

These wide-mouthed bowls are excellent for sponges.

- » Flat bottom and straight sides for extra stability
- » Available in a range of sizes

Item No.	Description	Pkg.
DYND0552Z	16-oz./500 ml	1/ea
DYND0575Z	20-oz./600 ml	1/ea
DYND05100Z	30-oz./900 ml	1/ea
DYND05106Z	36-oz./1100 ml	1/ea
DYND05148Z	48-oz./1400 ml	1/ea
DYND05107Z	64-oz./2000 ml	1/ea

Washbasins

These basins have a wide, rounded flange for easy pouring, handing and carrying.

- » Flat bottom and straight sides provide extra stability

Item No.	Description	Pkg.
DYND0530Z	4-qt./4 L	1/ea
DYND0530	4-qt./4 L	12/cs

Bedpans

Sturdy and autoclavable, these bedpans support up to 500 pounds.

- » Contour and fracture design options

Item No.	Description	Pkg.
DYND05T0015RZ	Contour Bedpan, Adult, 13.75 L x 11.5 W x 4" H (34.9 x 29.2 x 10.2 cm)	1/ea
DYND05TRZ	Contour Bedpan, Pediatric, 12.13 L x 9.75 W x 3" H (30.8 x 24.8 x 7.6 cm)	1/ea
DYND05T23RZ	Fracture Bedpan, Adult, 11.88 L x 9.5 W x 2.88" H (30.2 x 24.1 x 7.3 cm)	1/ea



DYND05T0015RZ



Emesis Basins

These easy-to-clean basins have a seamless design with rounded edges.

- » Suitable for surgical purposes, routine oral care and ear irrigation

Item No.	Description	Pkg.
DYND056Z	10-oz./300 ml	1/ea
DYND056	10-oz./300 ml	12/cs
DYND058Z	12-oz./350 ml	1/ea
DYND058	12-oz./350 ml	12/cs
DYND0510Z	26-oz./800 ml	1/ea
DYND0510	26-oz./800 ml	12/cs
DYND0512Z	50-oz./1500 ml	1/ea
DYND0512	50-oz./1500 ml	12/cs

Autoclavable Plastic Products

These plastic utensils are engineered and tested to withstand water temperatures up to 250-degrees Fahrenheit and 15 pounds of pressure for 15 to 20 minutes. Their blue color makes them easy to identify.

Item No.	Description	Pkg.
DYND88200	Contour Bedpan	12/cs
DYND88205	Pontoon Bedpan	12/cs
DYND88220	Fracture Bedpan	12/cs
DYND88552	Mega Fracture Bedpan	4/cs
DYND88235	Urinal w/Blue Lid	12/cs
DYND88310	5-qt. Round Washbasin	12/cs
DYND88320	500 ml Emesis Basin	12/cs
DYND88342	8-qt. Rectangular Washbasin	12/cs



Bedside Basins

For patients who need assistance with personal hygiene, our washbasins and other items help make bedside care quick and easy.




A Washbasins

These plastic basins accommodate bed baths and limb soaking.

- » Rolled rims make carrying easy
- » A soap dish can be clipped to rim for extra convenience
- » 6-qt. measures 9.5 L x 12.75 W x 4.5" H (24 x 32.4 x 11.4 cm)
- » 8-qt. measures 10 L x 13.5 W x 4.75" H (25.4 x 34.2 x 12 cm)

Item No.	Description	Pkg.
DYND80347	6-qt., Graphite	50/cs
DYND80306	6-qt., Gold	50/cs
DYNEC80347C	6-qt., Pigment Free	50/cs
DYND80342	8-qt., Graphite	50/cs
DYND80301	8-qt., Gold	50/cs
DYNEC80342C	8-qt., Pigment Free	50/cs



DYND80102

B Emesis Basins

Perfect for oral hygiene and ear irrigation, these sturdy basins feature molded graduation in milliliters and fluid ounces.

- » 500 ml (16-oz.) measures 8.5" (21.6 cm)
- » 700 ml (20-oz.) measures 9" (23 cm)

Item No.	Description	Pkg.
DYND80327	500 ml, Graphite	250/cs
DYND80321	500 ml, Gold	250/cs
DYNEC80327C	500 ml, Pigment Free	250/cs
DYND80367	700 ml, Graphite	250/cs
DYNEC80367C	700 ml, Pigment Free	250/cs

Sitz Baths

These localized baths provide soothing rectal and perineal soaking therapy. The graduated, 2000 ml solution bag has imprinted instructions for easy use and can be hung from any available fixture. The bowl fits over standard toilets.

- » Shut-off clamp on tubing controls solution flow
- » High-impact plastic bowl is contoured for patient comfort

Item No.	Description	Pkg.
DYND80102	Graphite	10/cs
DYND80101	Gold	10/cs
DYNEC80102C	Pigment free	10/cs

C Round Washbasins

These lightweight basins are excellent for sponge baths and limb soaking.

- » Plastic with rolled rims
- » 10.75 x 4.25" (27.3 x 10.8 cm)

Item No.	Description	Pkg.
DYND80317	5-qt., Graphite	50/cs
DYND80311	5-qt., Gold	50/cs
DYNEC80317C	5-qt., Pigment Free	50/cs

Soap Dishes

These ridged dishes keep soap handy and dry.

- » Clips onto rectangular washbasins
- » 3.75 L x 2.75 W x 0.5" H (9.5 x 7 x 1.3 cm)

Item No.	Description	Pkg.
DYND80491	Graphite	100/cs
DYND80492	Gold	100/cs
DYNEC80491C	Pigment free	100/cs

Round O.R. Basins

This non-sterile basin has molded graduations in quarts and centimeters.

- » Fits ring stands with an opening of 12–13.25" (30.5 x 33.7 cm)

Item No.	Description	Pkg.
DYNDX50346	6-qt., Light Blue	50/cs



DYND80491

Bedside Hydration

Our attractive hydration aids ensure patients always have water and other beverages within reach. We offer stocked and customizable carafes in sizes ranging from 16 to 32 ounces.





Insulated Carafes

These customizable carafes can be used for both hot and cold liquids. The 32-oz. holds enough water to last through the morning or afternoon. The smaller versions are perfect for patients with intake limitations.

- » Molded inner wall offers double-walled insulation
- » Lightweight with a secure-grip handle
- » Graduated in fluid ounces and cubic centimeters

Item No.	Description	Pkg.
DYC80540P	32-oz., Clear w/Blue Graduations and Blue Lid	48/cs
DYC80544	32-oz., Granite Body and Lid	48/cs
DYC80551	32-oz., White Body and Lid	48/cs
DYC80545	32-oz., Mauve Body and Lid	48/cs
DYND80558	32-oz., Clear w/Blue Graduations, Blue Lid and Handle	50/cs
DYND80557	24-oz., Clear w/Blue Graduations, Blue Lid and Handle	32/cs
DYC80541P	22-oz., Clear w/Blue Graduations and Blue Lid	72/cs
DYC80556	20-oz., Clear w/Blue Graduations, Blue Lid and Handle	50/cs

Accessories for Insulated Carafes

DYND80541	Jumbo Straw w/Safety Cap	104/cs
DYND80542	Poly Bag Liner	500/cs
DYND80540	Plastic Cup Liner for DYC80540P	500/cs

16-oz. Non-Insulated Carafe

These sporty carafes make excellent giveaways for children and seniors.

- » Clear with sparkles and a blue lid
- » Graduated in fluid ounces and cubic centimeters
- » Includes straw

Item No.	Description	Pkg.
DYND80552	16-oz., Sparkly w/Blue Graduations and Blue Lid	60/cs



You Design It. We Make It.

You want to attract patients and keep them loyal. Giving out branded and useful items like custom carafes from Medline can help your patients remember their stay and share their experience. Custom carafes also make great giveaways at community events. You supply the photo, graphic, logo or other artwork, and we take care of the rest.

Contact your Medline Representative to begin creating your custom carafe. Or use myownmug.com—it's quick and easy!

All About the Carafe

- » Your photo covers 93% of the carafe
- » In-mold labeling technology produces a high-quality, photo-like image that won't scratch off or fade
- » Made of rugged polypropylene, BPA-free
- » Dishwasher safe up to 165°F
- » Lid with hinge doors locks in place, facilitating straw placement, pouring and drinking directly from the lid
- » Recessed bottom prevents condensation so the carafe cannot form a seal with the moisture and "lock" itself to the tabletop
- » Double-wall thermal design keeps liquids cold without sweating
- » Capacity: 32-oz., 24-oz., 22-oz. and 20-oz.






Non-Insulated Plastic Pitchers

These pitchers have a detachable lid that pops off with one quick thumb movement.

- » New straw opening in lid lets patients drink directly from pitcher
- » Easy-pour spout holds back ice
- » Reinforced tabs under handle adds stability when pitcher is full of water or ice

Item No.	Description	Pkg.
DYND80535	32-oz., Graphite	100/cs
DYND80521	32-oz., Gold	100/cs
DYNEC80535C 	32-oz., Pigment Free	100/cs



Accessories for Pitchers

Designed to fit our plastic pitchers (DYND80535 and DYND80521).

- » Insulating foam and plastic liners help prevent surface sweating
- » White foam outer jacket provides added insulation for keeping liquids cold
- » Pitchers can be packaged with any inner or outer liner through our non-sterile kit product line

Item No.	Description	Pkg.
DYND80532	Foam Outer Jacket	200/cs
DYND80526A	Foam Inner Liner, White	500/cs
DYND80528	Graduated Plastic Inner Liner, Translucent	500/cs
DYND80541	Flexible Straw w/Safety Cap	104/cs
DYND80565	Replacement Lids	50/cs

Foam Pitchers

Offering excellent insulation, these thick, white foam pitchers keep beverages at the preferred temperature.

- » Cap rotates for effortless pouring or ice retention
- » Easy-grip, comfortable handle
- » Replacement plastic liners are available

Item No.	Description	Pkg.
DYND80412	32-oz. Foam Pitcher and Liner	50/cs
DYND80408	32-oz. Foam Pitcher, No Liner	50/cs
DYND80414	Liner for Inside Pitcher	500/cs



Carafe With Cup Cover

Great for patients with limited dexterity, these carafes are easy to open, close and hold.

- » Snap-on top is also the lid
- » Wide mouth simplifies pouring
- » Broad base offers added stability

Item No.	Description	Pkg.
DYND80403	32-oz., Graphite	20/cs

Service Trays

These trays hold different combinations of carafes, pitchers and tumblers.

- » Ridges on the Deluxe version help keep accessories dry and prevent sticking to the tray
- » Deluxe trays measure 10 L x 6 W x 0.5" H (25.4 x 15.2 x 1.27 cm)
- » Standard trays measure 4.6 L x 7.5 W x 0.25" H (11.6 x 19 x 0.6 cm)

Item No.	Description	Pkg.
DYND80445	Deluxe Tray w/Ridges, Graphite	100/cs
DYND80441	Deluxe Tray w/Ridges, Gold	100/cs
DYND80438	Standard Tray, Smooth, Graphite	500/cs
DYNEC80445C	Deluxe Tray w/Ridges, Pigment Free	100/cs
DYND87605	Pitcher, Tumbler and Tray Set, Graphite	12/cs



Tumblers

These 9-oz. tumblers have molded graduations in fluid ounces, cubic centimeters and milliliters.

- » Available in a variety of colors
- » 25 tumblers per sleeve, 20 sleeves per case

Item No.	Description	Pkg.
DYND80454	Molded Graduations, Graphite	500/cs
DYND80456	Molded Graduations, Gold	500/cs
DYND80450	Molded Graduations, Translucent	500/cs
DYND80460	Molded Graduations, Clear	500/cs
DYND80455	Imprinted Graduations in Pink, Translucent	500/cs



DYND87605

Bedpan & Urinal Holders

These wall-mountable and cage-style holders offer convenient storage for bedpans, urinals, graduates and specimen pans.



Urinal Holders with Vinyl Strap

Perfect for keeping urinals within easy reach of the patient, these holders feature a plastic cage and vinyl strap that fasten to bedrails, wheelchairs or walkers.

- » Adjusts to multiple angles using a “swiveling” action
- » 10 H x 4.38 W x 4.38” from exterior, (25.4 x 11.1 x 11.1 cm)
- » Standard options have a 6” strap, long “L” versions have an 8” strap

Item No.	Description	Pkg.
PGA201	EZP Urinal Holder w/Vinyl Strap	50/cs
PGA201-25	EZP Urinal Holder w/Vinyl Strap	25/cs
PGA201L	EZP Urinal Holder w/Long Vinyl Strap	50/cs
PGA201L-25	EZP Urinal Holder w/Vinyl Strap	25/cs
PGA201SHRT	EZP Urinal Holder, Armless w/Vinyl Strap	50/cs
PGA201SHRT-25	EZP Urinal Holder, Armless w/Vinyl Strap	25/cs



PGA201

Urinal Holders with Hook-and-Loop Strap

These holders have a hook-and-loop strap that attaches easily to bedrails, wheelchairs or walkers.

- » Adjusts to multiple angles using a “swiveling” action
- » 10 H x 4.38 W x 4.38” from exterior, (25.4 x 11.1 x 11.1 cm)

Item No.	Description	Pkg.
PGA101	EZP Urinal Holder w/Hook and Loop Strap	50/cs
PGA101-25	EZP Urinal Holder w/Hook and Loop Strap	25/cs
PGA101SHRT	EZP Urinal Holder, Armless w/Hook and Loop Strap	50/cs
PGA101SHRT-25	EZP Urinal Holder, Armless w/Hook and Loop Strap	25/cs



PGA101

Bedpan Holders

These wall-mounted steel frames are powder-coated to resist rusting.

- » Holds one bedpan and two urinals
or one bedpan, one urinal and one graduate

Item No.	Description	Pkg.
PCM1043	Bedpan/Urinal Holder, 13 W x 12.5 H x 5.5” From the Wall (33 x 31.8 x 14 cm)	6/cs
PCM1054	Bedpan/Graduate/Urinal Holder, 20.5 W x 13 H x 5.5” From the Wall (52.1 x 33 x 14 cm)	6/cs



PCM1054

Specimen Pan Holders

These wall-mounted enamel wire frames keep specimen pans and urine hats within easy reach.

- » PCM1052 holds one specimen pan
- » PCM1053 holds two specimen pans

Item No.	Description	Pkg.
PCM1052	14 W x 7 H x 2.25” From the Wall (36.6 x 17.8 x 5.7 cm)	6/cs
PCM1053	8.5 W x 2 H x 7.5” From the Wall (21.6 x 5.1 x 19.1 cm)	6/cs



Medline Industries, Inc. One Medline Place, Mundelein, IL 60060

FOLLOW US [f](#) [in](#) [t](#) [BLOG](#)

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

We reserve the right to correct any errors that may occur within this brochure. Some products may not be available for sale in Mexico or Canada.

© 2016 Medline Industries, Inc. All rights reserved. Durafit, greensmart and Medline are registered trademarks of Medline Industries, Inc. MKT1657154 / LIT118 / 2.5M / JPG / 3