



# SANI-CLOTH® AF3

## GERMICIDAL DISPOSABLE WIPE

Protecting your patients has never been easier. Fragrance Free - Ideal for use around patients and staff, especially those with respiratory sensitivities.

### Sani-Cloth® AF3 Benefits:

- The first alcohol free, quat based disinfectant wipe with a 3 minute kill time
- Perfect for equipment and other hard non-porous surfaces sensitive to alcohol
- Fragrance, alcohol, bleach, acid, phenol, acetone, and ammonia free
- Meets CDC, OSHA and CMS Tag F441 guidelines
- Bactericidal, Tuberculocidal, Virucidal



EPA Reg. No. 9480-9

Contact Time  
(Minutes)



# of Microorganisms  
Effective Against

44

Disinfection  
Formulation

Alcohol  
Free Quat

Effective against 44 microorganisms in 3 minutes including the following MDROs, bloodborne pathogens and viruses:<sup>1</sup>

- *Klebsiella pneumoniae* - NDM-1
- *Klebsiella pneumoniae* - Carbapenem Resistant
- ESBL producing *E. coli*
- CA-MRSA USA 300 and 400
- RSV
- HBV
- HCV
- HIV

**SANI-CLOTH®**

The #1 brand of disinfectant wipes in healthcare.<sup>2</sup>


<sup>1</sup>See Technical Data Bulletin for a complete list of microorganisms. <sup>2</sup>GHX Intelligence Data 2010.



**ACCESSIBILITY:** Deep well lid provides ample space to store and access next wipe.



**COMPLIANCE:** Prominent contact time for easy identification and improved compliance with accreditation standards.



**IDENTIFICATION:** Categorization of pathogens allows for quick ID of relevant efficacy claims.

## ORDER INFORMATION

SANI-CLOTH® AF3	Reorder No.	Wipe Size	Packaging	Case Wgt. (lbs)	Case Cube (ft.)	TI/HI
A. Large Canister	P13872	6" x 6.75"	160/can, 12 cans/cs	26.4	1.41	10/3
B. Extra-Large Canister	P63884	7.5" x 15"	65/can, 6 cans/cs	15.3	1.010	10/4
C. Large Individual Packets	H59200	5" x 8"	50/box, 10 boxes/cs	7.99	0.540	20/4
D. Extra-Large Individual Packets	U27500	11.5" x 11.75"	50/box, 3 boxes/cs	7.6	0.500	24/3



A.

B.

## SANI-CLOTH® EDUCATIONAL & COMPLIANCE ACCESSORIES

Developed to assist you in facility protocol compliance, as well as meet The Joint Commission and CMS requirements.<sup>3,4</sup>



**Infection Prevention Pak** for mobile infection prevention - hold packets, gloves and other PPE



**PDI 3-in-1 Universal Sani-Bracket** for wall/pole mounting of a Sani-Cloth® Canister



**Sani-Canister Caddy** for placement on flat surfaces (bilingual)



**Sani-TAG™ Equipment ID System** for identification of clean/soiled equipment



**Educational wall charts and signage** (bilingual)

www.pdipdi.com

1-800-999-6423

© 2012 PDI.

<sup>3</sup>Infection Prevention and Control (IC) Standard IC.02.02.01 and Elements of Performance (EP) 1 and 2, The Joint Commission, October 2009. <sup>4</sup>Guideline for Disinfection and Sterilization in Healthcare Facilities, Centers for Disease Control and Prevention, 2008. Sani-Cloth® and Sani-Bracket® are registered trademarks of Professional Disposables International, Inc. Sani-Canister Caddy™ and Sani-TAG™ are trademarks of Professional Disposables International, Inc. PDI04125798

**PDI**  
BE THE DIFFERENCE™