

Wallaby Pediatric Wheelchair

HCPCS: E1236



FIGURE A



- Standard Transit Tie Down Brackets
- WC19 Compliant when used with WB8021 for users between 50 and 80 lbs
- Height adjustable push handles
- Flip back padded desk arms
- Pelvic belt standard
- Push to lock wheel locks
- 22" Rear wheels with 6" front casters
- Swingaway footrests with tool free height adjustment and calf strap

PRODUCTS

Description

| | |
|-------------|--------------------------------|
| WB1200-2GFR | Wallaby, 12", Fire Truck Red |
| WB1400-2GJB | Wallaby, 14", Jet Fighter Blue |

SPECIFICATIONS

Description

| | |
|--|------------|
| Armrest Height | 7" |
| Backrest Height (FC 2000N FC 4000N) | 14" |
| Product Weight with footrests | 37 lbs. |
| Seat To Floor Height | 19" |
| Weight Capacity | 150 lbs |
| Weight Capacity (13046 13046-HD) | 150 lbs. |
| Weight Capacity when used in transport | 50-80 lbs. |

ACCESSORIES

Description

| | |
|---------|---------------------------------------|
| WB 8021 | Headrest Extension, Wallaby,1ea, CS |
| WB ELR | Elevating Leg Rests, Wallaby, 1pr, CS |
| STDS819 | Anti Tippers,Cruiser,1pr/bx, CS |