



## Mayo-Instrument Stand Double Post

- Double post provides added support for use in surgical procedures
- Removable stainless steel tray measures 19" x 12.5"
- Adjustable height from 32.5"– 52"
- Two hooded casters

PRODUCTS	Description
13045	1/cs

SPECIFICATIONS	Description
Tray Dimensions	19" (W) x 12.5" (D)
Height	32.5" - 52"
Base Dimensions	21" (W) x 19" (D)
Weight	11 lbs
Carton Shipping Weight	14 lbs.
Warranty	Limited Lifetime

ACCESSORIES	Description
13035T	1/cs