




Padded Transfer Bench
Tool-free Assembly Back, Legs and Arms



FIGURE A



FIGURE B



INSET D

- Comfortable cushioned seat and backrest
- Tool-free assembly of back, legs and arm (Figure A - C)
- New pinch-free cover allows for push pins to be depressed without pinching fingers
- A-frame construction provides stability
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs

PRODUCTS

Description

12005KD-2	Knocked Down, Tool-free Assembly Back, Legs and Arms, 2/cs
12005KD-1	Knocked Down, Tool-free Assembly Back, Legs and Arms, 1/cs

SPECIFICATIONS

Description

Outside Legs	29" (W) x 19.75" (D)
Seat Dimensions	24" (W) x 16" (D)
Seat Height	17.75" (L) x 21.75" (H)
Weight	11.75 lbs.
Carton Shipping Weight (-1 -2)	14 lbs. 26.5 lbs.
Weight Limit	400 lbs.
Warranty	Limited Lifetime

ACCESSORIES

Description

12011SS	Sm Suction Cups Pr., 1pr/bx
12011SL	Lg Suction Cups Pr., 1pr/bx