

Plastic Tub Rail



RETAIL BOX

- Attractive, ergonomically-designed tub rail provides a comfortable, safe grip
- Handle is perpendicular to tub
- 16.5" in height
- Tool-free, easy to install
- Width adjusts from 3" - 6.25"
- Retail packaged

PRODUCTS

Description

12044	White, 1/bx
-------	-------------

SPECIFICATIONS

Description

Adjustable Clamp	3" (smallest) 6.25" (widest)
Height	16.5"
Carton Shipping Weight	7 lbs.
Weight	5 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime