



## Rebel Wheelchair Single Axle



FIGURE A

INSET C

- Quick-release rear wheels allow for easy storage or transport (Figure A)
- Fold-down back and folding frame allows transport in almost any trunk (Figure B)
- Easily removable swing-away footrests are also convenient for storage
- Tool-free removable front casters
- Desk-style armrests make transfers easy and accessibility better
- Lightweight and extremely maneuverable, easy to navigate tight areas
- Large accommodating seating surface
- Carton size makes the Rebel perfect for retail

PRODUCTS	Description
RTLREB18DDA-SF	Rebel W/C 18 FlameRd DetDsk SwgFt, 1/cs

SPECIFICATIONS	Description
Seat Width	18"
Overall Width	25.5" (open)   10.5" (closed)
Overall Height	35"
Seat-To-Floor Height	19.5"
Back Height	17"
Seat Dimensions	18" (W) x 16" (D)
Overall Length (with riggings)	43"
Weight	37 lbs.
Carton Shipping Weight	42 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

ACCESSORIES	Description
STDS806	Anti-Fold Bar, 1/ea
STDS1040S	Universal Cup Holder, 1/ea
STDS831	Heel Loops, 1/bx
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS801GT	Wheel Lock Extension, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea

PARTS	Description
STDS806	Anti-Fold Bar, 1/ea
STDS1040S	Universal Cup Holder, 1/ea

PARTS	Description
STDS831	Heel Loops, 1/bx
STDS850	Seat Belts, Auto-Clasp Type Closure, 1/ea
STDS851	Seat Belts, Hook and Loop Fastener Type Closure, 1/ea
STDS855	Seat Belts, Auto-Clasp Type Closure, Bariatric, 1/ea
STDS856	Seat Belts, Hook and Loop Fastener Type Closure, Bariatric, 1/ea
STDS803	Universal Oxygen "E" Cylinder Holder, 1/ea
STDS804	Universal Oxygen "E" Cylinder Holder with IV Pole Attachment, 1/ea
STDS801GT	Wheel Lock Extension, 1/ea
STDS1034	Wheelchair Cane/Crutch Holder, 1/ea