

## Transfer Tub Bench

HCPCS: E0247



FIGURE A




FIGURE B



FIGURE C

- Tool-free assembly of back, legs and arm (Figure A - C)
- A-frame construction provides stability
- Durable blow-molded plastic bench and backrest
- Height adjusts in ½" increments with unique "Dual Column" extension legs
- Reversible to accommodate any bathroom
- Pinch-free lever allows for push pins to be depressed without pinching fingers

PRODUCTS	Description
12011KD-1	Transfer Tub Bench, CS
12011KD-2	Transfer Tub Bench, CS
SPECIFICATIONS	Description
Carton Shipping Weight	13 lbs.   13 lbs.   24 lbs.
Outside Legs	32" (W) x 23" (D)
Seat Dimensions	26" (W) x 18.5" - 19.5" (D)
Seat Height	17.5" - 22.5"
Warranty	Limited Lifetime
Weight	10 lbs.
Weight Capacity	400 lbs
ACCESSORIES	Description
12011SS	SuctionCups,Sm.12011,1pr/bx, BX
12011SL	SuctionCups,Lg.12011,1pr/bx, BX