



1705

UP-DRAFT II® OPTI-NEB® NEBULIZER

Designed to be a universal nebulizer, the Up-Draft II gives you added flexibility in flow rates, delivering flow rates ranging from 5 to 9 LPM and making it a popular choice for those looking to consolidate options into one high-performance device.

- Nebulizes at angles up to 45° to accommodate special treatment requirements
- Produces fine, dense particles for effective therapy
- Easy-seal, threaded cap and 8 cc capacity jar
- Anti-spill design prevents loss of medication in any position
- Jet stay-in-place design for intentional nebulization
- Available in universal and individualized versions
- Recommended flow rates: 5 to 9 LPM

Individualized

Updraft II® Opti-Neb® Nebulizer

ITEM CODE	NEBULIZER	TEE	MOUTHPIECE	RESERVOIR TUBE	7 FT STAR LUMEN TUBING	CONNECTOR
1705	•				•	Standard
1706	•				•	Standard

II® Opti-Neb® Nebulizer
quality aerosols at flows
choice for many facilities
performance product.

ays in place unless
onally removed
ble with standard or
sal tubing connector

recommended flow rate:
LPM

ally packaged, 50 to a case.

MASK	
	Adult, Elongated (Item 1083)
	Adult (Item 1084)