

Premium Raised Toilet Seat with Removable Arms, Standard Seat



Α	Overall Seat Width	13.75"
В	Inside Seat Width	9.75"
С	Overall Seat Depth/Length	16.5"
D	Inside Seat Depth/Length	11.25"

Е	Overall Seat Height	3.5"
F	Overall Width (with Handles)	21.75"
G	Overall Height (with Handles)	9.25"
Н	Mounting Hole Width	5.5"


Premium Raised Toilet Seat with Removable Arms, Elongated Seat



А	Overall Seat Width	14"
В	Inside Seat Width	10"
С	Overall Seat Depth/Length	18.5"
D	Inside Seat Depth/Length	13"

E	Overall Seat Height	3.5"
F	Overall Width (with Handles)	21"
G	Overall Height (with Handles)	10"
Н	Mounting Hole Width	5.5"


Premium Plastic, Raised, Regular/Elongated Toilet Seat with Lock



Α	Overall Seat Width	15.5"
В	Inside Seat Width	8"
С	Overall Seat Depth/Length	16.5"
D	Inside Seat Depth/Length	9"

Е	Overall Seat Height	5"
F	Overall Width (with Handles)	22.5"
G	Overall Height (with Handles)	9"
Н	Seat Lock Range	11.5"-15.5"

Premium Plastic Raised, Regular/Elongated Toilet Seat, with Lock



Α	Overall Seat Width	17"
В	Inside Seat Width	8"
С	Overall Seat Depth/Length	16.5"

D	Inside Seat Depth/Length	9"
Е	Overall Seat Height	4.5"
F	Seat Lock Range	11.5-15.5"