



Free-Standing Silver Vein Trapeze with Base and Wheels

- Designed to assist larger individuals to change positions while in bed and aid in the transfer from bed to chairs
- Easy-to-assemble (Requires only a drill, rubber mallet and socket with wrench)
- The overhead boom swivels from 0° to 180°
- Chain is adjustable from edge of headboard to the end of boom by sliding to desired position

PRODUCTS

Description

13039SV	For use by patients up to 1000 lbs, 1/cs
---------	--

SPECIFICATIONS

Description

Total Dimensions	73" (H) x 43" (W) x 48.5" (D)
Inside Width Between Legs	37.5"
Outside Width Between Legs	48"
Boom Length	36"
Weight	100 lbs.
Carton Shipping Weight	110 lbs.
User Weight Limit	1000 lbs.