

All foams soak it up, KerraFoam® *locks it up*

LOCKS IN MOISTURE AND HELPS TO REDUCE MACERATION

Older generation foam dressings claim superior fluid handling by virtue of their MVTR rates not absorption rates. However, MVTR is fine under laboratory conditions but in a clinical situation, with changes in humidity under bandages and temperature changes means that MVTR can reduce, substantially increasing risk of maceration!

EXCELLENT RETENTION CAPABILITIES

KerraFoam Non-adhesive has excellent fluid retention capabilities in comparison to other foam dressings – retaining more than the two leading foam brands.² Instantly locking harmful excess exudate away from the skin helps to prevent maceration and wound edge breakdown.

HIGHLY COMFORTABLE AND CONFORMABLE

Non-adherent wound contact layer for use on fragile and sensitive skin.

Conforms well to awkward areas without limiting movement.





Super Smart Foam Dressing

• OUTER POLYURETHANE FILM LAYER

A smooth, water-resistant and bacterial barrier, helping to create a moist protective healing environment.

ABSORBENT CORE AND FOAM PAD

Combined to both absorb and lock-in exudate, even under compression!

UNIQUE HORIZONTAL WICKING LAYER

Efficiently distributes exudate, reducing the risk of maceration. Increases capacity for absorption. Helps reduce the number of dressing changes. Helps to prevent strikethrough.

SOFT, NONWOVEN MATERIAL

Non-adherent contact layer to provide comfort for the patient



ORDER INFORMATION			
Product Code	Size	Dressings per box	HCPCS
CWL1116	4 x 5 in	10	A6210
CWL1117	6 x 6 in	10	A6210
CWL1118	4 x 8 in	10	A6210
CWL1119	8 x 8 in	10	A6211
CWL1131	8 x 12 in	10	A6211
CWL1120	8 x 17 in	5	A6211

To place an order please contact your medical supply distributor or call 855-522-2211



Crawford Healthcare Inc. | 2005 South Easton Rd. | Suite 203 | Doylestown, PA 1890
Tel 855-522-2211 | Email us.info@crawfordhealthcare.com

www.crawfordhealthcare.com/u