

Surgical Cautery Tips

McKesson Argent surgical cautery tips are ideal for procedures requiring precise homeostasis: cardiovascular, thoracic, occult-plastic, orthopedic and hand surgery procedures. These cautery tips are also excellent for evacuation of subungual hematoma. Available in a variety of tip configurations and are battery operated with 3V AA batteries.

Features and Benefits

- » Ideal for procedures requiring precise homeostasis: cardiovascular, thoracic, occult-plastic, orthopedic and hand surgery procedures
- » Available in a variety of tip configurations and are battery operated with 3V AA batteries
- » Temperature at activation can reach 1,800°F – 2,200°F with a constant burn time of 10 minutes
- » Each cautery includes a protective cover cap for storage and safe disposal
- » Sterile
- » Not made with natural rubber latex



MFR #	Description	QTY
231	High Temp Cautery, Fine Tip	1/pouch, 10 pouches/bx
243	High Temp Cautery, Loop Tip	1/pouch, 10 pouches/bx
244	High Temp Cautery, Elongated Tip	1/pouch, 10 pouches/bx
22-100	Low Temp Cautery, Ophthalmic Tip	1/pouch, 10 pouches/bx
22-102	Low Temp Cautery, Fine Tip	1/pouch, 10 pouches/bx

Satisfaction Guaranteed
 If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Surgical Cautery Tips

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908