


SPECIFICATIONS

| | Go-Go® Ultra X | Go-Go® LX w/ CTS Suspension | Go-Go Elite Traveller® | Go-Go® Sport | Go Chair® | Go-Go® Folding Scooter |
|--|--|--|--|--|---|--|
| Max. weight capacity: | 260 lbs. | 300 lbs. | 300 lbs. | 325 lbs. | 300 lbs. | 250 lbs. |
| Maximum speed: ¹ | Up to 4 mph | Up to 4.4 mph | Up to 4 mph | Up to 4.7 mph | Up to 3.7 mph | Up to 4 mph |
| Ground clearance: ² | 2.25" at midframe 1.25" at motor | 1.5" at motor | 2" | 3" at midframe 1.5" at motor | 1.6" (at anti-tip bracket) 1.9" (at frame) | 2.3" |
| Turning radius: ² | 3W: 33" / 4W: 44" | 3W: 33" / 4W: 44.75" | 3W: 33" / 4W: 44" | 3W: 36" / 4W: 50.5" | 25.4" without foot platform | 32" |
| Overall length: ² | 3W: 37" / 4W: 39.5" | 3W: 39.75" / 4W: 41" | 3W: 37" / 4W: 39.5" | 3W: 42" / 4W: 43" | 33" with foot platform | 37" |
| Overall width: ² | 19.5" | 20.5" | 19.5" | 21.25" | 22.6" | 19" |
| Battery charger: | Off-board, 2 amp | Off-board, 2 amp | Off-board, 2 amp | Off-board, 2 amp | Off-board, 2 amp | 2-amp, off-board (SLA) or 2.5-amp, off-board (lithium-ion) |
| Per-charge range: ^{1,5} | Up to 8 miles | 3W: Up to 9.44 miles 4W: Up to 10.5 miles | 3W: Up to 10.4 miles ⁸ 4W: Up to 9.7 miles ⁸ | 3W: Up to 12.7 miles 4W: Up to 12 miles | Up to 8.74 miles | Up to 12.5 miles (lithium-ion) Up to 9.3 miles (SLA) |
| Battery weight: ⁴ | 9 lbs. each ⁶ | 1.4 lbs. each ⁷ | 1.4 lbs. each ⁷ | 1.4 lbs. each ⁷ | 31 lbs. (battery pack) | 12 AH SLA - 9 lbs. ea. 11.5 AH lithium-ion - 5 lbs. ea. 11.5 AH SLA (2 required) 11.5 AH lithium-ion (1 required) |
| Batteries (2 required) ³ (12V, deep cycle) | 12 AH | 18 AH | 18 AH | 18 AH | 18 AH | 51.4 lbs. |
| Weight w/o battery pack ² | 3-wheel: 75 lbs. 4-wheel: 78.5 lbs. | 3-wheel: 85.4 lbs. 4-wheel: 87.4 lbs. | 3-wheel: 75 lbs. 4-wheel: 78.5 lbs. | 3-wheel: 83.2 lbs. 4-wheel: 99.2 lbs. | 36 lbs. (battery pack) | 51.4 lbs. (unit does not disassemble - unit weight without battery pack) |
| Heaviest piece: ² | 3-wheel: 29.5 lbs. (front section) 4-wheel: 34 lbs. (front section) | 3-wheel: 34 lbs. (front section) 4-wheel: 36.6 lbs. (front section) | 3-wheel: 29.5 lbs. (front section) 4-wheel: 34 lbs. (front section) | 3-wheel: 35 lbs. (front section) 4-wheel: 44 lbs. (front section) | 36 lbs. (rear section) | |



¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
²Due to manufacturing tolerances and continued product improvement, this specification is subject to a variance of +/- 3%.
³AGM or gel-cell type recommended.
⁴Battery weight may vary based on manufacturer.
⁵Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
⁶12 AH battery pack weight is 18 lbs.
⁷18 AH battery pack weight is 28 lbs.
⁸Range tested with 18AH batteries.

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.



GO GO TRAVEL MOBILITY

GO CHAIR

GO GO FOLDING SCOOTER



GO GO ULTRAX

GO GO LX WITH CTS SUSPENSION



GO GO ELITE TRAVELLER

GO GO SPORT



Pride
 Mobility Products Co.
 — LIVE YOUR BEST —

Pride
 Mobility Products Co.
 — LIVE YOUR BEST —

GO•GO®

ULTRA X



FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand
- Pride®'s exclusive black, non-scuffing tires
- Auto-connecting front to rear cable

GO•GO®

LX WITH CTS SUSPENSION



FEATURES

- CTS Suspension (Comfort-Trac front and rear suspension)
- Feather-touch disassembly permits simple frame separation with only one hand
- Includes 2 sets of easily changeable colored panels in Red and Blue
- Pride®'s exclusive black, non-scuffing tires
- Auto-connecting front to rear cable

GO•GO®

ELITE TRAVELLER®

&

GO•GO®

SPORT



FEATURES

- Feather-touch disassembly permits simple frame separation with only one hand
- Pride®'s exclusive black, non-scuffing tires
- Auto-connecting front to rear cable
- Includes two sets of easily changeable colored panels in two exciting colors - Red and Blue

GO•CHAIR®



FEATURES

- 18 amp batteries
- One hand feather-touch disassembly
- Fits in most vehicle trunks
- Top speed of 3.7 mph and enough power to get you where you want to go
- Choice of 6 colors

GO•GO®

FOLDING SCOOTER



FEATURES

- Unique folding design
- Folds in just three easy steps
- Height adjustable tiller
- Easy to transport at home, hotels, or airports
- Can be stowed upright or laid down