



**KAWASUMI**

KAWASUMI LABORATORIES AMERICA, INC.

[www.kawasumiamerica.com](http://www.kawasumiamerica.com)

**1-800-KAWASUMI (529-2786)**

P.O. Box 24355 • Tampa, FL 33623-4355

Phone: 813-630-5554 • Fax: 813-630-5033

e-mail: [orders@kawasumiamerica.com](mailto:orders@kawasumiamerica.com)



**KAWASUMI LABORATORIES, INC.**

Shinagawa Intercity Tower B, 2-15-2,

Konan, Minato-ku, Tokyo 108-6109, Japan

Kawasumi continuously improves and develops new products. Visit our website for more details!

PC-04





**KAWASUMI**  
KAWASUMI LABORATORIES AMERICA, INC.

U.S. Product Catalog

## K-Shield Advantage® Winged Blood Collection Sets

CE 0123

- Provides a safe collection environment for both clinician and patient
- Integrated safety device eliminates extra steps for disposal
- Fast and simple one-step activation with an “audible” click
- Needle is covered during removal, ensuring confidence in safety
- Superior needle sharpness\*
- Latex free - for patients, practitioners and end-users



\*NAMSA needle sharpness report available upon request.

For product codes, box and case sizes please see page 12





**K-Shield Advantage®**  
**Winged Blood Collection Sets**  
**with pre-attached Tube Holder**

CE 0123




- Provides a safe collection environment for both clinician and patient
- Integrated safety device eliminates extra steps for disposal
- Fast and simple one-step activation with an “audible” click
- Needle is covered during removal, ensuring confidence in safety
- Superior needle sharpness\*
- Convenience of a pre-attached vacuum tube holder
- Latex free - for patients, practitioners and end-users



\*NAMSA needle sharpness report available upon request.



For product codes, box and case sizes please see page 12






### K-Shield® Winged Blood Collection Sets $\text{CE 0123}$ with Multiple Sample Luer Adapter

- Provides a safe collection environment for both clinician and patient
- One handed operation
- Fast and easy to use integrated safety device reduces accidental needle sticks
- Needle point is immediately covered upon removal from venipuncture site, ensuring confidence in safety
- Superior needle sharpness\*
- Latex free - for patients, practitioners and end-users





### Standard Winged Blood Collection Sets $\text{CE 0123}$ with Multiple Sample Luer Adapter

- Low cost winged collection set with multiple sample luer adapter
- Self-sealing synthetic needle sheath
- Flexible wings with textured surface for secure grip
- Easy-grip hub fits tube holders with anti-stick devices
- Superior needle sharpness\*
- Latex free - for patients, practitioners and end-users



\*NAMSA needle sharpness report available upon request.

For product codes, box and case sizes please see page 12



### K-Shield Advantage® Small Vein Infusion Sets

- Provides a safe collection environment for both clinician and patient
- Integrated safety device eliminates extra steps for disposal
- Fast and simple one-step activation with an "audible" click
- Needle is covered during removal, ensuring confidence in safety
- Superior needle sharpness\*
- Latex free - for patients, practitioners and end-users

### Standard Small Vein Infusion Sets CE 0123

- Flexible wings with textured surface for secure grip which are color coded for easy identification
- Siliconized needle for easy insertion
- Kink resistant tubing
- Superior needle sharpness\*
- Latex free - for patients, practitioners and end-users



### Intermittent Injection Site

- Low cost injection site
- Male luer connector
- Allows clinician to perform multiple injections when attached to a winged infusion set
- Latex free - for patients, practitioners and end-users



Intermittent  
injection Site


\*NAMSA needle sharpness report available upon request.



### K-Shield Advantage® Blood Drawing Kits

- Low-cost therapeutic blood drawing
- Fast and simple one-step activation with an “audible” click
- Needle-free sampling port
- Does not contain unnecessary anticoagulant solution
- 600 mL blood bag
- Pre-attached 16G venipuncture needle
- Pinch clamp for easy prevention of backflow
- 48-inch tubing length
- Volume graduation marks
- Latex free - for patients, practitioners and end-users

BK-64A




For product codes, box and case sizes please see page 13







## K-Shield Advantage® Phlebotomy Sets

- Low-cost design incorporates pre-attached needles for accessing all vacuum bottles
- Fast and simple one-step activation with an “audible” click
- Pre-attached 16G venipuncture needle
- 16G needle with backeye optimizes blood flow
- Needle with wing for improved needle control
- Roller clamp to control blood flow
- 36-inch tubing length
- Latex free - for patients, practitioners and end-users



## Standard Phlebotomy Sets

- Sets are available with and without pre-attached 16G venipuncture needle
- No safety device
- Latex free - for patients, practitioners and end-users



For product codes, box and case sizes please see page 13




**K-Shield®  
Port Access Infusion Sets**

CE 0123

- **Low profile design** – streamlined design helps make dressing simple and efficient and improves patient comfort during long infusion times
- **Simple safety device activation** – clinician-friendly and easy to use; audible click lets you know the needle is locked and enclosed during activation
- **Non-coring siliconized needle** – minimizes irritation while the port is being accessed or when the needle is removed
- **Spacers on skin plate** – allows air flow and minimizes moisture between patient’s skin and underneath skin plate
- **Non-DEHP fluid path** – offers a safe, DEHP-free fluid path
- **Latex free** – for patients, practitioners and end-users

**Non-DEHP  
Fluid Path**




**Standard Port Access Infusion Sets**

CE 0123

- Low-cost port access sets available in a variety of models
- Non-DEHP fluid path
- Non-coring needle
- Needle hub and wing design offers secure and accurate needle placement
- Latex free - for patients, practitioners and end-users

For product codes, box and case sizes please see page 14



**Standard Port Access Infusion Sets** CE 0123

- Low-cost port access sets available in a variety of models
- Non-DEHP fluid path
- Non-coring needle
- Needle hub and wing design offers secure and accurate needle placement
- Latex free - for patients, practitioners and end-users

For product codes, box and case sizes please see page 14



**Non-DEHP  
Fluid Path**




### Gravity IV Administration Sets

- **Flexible** - extremely flexible tubing, does not kink
- **Long Length** - long length sets for accommodating a wide variety of surgical procedures where the patient is distant from the IV pole
- **Needle Safe** - Available with NAC® (Needleless Access Connectors), providing a needle-free environment
- **Available** with .22 micron, low protein binding in-line filter
- **Wide selection** of sets and components
- **Non-DEHP fluid path** - all sets feature a Non-DEHP fluid path including tubing and drip chamber
- **Latex Free** - for patients, practitioners and end-users
- **Economical** - high-quality, low-cost gravity IV administration sets for clinical needs as well as meeting materials management budget requirements


Code #	Description (50 sets per case)	Length	Drip Rate	Priming Volume
IV-009N	Standard set with NAC®	84 inches	10 drops/mL	18 mL
IV-009Y	Set with needle-pierceable Y site	84 inches	10 drops/mL	18 mL
IV-013Y	Micro-drip set with needle-pierceable Y site	84 inches	60 drops/mL	18 mL
IV-016Y	Set with two needle-pierceable Y sites	100 inches	10 drops/mL	23 mL
IV-019	Vented secondary set	40 inches	20 drops/mL	10 mL
IV-020Y	Vented standard set with needle-pierceable Y site	84 inches	20 drops/mL	18 mL
IV-025Y	Vented set with .22 micron filter and needle-pierceable Y site	92 inches	20 drops/mL	22 mL
IV-F02Y	Add-on .22 micron filter set and needle-pierceable Y site	21.5 inches	N/A	6.4 mL

For product codes, box and case sizes please see page 15


One Injection Site Sets




IV-009N Standard set with NAC®



IV-009Y Set with needle-pierceable Y-site




IV-013Y Micro-drip set with needle-pierceable Y-site



IV-020Y Vented standard set with needle-pierceable Y-site

Two Injection Site Sets




IV-016Y Set with two needle-pierceable Y-sites

Secondary Set



IV-019 Vented secondary set

.22 Micron Filter Sets



IV-025Y Vented set with .22 micron filter and needle-pierceable Y-site  
 € 0123



IV-F02Y Add-on .22 micron filter set and needle-pierceable Y-site  
 € 0123

NAC® (Needleless Access Connector) is a registered trademark of CareFusion

For product codes, box and case sizes please see page 15

## Ordering Information

## Winged Blood Collection Sets

(50 per box / 4 boxes per case)



Page 2



Page 4



Page 4

Code#			Needle Gauge	Wing Color	Needle Length	Tubing Length
DBMA-21G	DBMS-21G	DBM1-21G	21g	Green	3/4"	12"
DBMA-23G	DBMS-23G	DBM1-23G	23g	Light Blue	3/4"	12"
DBMA-25G	DBMS-25G	DBM1-25G	25g	Dark Blue	3/4"	12"

With Pre-attached Tube Holder <small>Page 3</small>						
(50 per box / 2 boxes per case)						
Code#			Needle Gauge	Wing Color	Needle Length	Tubing Length
DBMH-21G			21g	Green	3/4"	12"
DBMH-23G			23g	Light Blue	3/4"	12"
DBMH-25G			25g	Dark Blue	3/4"	12"



Page 5

## K-Shield Advantage® Small Vein Infusion Sets

(50 per box / 10 boxes per case)

Code#	Needle Gauge	Wing Color	Needle Length	Tubing Length	Priming Volume
D3KA-19G	19g	Brown	3/4"	12"	.34mL
D3KA-21G	21g	Green	3/4"	12"	.34mL
D3KA-23G	23g	Light Blue	3/4"	12"	.34mL
D3KA-25G	25g	Dark Blue	3/4"	12"	.34mL



Page 5

## Standard Small Vein Infusion Sets

(50 per box / 10 boxes per case)

Code#	Needle Gauge	Wing Color	Needle Length	Tubing Length	Priming Volume
D3K2-19G	19g	Brown	7/8"	12"	.34mL
D3K2-21G	21g	Green	3/4"	12"	.34mL
D3K2-23G	23g	Light Blue	3/4"	12"	.34mL
D3K2-25G	25g	Dark Blue	3/4"	12"	.34mL
D3K1-21G	21g	Green	3/4"	3.5"	.11mL
D3K1-23G	23g	Light Blue	3/4"	3.5"	.11mL
D3K1-25G	25g	Dark Blue	3/8"	3.5"	.11mL



# Ordering Information

## Intermittent Injection Site

Page 5

(50 per box / 4 boxes per case)

Code #	Description	Length
INT-01	Injection site with male luer connector	

## K-Shield Advantage® Blood Drawing Kits

Page 6

(25 per box / 2 boxes per case)

Code #	Description	Length
BK-64A	K-Shield Advantage® 600 mL blood bag with 48 inches 48" tubing, pinch clamp and pre-attached 16G venipuncture needle with safety device, and needle-free sampling port	48 inches

## Standard Phlebotomy Sets

Page 7

(25 per box / 2 boxes per case)

Code #	Description	Length
BO-1000	Phlebotomy set with 16G needle vacuum bottle	36 inches
BO-16S1	Phlebotomy set with 16G vacuum bottle needle and pre-attached 16G venipuncture needle with backeye	36 inches



## K-Shield Advantage® Phlebotomy Sets

Page 7

(25 per box / 2 boxes per case)

Code #	Description	Length
BA-16S1	K-Shield Advantage® Phlebotomy set with 16G vacuum bottle needle and pre-attached 16G venipuncture needle with backeye and safety device	36 inches

## Ordering Information

			Page 8  (20 per box / 5 boxes per case)		Page 9  (25 per box / 4 boxes per case)	
Needle Gauge	Needle Length	Pinch Clamp Color	Microbore Tubing		Microbore Tubing	
			Without NAC®	With NAC®	Without Y-SITE	With Y-SITE
19	3/4"	Brown	PS03090	PS01N90	PI03090	PI01Y90
19	1"	Brown	PS03091	PS01N91	PI03091	PI01Y91
20	3/4"	Yellow	PS03000	PS01N00	PI03000	PI01Y00
20	1"	Yellow	PS03001	PS01N01	PI03001	PI01Y01
22	3/4"	Black	PS03020	PS01N20	PI03020	PI01Y20
22	1"	Black	PS03021	PS01N21	PI03021	PI01Y21
20	1 1/2"	Yellow	-	-	-	PI01Y03
			Standard Tubing		Standard Tubing	
19	1"	Brown	PS04091	-	PI04091	-
20	3/4"	Yellow	PS04000	-	PI04000	-
20	1"	Yellow	PS04001	-	PI04001	-

NAC® (Needleless Access Connector) is a registered trademark of CareFusion

# Ordering Information



Gravity IV Administration Sets <span style="float: right;">Page 10, 11</span>				
(50 sets per case)				
Code #	Description	Length	Drip Rate	Priming Volume
IV-009N	Standard set with NAC®	84 inches	10 drops/mL	18 mL
IV-009Y	Set with needle-pierceable Y-site	84 inches	10 drops/mL	18 mL
IV-013Y	Micro-drip set with needle-pierceable Y-site	84 inches	60 drops/mL	18mL
IV-016Y	Set with two needle-pierceable Y-sites	100 inches	10 drops/mL	23 mL
IV-019	Vented secondary set	40 inches	20 drops/mL	10 mL
IV-020Y	Vented standard set with needle-pierceable Y-site	84 inches	20 drops/mL	18 mL
IV-025Y	Vented set with .22 micron filter and needle-pierceable Y-site	92 inches	20 drops/mL	22 mL
IV-F02Y	Add-on .22 micron filter set and needle-pierceable Y site	21.5 inches	N/A	6.4 mL

NAC® (Needleless Access Connector) is a registered trademark of CareFusion