

ENTERAL FEEDING IRRIGATION KITS AND TRAYS

AMSure® Enteral Feeding Irrigation Kits and Trays are designed for simplicity and efficiency of care.

Kits are non-sterile and include a patient label, 500 mL graduated container, either a poly pouch or re-sealable bag and a 60-mL syringe with or without an ENFit® tip.

Trays are sterile and include 500 mL graduated container, 1000 mL tray, alcohol prep pad and moisture-proof under pad, a 60-mL syringe with or without an ENFit® tip.

- ENFit® connector sets are compliant with ISO 80369-3 to help reduce risk of misconnections
- Latex-free, DEHP-free and BPA-free
- Includes syringe tip protector cap

AS121



AS136

KITS

ITEM CODE	ENFIT® SYRINGE TIP	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL RING PISTON SYRINGE	POLY POUCH	RE-SEALABLE POLE BAG	PATIENT LABEL	QUANTITY/ UNIT
AS121		•	•				•	30/Case
AS126		•		•			•	30/Case
AS127		•		•		•	•	30/Case
ENK126	•	•		•	•		•	30/Case
ENK127	•	•		•		•	•	30/Case

TRAYS

ITEM CODE	ENFIT® SYRINSGE TIP	1000 mL OUTER TRAY	500 mL GRADUATED CONTAINER	60 mL BULB SYRINGE	60 mL THUMB CONTROL RING PISTON SYRINGE	STERILE	ALCOHOL PAD	LARGE UNDERPAD	QUANTITY/ UNIT
AS130		•	•	•		•	•	•	20/Case
AS136		•	•		•	•	•	•	20/Case
ENK136	•	•	•		•	•	•	•	20/Case