

TENS - Transcutaneous Electrical Nerve Stimulator



BioStim® M7 w/ SD (Strength Duration)

- Best in Class, Digital lightweight and functional dual-channel TENS.
- 7 programmable modes including Strength Duration (SD1 & SD2).
- 1% tolerance means the unit holds the integrity of the settings chosen by therapist.
- Patient lock & compliance meter, timer.
- Unit comes with four electrodes, 48" lead wires, two "AA" batteries, carrying case.
- 5-year warranty (USA only)

Prescription Required	Part# KBSM7	HCPCS -E0730	EA
-----------------------	-------------	--------------	----



BioStim® SD (Strength Duration)

- Extremely durable and portable Digital lightweight, dual-channel TENS.
- 4 modes of operation including Strength Duration (SD1 & SD2).
- Clinician can lock the device in place to ensure the settings are retained.
- Patient lock & compliance meter, timer.
- Unit comes with four electrodes, 48" lead wires, two "AA" batteries, carrying case.
- 5-year warranty. (USA only)

Prescription Required	Part# KBSD	HCPCS -E0730	EA
-----------------------	------------	--------------	----



BioStim® Plus

- Extremely durable dual-channel TENS.
- Dual-channel TENS with a large, 2" LCD screen and digital keypad.
- 4 modes of operation: continuous, width modulation burst and cycled burst.
- Patient lock & compliance meter, timer.
- Unit comes with four electrodes, 48" lead wires, (2) "AA" batteries and carrying case.
- 5-year warranty. (USA only)

Prescription Required	Part# KBSP	HCPCS -E0730	EA
-----------------------	------------	--------------	----




IMPULSE® TENS D5

- Dual-channel with flip-top safety cover.
- 5 programmable modes.
- Patient lock & compliance meter, timer.
- Amplitude resets to zero when lead wires or electrodes detach from patient.
- Unit comes with four electrodes, 48" lead wires, two "AA" batteries, carrying case.
- 3-year warranty. (USA only)

Prescription Required	Part# KIMTD5	HCPCS -E0730	EA
-----------------------	--------------	--------------	----

BioMed® 2000 XL / 2000 XT



- Dual-channel, analog TENS device provides nerve stimulation to ease pain.
- Three modes of operation.
- Features covered amplitude knobs, hinged battery cover and horizontal belt clip for comfortable wear.
- Unit comes with four electrodes, lead wires, 9-volt battery, carrying case and instructions
- 5-year warranty. (USA only)

Prescription Required	# KBXL	HCPCS -E0730	EA
-----------------------	--------	--------------	----

IMPULSE® 3000 / 3000T



- Dual-channel, analog TENS device provides nerve stimulation to ease pain.
- Three modes of operation.
- Features covered amplitude knobs, hinged battery cover and horizontal belt clip for comfortable wear.
- Kit includes four pigtail electrodes, two pairs of lead wires, one 9-volt battery, carrying case and instructions.
- 3-year warranty. (USA only)

Prescription Required	# KIM3	HCPCS -E0730	EA
-----------------------	--------	--------------	----

Prescription Required	# KIM3T (timer)	HCPCS -E0730	EA
-----------------------	-----------------	--------------	----

	Impulse® 3000	BioMed® 2000	Impulse® D5	BioStim® Plus	BioStim® SD	BioStim® M7
Channels	2-Analog	2-Analog	2-Digital	2-Digital	2-Digital	2-Digital
Amplitude (mA)	1 to 98	1 to 98	1 to 98	1 to 98	1 to 98	1 to 98
Pulse Rate (Hz)	2 to 150	2 to 150	1 to 120	1 to 150	1 to 200	1 to 200
Pulse Width (µs)	50 to 250	50 to 250	25 to 250	10 to 250	10 to 250	10 to 250
MODES:	3	3	5	4	4	7
Constant	Yes	Yes	Yes	Yes	Yes	Yes
Burst 1 adjustable	Yes	Yes	Yes	Yes	Yes	Yes
Burst 2 preset	*	*	Yes	Yes	*	*
Width modulation	Yes	Yes	Yes	Yes	*	Yes
Rate modulation	*	*	Yes	*	*	Yes
Rate/Width modulation	*	*	*	*	*	Yes
Strength Duration 1	*	*	*	*	Yes	Yes
Strength Duration 2	*	*	*	*	Yes	Yes