

# PATIENT HANDLING

## Hydraulic Patient Lifter #0087

- Dual purpose, padded swivel bar for added safety
- Can be used with 4 or 6 point sling
- Rotating handle on hydraulic pump
- Adjustable width base
- Lifter only weighs 68 lbs and disassembles for easy one man delivery
- Heavy duty round tube steel construction
- Chip resistant brown vein finish
- Two (2) Locking Casters
- Warranty: 1 year on pump, 3 years on frame

### Technical Specifications

Base Closed	22" (can easily fit through small door openings)
Base Open	Up to 42.25" for patient safety
Max Height	64"
Min Height	22"
Base Height	6.5'
Castors	5"
Weight Capacity	450 lbs.

## Replacement Lifter Pump #0086P

- Warranty: 1 year

## Patient Lifter Chains #5075



## Patient Lifter Slings

- Four point strap hookup
- Mesh or nylon slings are washer safe and can be used in wet or dry environments
- Commode Opening: 8.5"x 11"

### Product Codes

SL110	Medium mesh full body
SL111	Large mesh full body
SL112	Medium nylon full body
SL113	Large nylon full body
SL117	X-large nylon full body
SL114	Medium mesh full body w/ commode opening
SL115	Large mesh full body w/ commode opening
SL116	X-large mesh full body w/ commode opening

