

PG II-Rated Large Volume Containers

The PG II-Rated containers allow for transport of sharps, cultures, and stocks. This rating eliminates the need for costly secondary packaging and provides a solution for safe disposal of sharps, chemotherapy waste, cultures, and stocks.



PG II Sharps Containers and PG II Chemotherapy D.O.T. Compliant

- Multiple sizes (8, 12, 18 gallon) for containment of disposable contaminated sharps and trace chemotherapy waste.
- Properly labeled ID# "UN3291", "PG II" and "Regulated Medical Waste". Eliminates the need for costly D.O.T. compliant secondary packaging during transport of sharps, cultures, and stocks when transported by a private or contract carrier.
- Absorbent pad helps contain excess fluid.
- Gasketed lid provides leak resistance.
- Label differentiates containers from those not D.O.T. compliant.
- Hinged lids or sliding lids allow for temporary and final opening and closing.
- Molded-in carrying handles allow for easier transport.
- SharpsCart™ foot pedal carts and floor brackets provide extra stability and easier transport during use.
- Sliding gasket lids have a round 6" opening to accommodate up to 1000mL bottles.

Code	Description	Dimensions	Ship Case
PG II Sharps Disposal Containers			
8997PG2	8 Gallon, Red, Hinged Lid	17.75"H x 11"D x 15.5"W	10
8932PG2	12 Gallon, Red, Hinged Lid	18.75"H x 12.75"D x 18.25"W	10
8998PG2	18 Gallon, Red, Hinged Lid	26"H x 12.75"D x 18.25"W	5
8997SPG2	8 Gallon, Red, Sliding Lid	17.75"H x 11"D x 15.5"W	10
8936PG2	12 Gallon, Red, Sliding Lid	18.75"H x 12.75"D x 18.25"W	10
8998SPG2	18 Gallon, Red, Sliding Lid	26"H x 12.75"D x 18.25"W	5
PG II Chemotherapy Containers			
8985PG2	8 Gallon, Yellow, Hinged Lid	17.75"H x 11"D x 15.5"W	10
8931PG2	12 Gallon, Yellow, Hinged Lid	18.75"H x 12.75"D x 18.25"W	10
8989PG2	18 Gallon, Yellow, Hinged Lid	26"H x 12.75"D x 18.25"W	5
8985SPG2	8 Gallon, Yellow, Sliding Lid	17.75"H x 11"D x 15.5"W	10
8934PG2	12 Gallon, Yellow, Sliding Lid	18.75"H x 12.75"D x 18.25"W	10
8939PG2	18 Gallon, Yellow, Sliding Lid	26"H x 12.75"D x 18.25"W	5