

Low Activity Liner

The Low Activity Liner provides superior skin conditioning and comfort for low to moderate activity level patients.

- Lined with a medium durometer gel
- Stabilizing Matrix reduces distal elongation and minimizes pistoning
- Liner easily conforms to limbs
- Excellent for sensitive skin
- Retrofittable to Silipos Explorer™ sockets
- Available in 3mm, 6mm, and 9mm gel thickness
- Packaged in a box



Pin Liner, 3mm Gel Thickness

Item #	Size	Style
18101	Medium	Uniform
18107	Medium +	Uniform
18113	Large	Uniform
18116	Large +	Uniform
18119	XLarge	Uniform

Item #	Size	Style
18102	Medium	Tapered
18108	Medium +	Tapered

Pin Liner, 6mm Gel Thickness

Item #	Size	Style
18103	Medium	Uniform
18109	Medium +	Uniform
18114	Large	Uniform
18117	Large +	Uniform
18120	XLarge	Uniform

Item #	Size	Style
18104	Medium	Tapered
18110	Medium +	Tapered

Pin Liner, 9mm Gel Thickness

Item #	Size	Style
18105	Medium	Uniform
18111	Medium +	Uniform
18115	Large	Uniform
18118	Large +	Uniform

Item #	Size	Style
18106	Medium	Tapered
18112	Medium +	Tapered

Low Activity Liner Sizing Guidelines

Size	Proximal Circumference		Distal Circumference	
	Relaxed	Fully Stretched	Relaxed	Fully Stretched
Medium	9" (23cm)	16" (41cm)	8" (20cm)	12" (30cm)
Medium+	13" (33cm)	21" (53cm)	9" (23cm)	12" (30cm)
Large	15" (38cm)	22" (56cm)	9" (23cm)	15" (38cm)
Large+	16" (41cm)	27" (68cm)	12" (30cm)	16" (41cm)
X-Large	16" (41cm)	30" (76cm)	12" (30cm)	20" (51cm)

All sizes are approximate.



Cushion Liner, 3mm Gel Thickness

Item #	Size	Style
18141	Medium	Uniform
18147	Medium +	Uniform
18153	Large	Uniform
18156	Large +	Uniform
18159	XLarge	Uniform

Item #	Size	Style
18142	Medium	Tapered
18148	Medium +	Tapered

Cushion Liner, 6mm Gel Thickness

Item #	Size	Style
18143	Medium	Uniform
18149	Medium +	Uniform
18154	Large	Uniform
18157	Large +	Uniform
18160	XLarge	Uniform

Item #	Size	Style
18144	Medium	Tapered
18150	Medium +	Tapered

Cushion Liner, 9mm Gel Thickness

Item #	Size	Style
18145	Medium	Uniform
18151	Medium +	Uniform
18155	Large	Uniform
18158	Large +	Uniform

Item #	Size	Style
18146	Medium	Tapered
18152	Medium +	Tapered