

solo®

DIGITAL CLINICAL SCALES

DIGITAL PHYSICIAN SCALE WITH HEIGHT ROD



- 550 lb x 0.2 lb / 250 kg x 0.1 kg Capacity
- Weight, Height, and Body Mass Index
- 14 in W x 15 in D x 2.5 in H Platform Size
- 6 AA Batteries (included) or Optional AC Adapter
- Inline Mechanical Height Rod

**THE MOST ECONOMICAL,
HIGH-QUALITY DIGITAL
PHYSICIAN SCALE
ON THE MARKET!**

*Detecto's solo® scale stands
alone in economic clinical
weighing*

The oversized 14 x 15
inch platform with rubber
mat is only 2.4 inches
high making it easy for
patients to step on/off.

*No power outlet available, no
problem! The solo® comes with 6 AA
batteries for weighing anywhere.*



The solo® indicator
features 8 simple buttons
that are easily identifiable
in any language.



*The AC port is
located in the base,
so no pesky wires
dangle from up on
the indicator (AC Adapter
optional: model PD-AC)*

 **DETECTO®**
www.Detecto.com