



Non-Tilt Overbed Table



FIGURE A

- Spring-loaded lift mechanism – positions top in 1" increments from 28" - 45"
- ¾" x 15" x 30" high pressure laminate top with T-molding
- Gunstock walnut is low pressure laminate
- Twin wheel 1.5" casters
- Available with “H” (Figure A) or “U” base in chrome or opal powder coat
- Optional dual-mirror vanity drawer
- Base and columns available in additional powder coat finishes
- Available on drop ship basis only

SPECIFICATIONS	Description
Height (adjustable)	28" - 45"
Weight	24 lbs.
Carton Shipping Weight	27 lbs.
Weight Limit	100 lbs.
Base Dimensions	14" x 24"
Tabletop Surface Dimensions	15" x 30"
Warranty	Limited Lifetime

[Click here to view detailed product descriptions and prices.](#) **Call 800-397-5899 to order today!**