



Urinary Drainage Bags

Latex-Free

Item

Qty/Unit

AS322

20/Case

2000mL, Anti-Reflux Chamber,
Pre-pierced Needle Free Sampling Port
(Luer Slip or Blunt Cannula Compatible),
Universal Double Hook and Rope Hanger,
T-Tap Drain Port, Sterile Fluid Pathway



AS32200

20/Case

2000mL, Low Profile
Anti-Reflux Chamber,
Pre-pierced Needle-Free Sampling Port
(Luer Slip or Blunt Cannula Compatible),
Single Hook and Rope Hanger,
T-Tap Drain Port, Sterile Fluid Pathway



AS32600

20/Case

2000mL, Low Profile,
Anti-Reflux Chamber,
AMSafe® Needle-Free Sampling Port
(Luer Lock or Luer Slip Compatible),
Single Hook and Rope Hanger,
T-Tap Drain Port, Sterile



AS302

20/Case

2000mL, Anti-Reflux Device
Pre-pierced Needle Free Sampling Port
(Luer Slip or Blunt Cannula Compatible),
Universal Double Hook and Rope Hanger,
T-Tap Drain Port, Sterile Fluid Pathway



AS332

20/Case

2000mL, Anti-Reflux Device,
Pre-pierced Needle-Free Sampling Port
(Luer Slip or Blunt Cannula Compatible),
Single Hook and Rope Hanger,
T-Tap Drain Port, Sterile Fluid Pathway