

Live the Longer
Life You Love[®]

5000 BED RAIL ADVANTAGE TRAVELER

Portable Bed Rail With Free Organizer





EASILY GET OUT OF BED



ORGANIZER KEEPS ITEMS CLOSE BY



FOLDS FOR TRAVEL

Features & Benefits

- **Ergonomic Cushion Handle** - Allows for easy transfer in and out of bed.
- **4-Pocket Organizer Included** - Prevents entrapment and provides storage space for handy items.
- **Portable** - Easily folds for travel or storage.
- **Safety Strap** - Secures rail to any home or hospital bed.

Specifications

Rail Height: 22"

Rail Width: 17"

Base Dimensions: 21" x 17"

Weight of product: 8 lbs

Weight Capacity: 400 lbs

Made for Bed Type: Attaches to any home or hospital bed with included Safety Strap

Package Dimensions: 23" x 18" x 2"

Case of 4 Dimensions: 23" x 18" x 8"; 34 lbs

How Do We Compare

Bed Rail Advantage Traveler

Other Bed Rails

No Assembly Required



Easily Folds for Travel



4-Pocket Organizer Pouch Included



Provides Assistance In & Out of Bed



Durable Powder Coat Finish



Easy Assembly (No Tools Required)



Open Bed Rail



Insert Safety Lock Pins



Secure to Bed Frame with Included Safety Strap