

Live the Longer
Life You Love®

2085 BEDCADDIE

Offers A Helping Hand To A Sitting Position In Bed





Features & Benefits

- **3 Hand Grips** – Ladder-like design makes sitting up easy.
- **Versatile** – Attaches to any bed frame.
- **Sturdy nylon construction** – Provides long lasting safety.
- **Adjustable**- Customize length to fit your preference.
- **Ergonomic** – Cushioned, no-slip rubber contour grip.
- **Hassle Free Assembly** – Installs and removes in seconds, no tools required.

Specifications

Strap Length: Adjustable 48" – 84"

Length between handles: 12"

Handle Grip: 5"

Weight Capacity: 300 lbs

Product Weight: 2 lbs

Case of 10 Dimensions: 8" L 8" W 11" H; 10 lbs

Easy Assembly (No Tools Required)

- Loop the lower BedCaddie strap around the bed frame
- Snap the easy connect buckle into place
- Adjust the strap length so that the first handle lies at waist level when lying in bed

