


Live the Longer
Life You Love®



200I COUCHCANE

Offers Support In And Out of a Couch or Recliner



INCLUDED
BONUS



HANDY 4-POCKET ORGANIZER



USE WITH YOUR FAVORITE COUCH OR CHAIR



POUCH KEEPS CONVENIENT ITEMS CLOSE



MAKES STANDING AND SITTING EASY

Features & Benefits

- **Ergonomic Safety Handle** - Makes standing and sitting easy and safe.
- **Saddle Bag Pouch** - Four pockets keep convenient items close by.
- **Adaptable** - Adjusts in length and height to fit either side of your favorite couch or chair. Works great with a lift chair.
- **Super Grip Rubber Pads** - Ensure stability and protect floor from damage.

Specifications

Height Adjustment Handle: 34" to 40"
 Dimensions of Handle: 6" x 6" 1.25 Diameter
 Base Length Adjustment: 20" to 36"
 Weight: 19.5 lbs
 Base Feet Measurements: 4" x 7"
 Weight Limit of Handle: 250 lbs
 Case of 3 Dimensions: 29" L x 19" W x 12" H; 60 lbs

Ways To Use Your CouchCane

Stand with confidence from your favorite couch



Get up easily from your recliner



Add additional support to your lift chair



Easy Assembly (No Additional Tools Required)



Slide Two Base Pieces Together



Insert Handle Into Base



Adjust Depth Of Base And Height Of Handle To Desired



Place Under Your Couch Or Chair