

# BSN medical Wound Care

Your  
Pathway to  
Innovation.

Your Partner  
in Wound  
Healing.

**Cutimed®**  
Closing wounds. Together.

## The assortment for every wound healing phase.

Discover Cutimed® advanced wound care from BSN medical, an innovative wound management range which reliably covers all wound healing phases.

Especially developed for the successful treatment of difficult and complex wounds, our state-of-the-art range of treatment options offers a complete selection of effective, patient- and user-friendly products, as well as wound therapy concepts for chronic wounds such as:

- ▶ **venous leg ulcers**
- ▶ **diabetic foot ulcers**
- ▶ **pressure ulcers**







## Cutimed® Gel

Gentle and effective hydration and debridement








- Pure formulation and preservative-free for less risk of allergic reactions
- Can be covered with many types of secondary dressings (e.g. film dressings) without risk of incompatibility
- Superior moisture donation allows for less frequent dressing changes and effective debridement







Size	REF No.	Pieces/Box	HCPCS
8 g / .3 oz.	7261000	10  + 10 	A6248
15 g / .5 oz.	7261001	10  + 10 	A6248
25 g / .8 oz.	7261002	10  + 10 	A6248

## Iconic wound management

- Wound depth**
-  Superficial
  -  Deep
  -  Superficial + deep

- Wound phase**
-  Necrotic
  -  Infected
  -  Slough
  -  Granulating
  -  Epithelializing

- Exudate level**
-  Low
  -  Low to Moderate
  -  Moderate to High
  -  High to Very High

## Cuticell®

Ointment dressing



- Prevents tissue growth with small mesh and does not disturb wound healing
- Minimizes pain at dressing change with smooth, acetate fabric with knit structures
- Permeable to air and allows exudate to drain freely into an absorbent secondary dressing
- Ointment impregnation helps maintain a moist healing environment


Size	REF No.	Dressings/Box	HCPCS
3 × 3 in.	7253900	10	N/A
3 × 3 in.	7253901	50	N/A
3 × 8 in.	7253902	10	N/A
3 × 8 in.	7253903	50	N/A

W O U N D B E D P R E P

# Cutimed® Sorbact®

## Physical attraction. The Sorbact® Way.

Bacteria and fungi are naturally hydrophobic (i.e., water repellent). The Sorbact® principle takes advantage of this.



## Easily combines with other dressings

Depending on the exudate levels, Cutimed® Sorbact® may be used with amorphous hydrogels, as well as secondary dressings including:

- Transparent film dressings
- Dressings for additional absorptions
- Hydroactive foam dressings

## Dressing Pads



- Absorbent core
- Exudate management

Size	REF No.	Dressings/Box	HCPCS
2.8 × 3.5 in. (7 × 8.9 cm)	7216101	5	A6251
2.8 × 3.5 in. (7 × 8.9 cm)	7216100	40	A6251
4 × 4 in. (10 × 10 cm)	7216201	5	A6251
4 × 4 in. (10 × 10 cm)	7216200	40	A6251
4 × 8 in. (10 × 20 cm)	7216300	20	A6252

## Ribbon gauzes



- Sorbact impregnated ribbon gauze
- Cavity wounds of all kinds

Size	REF No.	Dressings/Box	HCPCS
0.8 × 19.7 in. (2 × 50 cm)	7216600	20	A6266
2 × 78.7 in. (5 × 200 cm)	7216700	10	A6266

## Swabs



- Flat, folded dressing
- Deep or superficial wounds

Size	REF No.	Dressings/Box	HCPCS
1.6 × 2.4 in. (4 × 6.1 cm)	7216401	5	A6222
1.6 × 2.4 in. (4 × 6.1 cm)	7216400	40	A6222
2.8 × 3.5 in. (7 × 8.9 cm)	7216501	5	A6222
2.8 × 3.5 in. (7 × 8.9 cm)	7216500	40	A6222

## Round Swabs



- Sorbact impregnated acetate round swabs
- Deep or superficial wounds

Size	REF No.	Dressings/Box	HCPCS
walnut	7216800	70 (14 × 5)	A6222

## WCL



- Sorbact impregnated acetate fabric
- Shallow wounds and under compression

Size	REF No.	Dressings/Box	HCPCS
2 × 3 in. (5 × 7.6 cm)	7266200	10	A6206
4 × 4 in. (10 × 10 cm)	7266201	10	A6206
4 × 5 in. (10 × 12.7 cm)	7266202	10	A6207
4 × 8 in. (10 × 20 cm)	7266203	10	A6207
6 × 6 in. (15 × 15 cm)	7266204	10	A6207
8 × 8 in. (20 × 20 cm)	7266205	10	A6208

## Cutimed® Sorbact® at work in dry wounds.

The addition of a hydrogel to the Cutimed® Sorbact® acetate dressing allows bacterial management in dry to very lightly draining wounds.

## Gel




- Ready-to-use
- Combines Sorbact dressing with preservative-free hydrogel
- Promotes autolytic debridement
- Slough or partly necrotic wounds

Size	REF No.	Dressings/Box	HCPCS
3 × 3 in. (7.6 × 7.6 cm)	7261100	10	A6231
3 × 6 in. (7.6 × 15 cm)	7261101	10	A6232

# Cutimed® Siltec®

Super absorbent silicone foam dressing



- Allows for longer wear time with superior absorption
- Reduces risk of maceration with silicone wound contact layer, vertical absorption and super absorbent particles
- Promotes better quality of life for patient with pain free removal
- Can be used under compression

## Superior fluid management

### ▶ Smart Pore Structure allows vertical absorption

— reduces the risk of maceration to wound edges and peri-wound skin



Picture shows macerated edges before use




Picture shows true vertical absorption



Picture shows improvement of wound edges after 24 hours


### ▶ Super Absorbent Particles – SAP

— locks exudate away from wound for proven higher absorption and retention



### ▶ Highly Breathable Film for Dynamic MVTR

— transpires enough fluid to keep wound comfortable while maintaining a moist wound environment



**In-vitro tests confirm the outstanding fluid handling of Cutimed® Siltec®.** (SMTL test 2008, data on file)  
Laboratory testing results are intended to illustrate a product's performance under controlled conditions. Actual use results may vary.

## Painfree and atraumatic dressing changes

### ▶ Silicone Contact Layer

- ensures gentle adherence to peri-wound skin without any adhesion to the wound
- protects new layer of epidermal cells at dressing change
- allows high viscose exudate to pass through for absorption within foam structure

# Cutimed® Siltec®

Super absorbent silicone foam dressing



Size	REF No.	Dressings/Box	HCPCS
2 × 2½ in. (5 × 6 cm)	7263200	10	A6209
4 × 4 in. (10 × 10 cm)	7263201	10	A6209
4 × 8 in. (10 × 20 cm)	7263202	10	A6210
6 × 6 in. (15 × 15 cm)	7263203	10	A6210
8 × 8 in. (20 × 20 cm)	7263204	5	A6211



Size	REF No.	Dressings/Box	HCPCS
2 × 2½ in. (5 × 6 cm)	7263000	10	A6209
4 × 4 in. (10 × 10 cm)	7263001	10	A6209
6 × 6 in. (15 × 15 cm)	7263002	10	A6210



Size	REF No.	Dressings/Box	HCPCS
3 × 3 in. (7.5 × 7.5 cm)	7263100	10	A6212
5 × 5 in. (12.5 × 12.5 cm)	7263101	10	A6212
6 × 6 in. (15 × 15 cm)	7263102	10	A6212
7 × 7 in. (17.5 × 17.5 cm)	7263103	5	A6213
9 × 9 in. (22.5 × 22.5 cm)	7263104	5	A6213



Size	REF No.	Dressings/Box	HCPCS
7 × 7 in. (17.5 × 17.5 cm)	7264700	5	A6210
9 × 9 in. (22.5 × 22.5 cm)	7264701	5	A6210




Size	REF No.	Dressings/Box	HCPCS
6½ × 9½ in. (16 × 24 cm)	7264800	5	A6210



Size	REF No.	Dressings/Box	HCPCS
6½ × 9½ in. (16 × 24 cm)	7264801	5	A6210

# Cutisorb® Ultra

## Super absorbent wound dressing



- Absorbs more than 4x the exudate of other wound dressings\*
- Locks exudate away from the wound and surrounding skin even under pressure
- High fluid capacity results in less frequent dressing changes, lower cost of treatment, and allows longer periods of undisturbed healing

Size	REF No.	Dressings/Box	HCPCS
4 x 4 in. (10 x 10 cm)	7263900	20	A6251
4 x 8 in. (10 x 20 cm)	7263901	20	A6252
8 x 8 in. (20 x 20 cm)	7263902	20	A6252
8 x 12 in. (20 x 30 cm)	7263903	20	A6253

# Cutimed® Alginate

## Calcium alginate dressing



- Maintains a moist environment
- Fast gelling
- Highly absorbent
- Retains the exudate in a gel
- Excellent wet integrity

### Cutimed® Alginate compresses

Size	REF No.	Dressings/Box	HCPCS
2 x 2 in. (5 x 5 cm)	7263400	10	A6196
4 x 4 in. (10 x 10 cm)	7263401	10	A6196
4 x 8 in. (10 x 20 cm)	7263402	10	A6197

### Cutimed® Alginate rope

Size	REF No.	Dressings/Box	HCPCS
1 x 11¼ in. (2.5 x 30 cm)	7263403	5	A6199

# Cutimed® Hydro

## Hydrocolloid dressing



- Maintains a moist environment
- Stimulates granulation
- Retains the exudate in a gel
- Waterproof and breathable



### Cutimed® Hydro L



Size	REF No.	Dressings/Box	HCPCS
2 x 10 in. (5 x 25 cm)	7263600	12	A6235
3 x 3 in. (7.5 x 7.5 cm)	7263601	10	A6234
4 x 4 in. (10 x 10 cm)	7263602	10	A6234
6 x 6 in. (15 x 15 cm)	7263603	5	A6235

### Cutimed® Hydro B



Size	REF No.	Dressings/Box	HCPCS
3 x 3 in. (7.5 x 7.5 cm)	7263500	5	A6237
4 x 4 in. (10 x 10 cm)	7263501	5	A6237
6 x 6 in. (15 x 15 cm)	7263502	5	A6238
Sacral 6 x 7 in. (15 x 18 cm)	7263503	5	A6238

# Cutimed® Cavity

## Foam dressing



- Atraumatic
- Non-adhesive
- Soft and conformable
- Easy to apply
- Highly absorbent\*
- Easy to remove




Size	REF No.	Dressings/Box	HCPCS
2 x 2¾ in. (5 x 6 cm)	7266300	10	A6209
4 x 4 in. (10 x 10 cm)	7266301	10	A6209
6 x 6 in. (15 x 15 cm)	7266302	10	A6210

\* Data on file

# Venous leg ulcer management: Therapy approach

Selecting the appropriate compression system



## Comprilan®

Short stretch compression bandage



- Provides high working pressure and low resting pressure
- Can be worn comfortably, safely and effectively during rest and activity
- Can be washed and reused for economical treatment
- Cool and light, the open weave 100% cotton fabric allows the skin to breathe
- Better tolerated if pain is an issue
- Safe for mixed ulcers (0.6 < ABPI < 0.8)

Size	REF No.	Rolls/Box	HCPCS
4 cm x 5 m	7718700	1	N/A
6 cm x 5 m	0102600	1	N/A
8 cm x 5 m	0102700	1	N/A
10 cm x 5 m	0102800	1	N/A
12 cm x 5 m	0102900	1	N/A
10 cm x 10 m	7718800	1	N/A
12 cm x 10 m	7718900	1	N/A

Prevent, Treat and Manage chronic wound vascular conditions

JOBST® and Cutimed®

# JOBST® Comprifore®

## Multilayer compression system



- Provides effective levels of sustained graduated compression
- Provides built in safety and ease of application
- Ensures compliance and maximum healing for cost effective care
- Includes Cutimed® Sorbact® Wound Contact layer to reduce bioburden that may impede healing



### JOBST® Comprifore®

REF No. Quantity  
7266100 1 Set/Box  
4-layer compression bandaging system  
(includes microbe binding wound contact layer)



### JOBST® Comprifore® LF

REF No. Quantity  
7266101 1 Set/Box  
Latex-free, 4-layer compression bandaging system  
(includes microbe binding wound contact layer)



### JOBST® Comprifore® lite

REF No. Quantity  
7266102 1 Set/Box  
3-layer compression bandaging system for reduced compression  
(includes microbe binding wound contact layer)



### JOBST® Comprifore® lite LF

REF No. Quantity  
7266103 1 Set/Box  
Latex-free, 3-layer compression bandaging system for reduced compression  
(includes microbe binding wound contact layer)

Ankle circumference: **Comprifore® and Comprifore® LF** 7 - 10 in. (18 - 25 cm)  
**Comprifore® lite and Comprifore® lite LF** all

Product	#1	#2	#3	#4	Sorbact WCL
JOBST® Comprifore®	•	•	•	•	•
JOBST® Comprifore® LF	•	•	•(LF)	•(LF)	•
JOBST® Comprifore® lite	•	•		•	•
JOBST® Comprifore® lite LF	•	•		•(LF)	•
<b>HCPSCS</b>	A6441	A6443	A6449	A6454	A6207
<b>CPT Code</b>	29581-Application of multi-layer venous wound compression system, below the knee				

# Gelocast®

## The original Unna's Boot

### Gelocast®



- Dries quickly providing a semi-rigid cast
- Soothing calamine reduces skin irritation
- 100% cotton base does not unravel reducing waste
- Provides high working pressure and lower resting pressure

Size	REF No.	Dressings/Box	HCPSCS
3 in. x 10 yds. (8 cm x 9 m)	01052	1	A6456
4 in. x 10 yds. (10 cm x 9 m)	01053	1	A6456

# UlcerCARE™

## Two part compression system



- Easy to don
- Offers handling by patients or caregivers
- Provides high wearing comfort
- Delivers effective compression
- Available as zippered or non-zippered style



Size	Zippered Stocking and Liner		Stocking and Liner		Liner
	Beige Left	Beige Right	Black	Beige	White
S	114485	114520	114510	114479	114500
M	114486	114521	114511	114480	114501
L	114487	114522	114512	114481	114502
XL	114488	114523	114513	114482	114503
2XL	114489	114524	114514	114483	114504
3XL	114532	114530	114506	114484	114505
4XL	114536	114534	114508	114507	114509

**Zippered Stocking/Stocking Package** contains 2 liners per package

**Liner Package** contains 3 liners per package

**Zipper** is placed on the opposite side of the ulcer

**UlcerCARE™ System** - HCPCS Code A6532

**UlcerCARE™ Liner** - HCPCS Code A6530

Size	b	c
S	7" - 8¼"	11½" - 13¾"
M	8¼" - 9½"	13¼" - 16"
L	9½" - 10¾"	15¾" - 18½"
XL	10¾" - 12"	18" - 21¼"
2XL	12" - 13"	19" - 22"
3XL	13" - 14"	19¾" - 23"
4XL	14" - 15¼"	20½" - 23¾"

# Cutimed® Off-Loader Select

Total contact cast system

- 1st choice for wounds not showing significant healing at 4 weeks
- Cost effective
- Reduces inventory
- Complete treatment kit for diabetic foot ulcers
- Reduces shear and friction



## Cutimed® Off-Loader Plaster and Synthetic Combination

REF No.	Description	Qty
7800900	Cutimed® Foam	1 each
	2 × 2 in. Cotton Gauze	4 each
	4 in. Ortho-Flex® Elastic Plaster***	2 rolls
	4 in. Extra Fast Gypsona® S	3 rolls
	3 in. Delta-Lite® Conformable	1 rolls
	4 in. Delta-Lite® Conformable	2 rolls
	3 in. Delta Terry-Net™ Stockinette	1 each
	4 in. Specialist® 100 Cotton Cast Padding	2 rolls
	Delta Terry-Net™ Adhesive Felt	1 each
	Open Cell Foam	1 each

\*\*\*CAUTION: This product contains natural rubber latex which may cause allergic reactions.

## Cutimed® Off-Loader Select Synthetic

REF No.	Description	Qty
7800901	Cutimed® Foam	1 each
	2 × 2 in. Cotton Gauze	4 each
	3 in. Delta-Lite® Conformable	2 rolls
	4 in. Delta-Lite® Conformable	3 rolls
	3 in. Delta Terry-Net™ Stockinette	1 each
	4 in. Specialist® 100 Cotton Cast Padding	2 each
	Delta Terry-Net™ Adhesive Felt	1 each
	Open Cell Foam	1 each

Canvas Rocker Bottom Cast Shoe (not included)

REF No.	Description	Qty
43101408	Canvas Rocker Bottom Cast Shoe – Small	1 each
43101505	Canvas Rocker Bottom Cast Shoe – Medium	1 each
43101602	Canvas Rocker Bottom Cast Shoe – Large	1 each
43101709	Canvas Rocker Bottom Cast Shoe – X-Large	1 each

# Cutimed® PROTECT

Medical skin protection



- Protects wound margins and peri-wound skin
- Preserves and maintains intact skin
- Prevents maceration and peri-wound skin breakdown
- Promotes dressing adhesion and pain free dressing removal
- Available in foam applicator and spray

Size	REF No.	Pieces/Box	HCPCS
<b>Spray</b>			
28 ml / .95 oz.	7265300	1	A6250
28 ml / .95 oz.	7265301	12	A6250
<b>Foam Applicator</b>			
1 ml / .03 oz.	7265400	5	A6250
1 ml / .03 oz.	7265401	25	A6250
3 ml / 1.1 oz.	7265402	5	A6250
3 ml / 1.1 oz.	7265403	25	A6250

# Cutimed® ACUTE

Urea-based skin care



- Absorbs rapidly, penetrates deeply
- Rebalances skin moisture levels
- Re-establishes epidermal barrier function
- Foam mousse is gentle to apply
- Leaves skin smooth and supple
- Does not leave greasy residue
- Matches urea concentration to patient's needs

	Size	REF No.	Pieces/Box	HCPCS
5%	4.23 oz. (125 ml)	7264122	1	A6250
10%	4.23 oz. (125 ml)	7264123	1	A6250



# Leukomed®

## Composite dressing



- Absorbs 12x its weight
- Non-adherent wound pad
- Permeable to air



Dressing Size	Wound Pad Size	REF No.	Dressings/Box	HPCPS
3 1/8 x 4 in.	1 1/2 x 2 3/8 in.	7238001	50	A6254
3 1/8 x 6 in.	1 1/2 x 4 3/8 in.	7238002	50	A6254
4 x 7 1/8 in.	2 x 6 3/8 in.	7238003	50	A6254
4 x 10 in.	2 x 8 3/8 in.	7238004	50	A6254
4 x 11 1/8 in.	2 x 9 3/8 in.	7238005	50	A6255
4 x 13 3/4 in.	2 x 11 1/4 in.	7238006	50	A6255

# Leukomed® T Plus

## Post-op dressing



- Absorbs 12x its weight
- Waterproof and bacteria proof
- Non-adherent wound pad



**Pattern coat technology**  
– ensures high moisture vapor permeability

Dressing Size	Wound Pad Size	REF No.	Dressings/Box	HPCPS
2 x 3 in.	1 1/8 x 1 1/2 in.	7238200	50	A6203
3 1/8 x 4 in.	1 1/2 x 2 3/8 in.	7238201	50	A6203
3 1/8 x 6 in.	1 1/2 x 4 3/8 in.	7238202	50	A6203
4 x 7 1/8 in.	2 x 6 3/8 in.	7238214	50	A6203
4 x 10 in.	2 x 8 3/8 in.	7238203	50	A6204
4 x 11 1/8 in.	2 x 9 3/8 in.	7238204	50	A6204
4 x 13 3/4 in.	2 x 11 1/4 in.	7238205	50	A6204

# Coverlet®

## Extremely absorbent bandages



- Absorbs liquid faster than competitors\*\*
- Absorbs 10x its own weight in water
- 360-degree adhesive keeps wounds clean

Size		REF No.	Dressings/Box	HPCPS
3/4 x 3 in.	Strips	00230	100	A6413
1 x 3 in.		00231	100	A6413
1 x 3 in.		76448000	300	A6413
7/8 in. round	Spots	00301	100	A6413
1 1/4 in. oval		00303	100	A6413
1 1/2 x 2 in.	Patches	00330	100	A6413
2 x 3 in.		00340000	50	A6413
4 x 2 3/4 in.		7801000	50	A6413
1 1/2 x 3 in.	Knuckles	01390	100	A6413
2 1/2 x 1 3/8 in.	Toe Shield (non-sterile)	00801	100	A6413
Sm. finger tip	Digits	01306	100	A6413
Lg. finger tip		01307	50	A6413
3 x 3 in.	4-Wing	00385	50	A6413

\*\*Laboratory testing results are intended to illustrate a product's performance under controlled conditions. Actual use results may vary.



# Leukomed® T

## Transparent film dressing – sterile



- Waterproof and bacteria proof
- Permeable to water vapor and air
- Easy-to-use




**Pattern coat technology**  
– ensures high moisture vapor permeability

Size	REF No.	Dressings/Box	HCPCS
3 × 2 in.	7238100	50	A6257
3 7/8 × 4 in.	7238101	50	A6257
4 × 4 3/4 in.	7238108	50	A6258
4 1/4 × 5 1/2 in.	7238102	50	A6258
6 × 7 7/8 in.	7238109	50	A6258
6 × 10 in.	7238110	50	A6259

# Fixomull® transparent

## Waterproof wide-area fixation dressing – non-sterile



- Used in a wide variety of applications
- Protects from both internal and external factors
- Soil-resistant

Size	REF No.	Dressings/Box	HCPCS
2 in. × 10.9 yds.	7221600	1 roll	A6257
4 in. × 10.9 yds.	7221601	1 roll	A6258
6 in. × 10.9 yds.	7221602	1 roll	A6259

Fixomull® transparent is a non-sterile transparent film. If covering an open wound, a primary dressing needs to be applied first.

# Cover-Roll®

## Adhesive fixation dressing



- Acrylic adhesive dissolves in the presence of alcohol ensuring reduced skin irritation
- Air and moisture permeable reducing the risk of maceration
- Efficient for large wounds as the dressing is firmly secured on all sides
- Radio-transparent

Size	REF No.	Dressings/Box	HCPCS
2 in. × 10 yds.	02034	1 roll	A6448
4 in. × 10 yds.	02035	1 roll	A6449
6 in. × 10 yds.	02041	1 roll	A6450
8 in. × 10 yds.	02042	1 roll	A6450

# Cover-Roll® Stretch

## Conforming, adhesive fixation dressing



- Soft, stretchable non-woven polyester material adapts well to body contours
- Provides light compression without causing constriction or congestion
- The polyacrylate adhesive holds dressing safely and reliably in place
- Radio-transparent

Size	REF No.	Dressings/Box	HCPCS
2 in. × 2 yds.	45547	1 roll	A6448
4 in. × 2 yds.	45548	1 roll	A6449
6 in. × 2 yds.	45549	1 roll	A6450
8 in. × 2 yds.	45550	1 roll	A6450
12 in. × 2 yds.	45551	1 roll	A6450
2 in. × 10 yds.	45552	1 roll	A6448
4 in. × 10 yds.	45553	1 roll	A6449
6 in. × 10 yds.	45554	1 roll	A6450
8 in. × 10 yds.	45555	1 roll	A6450
12 in. × 10 yds.	45556	1 roll	A6450

# Elastomull®

## Elastic gauze bandage



- High density fabric for a soft, rich texture
- Superior performance and stretchability for controlling bleeding and edema
- Precision fit in areas with constant movement
- Helps ensure no slippage

Size	REF No.	Dressings/Bag	HCPCS
<b>Non-sterile</b>			
1 in. × 4.1 yds.	02088000	24 rolls	A4649
2 in. × 4.1 yds.	02089000	12 rolls	A4649
3 in. × 4.1 yds.	02101000	12 rolls	A4649
4 in. × 4.1 yds.	02102000	12 rolls	A4649
6 in. × 4.1 yds.	02103000	6 rolls	A4649
<b>Sterile</b>			
1 in. × 4.1 yds.	02075001	12 rolls	A4649
2 in. × 4.1 yds.	02076001	12 rolls	A4649
3 in. × 4.1 yds.	02070001	12 rolls	A4649
4 in. × 4.1 yds.	02071001	12 rolls	A4649
6 in. × 4.1 yds.	02072001	6 rolls	A4649

TRANSPARENT FILM

FIXATION

# Tensogrip®

## Elastic tubular stockinette



- Easy to apply
- Exerts even pressure and compression
- Washable and reusable
- Wide range of sizes

Size	REF No.	Rolls/Box	HCPCS
1.75 in x 11 yds	7579	1 roll	A6457
2.5 in x 11 yds	7580	1 roll	A6457
2.75 in x 11 yds	7581	1 roll	A6457
3 in x 11 yds	7582	1 roll	A6457
3.5 in x 11 yds	7583	1 roll	A6457
4 in x 11 yds	7584	1 roll	A6457
4.5 in x 11 yds	7585	1 roll	A6457
7 in x 11 yds	7586	1 roll	A6457
8.5 in x 11 yds	7587	1 roll	A6457
12.75 in x 11 yds	7588	1 roll	A6457
2.5 in x 11 yds	7580FL	1 roll	A6457
2.75 in x 11 yds	7581FL	1 roll	A6457
3 in x 11 yds	7582FL	1 roll	A6457
3.5 in x 11 yds	7583FL	1 roll	A6457
4 in x 11 yds	7584FL	1 roll	A6457
4.5 in x 11 yds	7585FL	1 roll	A6457

# Tensoshape®

## Elastic tubular bandage



- Easy to apply
- Anatomical design exerts even pressure and compression
- General light support for treatment of minor soft tissue and injuries

### Tensoshape® – Below the Knee

Size	Limb Measurement	REF No.	Unit of Measure
Small (B/C)	12.5 in x 14 ft	7589	10/case
Medium (C/D)	13.7 in x 15.2 ft	7590	10/case
Large (D/E)	14.8 in x 16.4 ft	7591	10/case

### Tensoshape® – Full Leg

Size	Limb Measurement	REF No.	Unit of Measure
Small (B/C)	12.5 in x 14 ft	7592	10/case
Medium (C/D)	13.7 in x 15.2 ft	7593	10/case
Large (D/E)	14.8 in x 16.4 ft	7594	10/case
Extra Large (D/E)	16 in x 17.6 ft	7595	10/case

Caution: This product contains latex rubber which may cause an allergic reaction.

# SplashCap™

## Wound Irrigation Shield



- Faster Wound Irrigation: Improved productivity and faster turnaround times.
- Safer Built-in Shield: Reduces risk of needle sticks and splattering when using bottle caps without shields.
- Proven Pressure (PSI): 4–15 PSI.
- Universal Fitting: Fits all major standard saline bottles.
- Better Wound Care: Standardize wound irrigation procedures in your department.
- High Volume: Deliver over 250 ml of irrigation fluid in seconds.
- Kid Friendly: No scary syringes.

REF No.	Description	Qty
7446800	SplashCap™ Wound Irrigation Shield	25/case
7446801	SplashCap™ Wound Irrigation Shield	100/case

# AbscessCap™

## Abscess and Deep Wound Irrigation



- Protruding Nozzle – allows penetration below the surface
  - one-handed probing
  - break up loculations
  - irrigate tunneling wounds
- Fits any standard saline bottle – delivers high volumes in a shorter amount of time
- Fits any standard syringe – versatility between preferred irrigation methods
- Built In Splash Shield – permits safe visualization of wound target while reducing splatter
- Large Nozzle Opening – creates a low pressure stream of fluid to safely flush the wound of pus and debris
- Sterile

REF No.	Description	Qty
7446806	AbscessCap™ Irrigation Shield	50 each

SplashCap™ and AbscessCap™ are trademarks of Splash Medical Devices, LLC.



## Training & Education

In addition to the highest quality products and service, BSN medical has an ongoing commitment to ensure a comprehensive range of training and education options based on customer needs. Our programs are designed to optimize educational support so that resources supporting our products are equipped with the latest clinical techniques and cost-effective use. Training is offered in a variety of forms including hands-on workshops, skill labs, manuals, and videos.

Please contact Customer Service or your local BSN medical Sales Rep for more information.

BSN medical, Inc.  
Toll-Free: (800) 552 1157  
Fax: (704) 358 4558  
Email: [BSN.Orders@BSNmedical.com](mailto:BSN.Orders@BSNmedical.com)

## Earn & Learns

BSN medical offers Advanced Wound Care learning centers across the US.

- **Chronic Wound Healing**
  - Cutimed® Sorbact® – Hydrophobic Dressings
- **Venous Leg Ulcers - Hands-On**
  - JOBST® Comprifore® Multi-Layer Compression Bandaging
  - Gelocast® – The Original Unna's Boot
  - JOBST® **UlcerCARE™** (2-in-1 Compression System)
- **Diabetic Foot Ulcers – Hands-On**
  - Cutimed® Off-Loader (Total Contact Casting)

### SUGGESTED HCPCS CODES

This document includes codes that might be used to bill for BSN medical products and related services. Each provider will have to verify the appropriate codes for each patient. It is the provider's sole responsibility to determine and submit appropriate codes, charges and modifiers for services rendered. Providers should contact insurers to verify correct coding procedures prior to submitting claims related to any product or service. BSN medical Inc. cannot guarantee coverage or reimbursement with the codes listed in this billing guide. In all cases, providers will need to follow local payer policies for billing and reimbursement.

**NOTE:** Please refer to the product label and/or package insert for full instructions on safe use of these products.

BSN medical Inc. 5825 Carnegie Blvd. Charlotte, NC 28209-4633  
Tel. (+1) 704 554 9933 Fax (+1) 704 358 4558 [www.bsnmedical.com](http://www.bsnmedical.com)  
To order toll-free: BSN medical (+1) 800 552 1157

