











# Catheter Sizing Continuum



- Catheters sizes measured in French scale.
- Pediatric sizes 6 to 10 French.
- Women sizes 12 to 16 French.
- Men sizes 14 to 22 French.
- 14 French is the most common size.

|   |                                  |
|---|----------------------------------|
|    | ● 6 French, 1.98 mm / 0.078 inch |
|  | ● 8 French, 2.64 mm / 0.104 inch |
|  | ● 10 French, 3.34 mm / 0.13 inch |
|  | ○ 12 French, 4 mm / 0.16 inch    |
|  | ● 14 French, 4.7 mm / 0.185 inch |
|  | ● 16 French, 5.3 mm / 0.21 inch  |
|  | ● 18 French, 6 mm / 0.24 inch    |
|  | ● 20 French, 6.7 mm / 0.265 inch |
|  | ● 22 French, 7.4 mm / 0.295 inch |

French Units are used to measure the outside diameter of a catheter.

One French Unit / 1 French (Fr)

1 Fr = 0.33 mm  
1 Fr = 0.013 inch