


Circumference



HYDROVEN GARMENT SPECIFICATIONS

Hydroven 1, Hydroven 3 and Hydroven FPR Arm Garments

Product #	Product Description	Length	Maximum Circumference	
			Upper Arm	Wrist
510003	Hydroven 3 Pump			
5101A51	Hydroven 1 Half Arm Garment	20"	–	20"
5101A68	Hydroven 1 Full Arm Garment	27"	24"	20"
5101A78	Hydroven 1 Full Arm Garment	31"	24"	20"
5103A68	Hydroven 3 Full Arm Garment	27"	24"	20"
5103A78	Hydroven 3 Full Arm Garment	31"	24"	20"
516003	Hydroven FPR Pump			
5163A68	Hydroven FPR Full Arm Garment	27"	24"	20"
5163A78	Hydroven FPR Full Arm Garment	31"	24"	20"

Hydroven 1, Hydroven 3 and Hydroven FPR Leg Garments

Product #	Product Description	Length	Maximum Circumference		
			Upper Thigh	Calf	Ankle
510003	Hydroven 3 Pump				
5101L50	Hydroven 1 Half Leg Garment	20"	–	20"	20"
5101L66	Hydroven 1 Full Leg Garment	26"	25"	20"	20"
5101L71	Hydroven 1 Full Leg Garment	28"	26"	20"	20"
5101L76	Hydroven 1 Full Leg Garment	30"	28"	20"	20"
5101L84	Hydroven 1 Full Leg Garment	33"	28"	20"	20"
5101L92	Hydroven 1 Full Leg Garment	36"	28"	21"	20"
5103L50	Hydroven 3 Half Leg Garment	20"	–	20"	20"
5103L66	Hydroven 3 Full Leg Garment	26"	25"	20"	20"
5103L71	Hydroven 3 Full Leg Garment	28"	26"	20"	20"
5103L76	Hydroven 3 Full Leg Garment	30"	28"	20"	20"
5103L84	Hydroven 3 Full Leg Garment	33"	28"	21"	20"
5103L92	Hydroven 3 Full Leg Garment	36"	28"	21"	20"
516003	Hydroven FPR Pump				
5163L50	Hydroven FPR Half Leg Garment	20"	–	20"	20"
5163L66	Hydroven FPR Full Leg Garment	26"	25"	20"	20"
5163L71	Hydroven FPR Full Leg Garment	28"	26"	20"	20"
5163L76	Hydroven FPR Full Leg Garment	30"	28"	20"	20"
5163L84	Hydroven FPR Full Leg Garment	33"	28"	21"	20"
5163L92	Hydroven FPR Full Leg Garment	36"	28"	21"	20"

Hydroven Insert Pieces (for all Hydroven garments / max of 2 inserts per garment)


Product #	Product Description	Length	Inches Added to Garment Circumference per Insert				
			Upper Thigh	Calf	Ankle	Upper Arm	Wrist
510LI50	Hydroven Half Leg Insert	20"	–	6"	6"	–	–
510LI66	Hydroven Full Leg Insert	26"	7"	6"	6"	–	–
510LI71	Hydroven Full Leg Insert	28"	7"	6"	6"	–	–
510LI76	Hydroven Full Leg Insert	30"	7"	6"	6"	–	–
510LI84	Hydroven Full Leg Insert	33"	7"	6"	6"	–	–
510LI92	Hydroven Full Leg Insert	36"	7"	6"	6"	–	–
510AI68	Hydroven Full Arm Insert	27"	–	–	–	7"	6"
510AI78	Hydroven Full Arm Insert	31"	–	–	–	7"	6"

Hydroven 1 and Hydroven 3 Arm Garments

**Half arm
measurements;**




**Full arm
measurements:**




- L = Half Arm: Length from fingertips to elbow
- Full Arm: Length from fingertips to 5 cm below the armpit with arm straight
- C = Area of the limb with the greatest circumference at points illustrated
- H = Circumference at hand level

Hydroven 1 and Hydroven 3 Leg Garments

**Half leg
measurements;**



**Full leg
measurements:**



- L = Half Leg: Length from heel to fold in knee
- Full Leg: Length from heel to 5 cm below groin
- C = Area of the limb with the greatest circumference at points illustrated

Hydroven Insert Pieces

Single chamber inflatable expansion inserts increase the circumference of leg or arm garments while providing true circumferential compression. Insert pieces increase the circumference by A cm at the distal position (bottom), by B cm at the proximal position (top), as well as circumferences of the mid thigh, knee, calf and ankle as shown in the table on page 2.

