

## Specifications

### MOISTURE OUTPUT AND TEMPERATURE AT 2 HOURS, 15 BREATHS / MINUTE (ISO 9360)

Vt 250 ml	20.5 mg H <sub>2</sub> O/L at 30.0oC
Vt 500 ml	31 mg H <sub>2</sub> O/L at 30.0oC
Vt 750 ml	29.7 mg H <sub>2</sub> O/L at 30.0oC
Vt 1000 ml	29.6 mg H <sub>2</sub> O/L at 30.0oC

### RESISTANCE TO FLOW (ISO 9360)

30 L/min	0.5 cm H <sub>2</sub> O
60 L/min	1.1 cm H <sub>2</sub> O
90 L/min	2.0 cm H <sub>2</sub> O
Dead space (connectors included - ISO 9360)	60 ml (compressed) 75 ml (extended)
Weight	12 g